



# Sabra

**dental products**

**SabraDent Surgical**

**Technology, Equipment, Supplies**

**PRODUCT CATALOG**



**Warranty Policy:**

Sabra Dental Warrants its products to the original purchaser against defects in materials and workmanship for periods in accordance with those stated throughout the product catalog, provided that the product has been operated and maintained in accordance with the product instruction manual, as recognized by a certified Sabra technician, and has not been subjected to apparent misuse, abuse, or accident.

**Payment Options:**

We accept Visa and MasterCard for payment of invoices within 10 days of the invoice date. After initial purchases with credit card/COD payment, if you would prefer to have open billing terms, please call the office at 1-800-888-4435 to set terms. For accounts with open billing, a 2% monthly interest charge will be applied to all amounts that are past the due date by 30 days. COD shipments have an additional \$11 service charge.

**Orders and Shipping:**

UPS Ground service is our standard shipping method. You may request another method of shipping or service. ie... UPS 2nd Day Air, FedEx Overnight etc. All such requests will be met within reason. We will make every effort to ship orders of any stocking item within two business days from the receipt of the order. At your request, we will ship partial orders if part of your shipment is on back order. If there is a special need, arrangements can be made for express delivery of back order fulfillment at the customers' expense. Orders received after 3PM on any day will not be processed until the next business day.

**Return Policy:**

All products for return must be pre-authorized with a Return Authorization number (RA#). Return packages should be clearly marked with the RA#. No returns will be accepted without an RA number. No returns on Handpieces that have been used - No Exceptions!

**Restocking Fees:**

All returns are subject to a 15% restocking fee, except when a return is accompanied by an order of at least equal value. All products must be returned in their original packaging, in clean condition, with a copy of the invoice and a description of the reason for return or an additional 15% re-packaging and handling fee may be applied. Products must be returned within 30 days of the ship date or will not be accepted for return. Shipping charges for any returned products are the customers' responsibility. No refunds, Company credit only. No returns on custom orders.

**Ordering Information:**

To place an order or for more information: Call Toll Free and speak to our customer service representatives.

Toll Free: 1-800-888-4435

Phone: 631-243-6575

Fax: 631-206-9140

Sabra Dental Products

289 Suburban Avenue

Suite E

Deer Park, NY 11729

Email: sales@sabradent.com

**TABLE OF CONTENTS**

|   |       |
|---|-------|
| Sabra Handpieces . . . . .                              | 2     |
| Implant and Micro Surgery Handpieces . . . . .          | 6     |
| Surgical Units and Electrical Motors . . . . .          | 8     |
| PNEU Micro / Hall Type Drills and Accessories . . . . . | 10    |
| Handpiece Optic System . . . . .                        | 11    |
| LED Lighting System . . . . .                           | 13    |
| Patient Monitoring/Oximetry . . . . .                   | 14    |
| TMJ Therapy Device . . . . .                            | 16    |
| Delivery Units / Parts and Accessories . . . . .        | 17-18 |
| Carbide Burs . . . . .                                  | 19    |
| Sabra Diamonds . . . . .                                | 20    |
| Sabra Polishers . . . . .                               | 23    |
| Diamond Strips . . . . .                                | 24    |
| Specialty Surgical Burs and Kits . . . . .              | 25    |
| Trephines and Mucosa Membrane Punches . . . . .         | 27    |
| Implant Burs and Specialty Parts . . . . .              | 28    |
| Implant Drivers and Tools . . . . .                     | 29    |
| Universal Drivers . . . . .                             | 32    |
| Surgical / Implant Accessories . . . . .                | 33    |
| Bone Mill Instruments and Kits . . . . .                | 34    |
| Extraction System / Expansion System . . . . .          | 35    |
| Bone Expansion and Bone Fixation Kits . . . . .         | 36    |
| Grafting / Allograft Materials . . . . .                | 38    |
| Resorbable Collagen Membranes . . . . .                 | 46    |
| Barrier Membranes / Mesh . . . . .                      | 48    |
| Hemostatic Agents and Wound Care . . . . .              | 49    |
| Dental Models / Case Presentation . . . . .             | 51    |
| Sutures . . . . .                                       | 56    |
| Scalpels . . . . .                                      | 58    |
| Instruments . . . . .                                   | 59    |
| Infection Control . . . . .                             | 61    |
| Impression Products . . . . .                           | 67    |
| Barrier Products . . . . .                              | 69    |
| Consumables . . . . .                                   | 70    |
| Gloves . . . . .  | 71    |
| Apparel and Eyewear . . . . .                           | 72    |
| Hygiene . . . . .                                       | 74    |
| Impression Products and Organizers . . . . .            | 75    |
| Endodontics . . . . .                                   | 78    |
| Amalgams . . . . .                                      | 79    |
| Materials and Cosmetic Dentistry . . . . .              | 80    |
| Turbines and Replacement Parts . . . . .                | 82    |
| Repair Services . . . . .                               | 83    |
| OSHA Compliance . . . . .                               | 84    |




Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Sabra Handpieces

## High Speed



### The Sabra 105° Series Fiber Optic Handpiece

#### ACCESS, VISION, TORQUE

The only handpiece available today with the new 105° angle head for better access, visibility and the large head for torque. The quietest handpiece in dentistry at 54.3db

- Large-head turbine for high-torque, smooth, concentric cutting
- Push button bur changing
- Straight body, no contra angle
- 105° head angle
- Operates at a quiet 54.3 decibels
- 360° swivel on a quick disconnect coupling

ORDER CODE SH00-H9 SABRA 105° FO - \$789.00

### Sabra Swivel Couplers



C-4



C-5



C-6

ORDER CODE C-4 Four hole non-retraction quick connect swivel coupler - \$139.00

ORDER CODE C-5 Five hole non-retraction ISO B optic quick connect swivel coupler - \$149.00

ORDER CODE C-6 Six pin lux type non-retraction ISO C optic quick connect swivel coupler - \$159.00

ORDER CODE C-6 LED Six pin lux type non-retraction ISO C optic quick connect swivel coupler with 1000 hour LED lamp - \$249.00

ORDER CODE BLB-016 Replacement 6 PIN BULB, for Sabra C-6 coupler - \$39.00

ORDER CODE BLB-LED Replacement L.E.D. 6 PIN BULB, for Sabra C-6 coupler - \$119.00



## 45 Degree Angle Handpieces




Photo courtesy of Dr. Jin Kim, DMD  
Garden Grove, CA

### Sabra OMS-45 Series Handpieces

*For the best access to 3rd molars.*

**The gold standard for high speed surgical drills**

#### Sabra OMS 45-II

High speed, quick connect, rear exhaust handpiece, non optic, torque head version

ORDER CODE SH00-H5 SABRA OMS 45-II - \$549.00

#### Sabra OMS 45-F0

High speed, quick connect, rear exhaust handpiece, optic, torque head version

ORDER CODE SH00-H6 SABRA OMS 45-F0 - \$689.00

#### Sabra 45-LS

High speed, quick connect, rear exhaust handpiece, optic, mini head version

ORDER CODE SH00-H7 SABRA 45-LS - \$689.00

## Electric E-Type 45 Degree Angle Handpiece Attachment



### Sabra Electric 45 degree angle E type handpiece

**The worlds first 45 degree electric contra angle handpiece**

- Power-unsurpassed
- Superior accessibility
- Versatile - worlds first two-way water spray function
- Can be used with quatro port mist spray for traditional dental procedures Or can be used with water jet for oral surgery procedures. Takes either 19mm or 25mm surgical bur
- Cellular optic
- Clean head feature prevents contaminant build up. Ceramic bearings 25% harder, and 50% lighter than stainless bearings for longer life, and cooler running. Dura coat titanium coating for longer lasting scratch finish, and comfort Solid titanium body. 30% lighter than stainless for reduced weight and better balance in the hand. Anti heat system in the cap to prevent heat build associated with other brand electric handpieces

ORDER CODE - C1064 - \$1499.00



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.



## High Speed



### Sabra Econotorque

- 350,000 - 400,000 RPM
- Autoclavable
- Extremely high torque
- Smooth, quiet operation
- Light weight

#### ORDER CODE

SH00-H10 Econotorque - \$189.00

SH00-H11 Econotorque Push Button - \$219.00

## Low Speed



### Sabra V Low Speed

#### Universal Straight Low Speed Handpiece

- Smooth vibration free operation, up to 5500 RPM
- 360 swivel superb maneuverability
- Compatible with all standard doriet style disposable or metal angles, ISO stones, cutting instruments, contra angles
- Lightweight @ 4.3 oz
- It is ideal for the Hygienist, but is applicable for any clinical low speed procedure

ORDER CODE SH00L5 - \$389.00

( burs and attachments not included )

### Sabra Prophy Gem

#### The Lightest Swivel Handpiece Today 2.3 oz (60 grams)

- Friction grip chucking mechanism
- Smooth, vibration free operation up to 4000 rpm
- Ergonomic design provides perfect balance
- 360 swivel for superb maneuverability
- Compatible with all standard doriet style disposable or metal angles
- Made in the U.S.A.

ORDER CODE SH00L6 - \$379.00

#### Additional Nose Cone

ORDER CODE SH00L6A - \$99.00



| Brand            | Weight In oz. | Length In Inches |
|------------------|---------------|------------------|
| Sabra Prophy Gem | 2.3           | 3.9              |
| NSK Hygiene Pro  | 2.5           | 4.0              |
| Kavo Prophy Wiz  | 2.7           | 4.3              |
| Midwest RDH      | 3.0           | 4.4              |

\*The NSK, Kavo and Midwest names are not registered Trademarks of Sabra Dental.



## E Type Low Speed



### Straight Handpiece

U-Type Nose 40,000 rpm, Ball Bearing  
ORDER CODE SH00-SS1 - \$181.00



### Contra Angle Handpiece

1:1 20,000 rpm for CA Bur  
ORDER CODE SH00-SC2 - \$119.00



### Contra Angle Handpiece

4:1 For Endodontic Twist (Quarter Turn)  
ORDER CODE SH00-C4T - \$179.00



### SC Contra Power Contra Angle

1:1 20,000 rpm for CA Bur, Ball Bearing  
ORDER CODE SH00-C1 - \$139.00



### SRC-E Endo Reduction Contra Angle

20:1  
ORDER CODE SH00-C3 - \$349.00



### Air Motor

20,000 rpm (0.2MPa), ISO 4 Hole (Midwest Type)  
ORDER CODE SH00-L8 - \$329.00



ISO E-Type



### Air Motor

5,000 rpm (0.2MPa), ISO 4 Hole (Midwest Type)  
ORDER CODE SH00-L7 - \$379.00



### Contra Angle Sheath

Sabra SC Power Contra Sheath (No Head)  
ORDER CODE SH00-CH3 - \$69.00



### External Waterclip for SABRA Contra Angles

ORDER CODE SHP-100 - \$36.00

## Spare Heads



20,000 rpm for CA Bur  
ORDER CODE SH00-CH2 - \$49.00



20,000 rpm Miniature Head for CA Bur  
ORDER CODE SH00-CHM - \$49.00



30,000 rpm Ball Bearing for CA Bur  
ORDER CODE SH00-CH1 - \$79.00



30,000 rpm Ball Bearing Slide Latch Type  
for CA Bur  
ORDER CODE SH00-CH1A - \$99.00



30,000 rpm For Endodontic Twist  
(Quarter Turn)  
ORDER CODE SH00-CHT - \$79.00

## Doriet Style Contra Angles



### U-Type Contra Angle

15,000 rpm for CA Bur  
ORDER CODE DC-10L - \$34.00



### U-Type Contra Angle

20,000 rpm for CA Bur Spare Head  
ORDER CODE DC-20L - \$69.00



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.

# Implant Handpieces



### SRC Power Reduction Contra Angle

#### 20:1 Reduction Contra Angle Sheath

High Torque Implant Contra Angle with built in external water clip, and internal water capability.  
**ORDER CODE SH00-C4 - \$549.00**



### SGM Series

- 16:1 MODEL SGM-ER16i ORDER CODE Y110126 - \$779.00
- 20:1 MODEL SGM-ER20i ORDER CODE Y110127 - \$779.00
- 32:1 MODEL SGM-ER32i ORDER CODE Y110128 - \$829.00
- 64:1 MODEL SGM-ER64i ORDER CODE Y110129 - \$829.00
- 256:1 MODEL SGM-ER256i ORDER CODE Y110130 - \$899.00
- 20:1 MODEL SGM-ER20i (with depth indicator) ORDER CODE Y110147 - \$798.00



# Micro Surgery Handpieces



### SGS-ES

MODEL SGS-ES ORDER CODE H264 - \$569.00

- For surgical burs (ø2.35)
- Twist chuck
- Max speed: 40,000 min<sup>-1</sup>

1:1 Direct drive



### SGS-E2S

MODEL SGS-E2S ORDER CODE H266 - \$969.00

- For surgical burs (ø2.35)
- Twist chuck
- Max speed: 80,000 min<sup>-1</sup>

1:2 Increasing



### SGA-ES

MODEL SGA-ES ORDER CODE H263 - \$669.00

- For surgical burs (ø2.35)
- Twist chuck
- Max speed: 40,000 min<sup>-1</sup>

1:1 Direct drive



### SGA-E2S

MODEL SGA-E2S ORDER CODE H265 - \$1069.00

- For surgical burs (ø2.35)
- Twist chuck
- Max speed: 80,000 min<sup>-1</sup>

1:2 Increasing



# Micro Saw Handpieces



1.8mm Reciprocating



**SGR2-E** 3:1 Reduction  
**MODEL SGR2-E ORDER CODE SH162 - \$1199.00**  
 • With external spray nozzle • Max speed:12,600 strokes/min



17° Oscillating



**SG02-E** 3:1 Reduction  
**MODEL SG02-E ORDER CODE SH164 - \$1399.00**  
 • With external spray nozzle • Max speed:12,800 strokes/min



7° sagittal



**SGT2-E** 3:1 Reduction  
**MODEL SGT2-E ORDER CODE SH163 - \$1299.00**  
 • With external spray nozzle • Max speed:12,600 strokes/min

# Blades

\$225.00 **MODEL SGR-1 ORDER CODE Y900072**  
 • Blade length: 1.75" • Pack of 10 blades

\$245.00 **MODEL SGR-2 ORDER CODE Y900073**  
 • Blade length: 1.75" • Pack of 10 blades

\$260.00 **MODEL SGR-3 ORDER CODE Y900074**  
 • Blade length: 1.75" • Pack of 10 blades

**MODEL SGO-1 ORDER CODE H174034**  
 • Blade length: 2.0" • Pack of 10 blades

**MODEL SGO-2 ORDER CODE H174035**  
 • Blade length: 2.0" • Pack of 10 blades

**MODEL SGO-3 ORDER CODE H174036**  
 • Blade length: 2.0" • Pack of 10 blades

**MODEL SGO-4 ORDER CODE H174037**  
 • Blade length: 2.0" • Pack of 10 blades

**MODEL SGO-5 ORDER CODE H174038**  
 • Blade length: 2.0" • Pack of 10 blades

**MODEL SGO-6 ORDER CODE H174039**  
 • Blade length: 2.0" • Pack of 10 blades

All Oscillating Blades - \$40.00 (Except SGO-5, SGO-6 - \$45.00)

\$235.00 **MODEL SGT-1 ORDER CODE Y900075**  
 • Blade length: 1.75" • Pack of 10 blades

\$245.00 **MODEL SGT-2 ORDER CODE Y900076**  
 • Blade length: 1.75" • Pack of 10 blades



**Internal Spray Nozzle**  
**ORDER CODE 5010456U0**  
 For SGM - \$35.00



**Internal Spray Nozzle**  
**ORDER CODE 5010457U0**  
 For SGQ - \$24.00



**External Spray Nozzle**  
**ORDER CODE 5010458U0**  
 For SGP - \$36.00



**External Spray Nozzle**  
**ORDER CODE 5010459U0**  
 For SGP - \$36.00



**Spray Nozzle**  
**ORDER CODE 5010460U0**  
 For SGS, SGA - \$29.00




**Spray Nozzle**  
**ORDER CODE 5010461U0**  
 For SGR, SGT, SGO - \$29.00




**Internal Spray Nozzle Clamp**  
**ORDER CODE 5010462U0 - \$19.00**



**Irrigation Tube Clamp**  
**ORDER CODE 5010463U0 - \$15.00**



**Y-Connector**  
**ORDER CODE 5010464U0 - \$24.00**



**Bur Stopper**  
**ORDER CODE 5013419U0 - \$5.00**



**Bur Stopper Set**  
**ORDER CODE 501232U0 - \$30.00**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
 ©2013 Sabra Dental Products, All Rights Reserved.



# Surgic Pro Complete Sets

Optic / Non-Optic



## Surgic Pro Complete Sets

**Optic Model:** Surgic Pro OPT(120V) - ORDER CODE Y1001949 - \$5399.00

**Optic Model:** Surgic pro OPT-D(120V) - ORDER CODE Y10020953 - \$6499.00

### Contents:

- Control Unit without data storage facility
- SGL70M LED Micromotor
- FC-78 Foot Control
- X-SG20L Optic Handpiece (20:1 Reduction)
- Irrigation tube (5 pcs.) and other accessories

## Surgic Pro Complete Sets

**Non-Optic Model:** Surgic Pro NON-OPT(120V) - ORDER CODE Y1001950 with SG20 Handpiece - \$4399.00

### Contents:


- Control Unit without data storage facility
- SG70M Non-Optic Micromotor
- FC-78 Foot Control
- SG20 Handpiece (20:1 Reduction)
- Irrigation tube (5 pcs.) and other accessories

### Specifications

| Control Unit with AHC  | Micromotor   | Foot Control  |
|--|--|---|
| <ul style="list-style-type: none"> <li>• Power Supply : AC120V or 230V 50/60 Hz</li> <li>• Max. Pump Output : 75mL/min</li> <li>• Programs : 8 Program / Implant System</li> <li>• Dimensions : W265 x D220 x H100 (mm)</li> </ul> | <ul style="list-style-type: none"> <li>• Torque : 5-80Ncm</li> <li>• Motor Speed : 200~40,000 min-1</li> <li>• Light Power : over 32,000 LUX (LED Micromotor)</li> </ul> | <ul style="list-style-type: none"> <li>• Foot Control Functions : Program Button, Speed Control Pedal, Coolant Solution Flow Volume Button, Forward / Reverse Button</li> </ul> |

## Foot Control

### Handle Set (Optional)



### Stay in Control

The Foot Control is user friendly and allows operation of all functions within the preset parameters without touching the control panel to avoid accidental activation of the micromotor outside the preset limits. The Surgic Pro/Surgic Pro+ is certified according to IPX8.

**Handle Set (Optional)** Handle Set can be easily attached to Foot control (FC-78).

**Handle Set - ORDER CODE Z1027001 - \$138.00**

**MODEL: FC-78 - ORDER CODE Z1102001 - \$1100.00**

- With 2m cord

## Sterilization Cassette (Optional)



The NSK Surgic Sterilization Cassette can accommodate the micromotor with cable, straight and contra-angles

**MODEL: SG-CASE**

**ORDER CODE 20001348 - \$250.00**

- Dimensions : W281 x D171.5 x H47 (mm)

## Carrying Case (Optional)



The NSK Carrying Case can accommodate all Surgic Pro components as well as the optional sterilization cassette.

**MODEL: Carrying Case (Surgic Pro)**

**ORDER CODE Y1001952 - \$350.00**

- Dimensions : W469 x D344 x H180 (mm)



# VarioSurge



## Applications

- Bone surgery
- Sinus lift
- Periodontal surgery
- Endodontic surgery
- Orthodontics

## Advantages

- Power:** Superior cutting capability shortens treatment time.
- Visibility:** Optic handpiece assures safe, easy surgical procedures.
- Minimally Invasive:** Eliminates heat generation during bone dissection.
- Safe:** Avoids injury to soft tissue even during accidental contact.
- Precision:** Cavitation effect of irrigation solution keeps surgical site blood free.

NSK's VarioSurg offers numerous compelling advantages, the first of which is strong, precise cutting power. It not only employs ultrasonic cavitation of the irrigation solution but also features an industry first optic handpiece to provide enhanced visibility of the surgical field. This ultrasonic approach to surgery reduces heat generation to minimize osteonecrosis and avoid damage to surrounding soft tissue. VarioSurg does this through ultrasonic micro vibrations which selectively cut only mineralized structures, even in the case of accidental contact. When an area requires special attention, precise cutting is assured, rapidly and without stress, while bone can be treated mechanically with VarioSurg's ultrasonic surgical techniques. NSK VarioSurg — new solution of osteotomy and osteoplasty.



## VarioSurg Ultrasonic Tips

**Optic ORDER CODE Y141395 - \$8999.00**    **Non Optic ORDER CODE Y1001195 - \$8399.00**

| Bone Surgery        |   |
|---------------------|---|
|                     | <b>\$170.00</b><br>MODEL SG1    ORDER CODE Z305-101<br>• Under mark 3mm and 5mm from tip    • Flat bevel, 0.5mm thick |
|                     | <b>\$170.00</b><br>MODEL SG2R    ORDER CODE Z305-102<br>• Right curved tip    • Flat bevel, 0.5mm thick               |
|                     | <b>\$170.00</b><br>MODEL SG2L    ORDER CODE Z305-103<br>• Left curved tip    • Flat bevel, 0.5mm thick                |
| Scraper             |   |
|                     | <b>\$122.00</b><br>MODEL SG3    ORDER CODE Z305-104<br>• Tip: sanded edge available to                                |
|                     | <b>\$122.00</b><br>MODEL SG4    ORDER CODE Z305-105<br>• Edge available to  |
|                     | <b>\$122.00</b><br>MODEL SG5    ORDER CODE Z305-106<br>• Rounded edge available to                                    |
| Sinus Membrane Tool |   |
|                     | <b>\$95.00</b><br>MODEL SG9    ORDER CODE Z305-110<br>• Flat circular convex blade    • Angled at 90°                 |
|                     | <b>\$95.00</b><br>MODEL SG10    ORDER CODE Z305-111<br>• Flat circular convex blade    • Angled at 135°               |
|                     | <b>\$165.00</b><br>MODEL SG11    ORDER CODE Z305-112<br>• Cone compressor   |

| Sinus Lift           |  |
|----------------------|--|
|                      | <b>\$150.00</b><br>MODEL SG6D    ORDER CODE Z305-107<br>• Diamond coating    • Average 3mm from tip                    |
|                      | <b>\$150.00</b><br>MODEL SG7D    ORDER CODE Z305-108<br>• Diamond coating    • Average 3mm from tip                    |
| Peric (Root Planing) |  |
|                      | <b>\$135.00</b><br>MODEL P20-S    ORDER CODE Z305-114<br>• Straight tip  |
|                      | <b>\$135.00</b><br>MODEL P25R-S    ORDER CODE Z305-115<br>(From above)    • Right curved tip                           |
|                      | <b>\$135.00</b><br>MODEL P25L-S    ORDER CODE Z305-116<br>(From above)    • Left curved tip                            |
| Retrograde Endo      |  |
|                      | <b>\$150.00</b><br>MODEL E30RD-S    ORDER CODE Z305-118<br>• Diamond coated    • For posterior teeth (right angled)    |
|                      | <b>\$150.00</b><br>MODEL E30LD-S    ORDER CODE Z305-119<br>• Diamond coated    • For posterior teeth (left angled)     |
|                      | <b>\$150.00</b><br>MODEL E31D-S    ORDER CODE Z305-120<br>• Diamond coated    • For anterior and posterior teeth (70°) |
|                      | <b>\$150.00</b><br>MODEL E32D-S    ORDER CODE Z305-121<br>• Diamond coated    • For anterior teeth (90°)               |

| Scaling   |   |
|---|---|
|   | <b>\$85.00</b><br>MODEL G1-S    ORDER CODE Z305-113   |
| Maintenance (V-Tip)   |   |
|   | <b>\$99.00</b><br>MODEL V10-S    ORDER CODE Z305-117<br>• Plastic Tip (not included)                                      |
|   | <b>\$42.00</b><br>MODEL V-P10    ORDER CODE Y900-184<br>• Pack of 3   |
| Tip Kits  |   |
| <b>\$745.00</b><br>MODEL Basic-S Kit    ORDER CODE Y900-687<br>Contents:<br>• SG1, SG3, SG5, SG6D, SG7D, SG11<br>• Tip holder       | <b>\$715.00</b><br>MODEL Bone Cut Kit    ORDER CODE Y900-688<br>Contents:<br>• SG1, SG3R, SG4, SG4L, SG6D<br>• Tip holder |
| <b>\$660.00</b><br>MODEL Sinus Lift Kit    ORDER CODE Y900-689<br>Contents:<br>• SG1, SG3, SG6D, SG2L, SG10, SG11<br>• Tip holder   |   |
| <b>\$610.00</b><br>MODEL Endo-S Kit    ORDER CODE Y900-691<br>Contents:<br>• G1-S, E30RD-S, E30LD-S, E31D-S, E32D-S<br>• Tip holder |   |



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
 ©2013 Sabra Dental Products, All Rights Reserved.

# PneuMicro Power Systems and Accessories



The Pneu-Micro Small Bone Power System available from Sabra offers a wide range of drills and saws for a variety of surgical procedures. Fully compatible with Linvatec® Hall® pneumatic system, Micro 100™ and Surgairtome Two®, the PneuMicro® is engineered and tested to Sabra standards of high-quality and will be an outstanding addition to your facility.



### Bur Guards

Medium

**Recommended** Bur Length 44.5mm - 69mm

**ORDER CODE 5010873M0 - \$359.00**



### Bur Guards

Long

**Recommended** Bur Length 70mm - 94mm

**ORDER CODE 5010874M0 - \$450.00**



### Cleaning Cap

**Connection** Fits PneuMicro™ Drills & Saws and Wiredrivers (-200)

**ORDER CODE 5007878M0 - \$78.52**



### PneuMicro™ High Speed Drill PM-M09-200/1

**ORDER CODE 5010867M0 - \$2950.00**



### Air/Nitrogen Hose

|                   |                      |                      |
|-------------------|----------------------|----------------------|
| <b>ORDER CODE</b> | 5007877M0 - \$563.00 | 5011371M0 - \$675.00 |
| <b>Length</b>     | 10 ft (3m)           | 15 ft (5m)           |

***We Buy, Refurbish and Sell Used Hall and Stryker Drills, and systems. If you have a system that you want to sell -please call us, or if you want to buy please look at [www.sabradent.com](http://www.sabradent.com) for inventory clearance items or call a professional service technician at 1-800-888-4435***

***Send in your broken old #296-4 Stryker Cord and we will exchange it with a new refurbished one with a 6 month warranty - \$665.00***

Hall is not a registered trademark of SABRA.



Call Toll Free at 1-800-888-4435



# Sabra Brite-Lite Plus Optic System with Integrated Curing Light System

## Sabra Dental Brite Lite Plus Sabra Brite 4 Circuit Universal Handpiece Light Source

Plus - High Intensity Integrated Rapid Cure L.E.D. Curing circuit All in one system that integrates into any dental unit to run any type/brand optic handpiece.



### BRITE-LITE- Four circuit universal handpiece light source

- Voltage controllable for use with any brand hoses / couplers
- Use any brand handpiece / coupler with 3 styles of Sabra hoses:
  - Sabra-Quick 360 quick connect hose for use with Sabra handpieces
  - Sabra ISO-B 5 hole tubing for use with any 5/6 hole handpiece
  - Sabra ISO-C 6 pin tubing for use with any 6 pin handpiece
- Can be retro-fit in most units
- Crisp, bright light
- Long lasting, inexpensive lamp module

### PLUS- High Intensity Integrated Rapid Cure L.E.D. Curing Circuit

- Built into dental unit, integrated into optic system
- Activate with existing foot pedal
- 468-472 Nm at 1000+ mw/cm2 intensity
- "0" heat emission
- Swivel curing 8mm or 11mm probes
- 10 second audible timing tone
- Autoclavable probes
- No recharging - always there and ready!

ORDER CODE BLP-1001 - \$389.00



Circuit 1  
Circuit 2  
Circuit 3  
Circuit 4  
Curing Circuit

### Available Tubing Colors



Light Gray    Gray    Black



Sabra Quick 360 Hose Kit  
ORDER CODE BLHK - 01,02,03 - \$279.00  
( Hose kit includes hose, bulb, and Sabra-Quick coupler )



LED 5 Hole Hose  
ORDER CODE BLH5LED - 10,11,12 - \$219.00



C-5 Coupler



ISO-C 6 Pin  
ORDER CODE BLH6 - 13,14,15 - \$179.00



C-6 Coupler

### Optic Bulb

ORDER CODE BLB-01N - \$32.00



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Sabra Brite-Lite Plus System Parts



## COMPONENTS

**ORDER CODE BLP-1001** Sabra Brite Lite Plus. Led Curing Light And Handpiece Light Source Illumination System (Includes, 1 variable control 4 circuit power pack, 1 low voltage transformer) - **\$389.00**

**ORDER CODE BLC-00** Sabra-Quick 360. Quick connect coupling attachment for Sabra-Quick swivel fiber optic tubing - **\$159.00**

**ORDER CODE BLB-01** Standard fiber optic replacement bulb for Sabra-Quick swivel hoses. - **\$32.00**

## CURING LIGHT PROBE AND ASSEMBLY

**ORDER CODE BLP-01** Sabra Brite Lite Plus Curing/Optics System Led Optics & Probe Assembly - **\$729.00** ( Includes 1x 8mm probe and Lower power assembly )

**ORDER CODE BLP-02** Standard 8mm Curing probe - **\$229.00**

**ORDER CODE BLP-03** Standard 11mm Curing probe - **\$229.00**

# Sabra-Quick Swivel Optic Hoses

## SABRA-QUICK SWIVEL HOSE KITS

**ORDER CODE BLHK-01** Sabra-Quick swivel lite grey hose kit- Includes 1 lite grey swivel optic tubing, 1 Sabra-Quick 360 coupler, 1 BLB-01 swivel optic bulb - **\$279.00**

**ORDER CODE BLHK-02** Sabra-Quick swivel black hose kit- Includes 1 black swivel optic tubing, 1 Sabra-Quick 360 coupler, 1 BLB-01 swivel optic bulb - **\$279.00**

**ORDER CODE BLHK-03** Sabra-Quick swivel grey hose kit- Includes 1 grey swivel optic tubing, 1 Sabra-Quick 360 coupler, 1 BLB-01 swivel optic bulb - **\$279.00**

## SABRA-QUICK SWIVEL HOSES

**ORDER CODE BLH-04** Sabra-Quick swivel optic tubing, lite grey - **\$179.00**

**ORDER CODE BLH-05** Sabra-Quick swivel optic tubing, black - **\$179.00**

**ORDER CODE BLH-06** Sabra-Quick swivel optic tubing, grey - **\$179.00**

## ISO B LED 5 HOLE HOSES:

**ORDER CODE BLH5LED-10** Replacement ISO B, 5 hole LED tubing, lite grey - **\$219.00**

**ORDER CODE BLH5LED-11** Replacement ISO B, 5 hole LED tubing, black - **\$219.00**

**ORDER CODE BLH5LED-12** Replacement ISO B, 5 hole LED tubing, grey - **\$219.00**

## ISO C 6 PIN HOSES

**ORDER CODE BLH6-13** Replacement ISO C, 6 pin optic tubing, lite grey - **\$179.00**

**ORDER CODE BLH6-14** Replacement ISO C, 6 pin optic tubing, black - **\$179.00**

**ORDER CODE BLH6-15** Replacement ISO C, 6 pin optic tubing, grey - **\$179.00**

## L.E.D. CURING LIGHT HOSES

**ORDER CODE BLLH-16** LED Curing Tubing/Replacement hose - **\$179.00**

**ORDER CODE BLLH-17** LED Curing Tubing/Replacement hose - **\$179.00**

**ORDER CODE BLLH-18** LED Curing Tubing/Replacement hose - **\$179.00**



# Pro-Lite Headlight

## Pro-Lite the Ultimate Headlight

Designed with over 28 years of dentist experience, Pro-Lite is the most cost-effective, ultra-lightweight, and ideal universal head light for all dental/medical treatments and applications. The future of dental/medical lighting, Pro-Lite, replaces the need for clumsy overhead lighting.

### Features and Benefits:

- The first only Hands Free Operation Headlight, minimizes cross contamination
- Ultra light weight, for maximum comfort (4 gm only)
- Mobility, time saving and stress reduction
- Ergonomics, reduce arm and shoulder fatigue
- Instant lighting on demand, work efficiently = increase productivity
- Long battery life, last all day
- Reduce patient anxiety, less stress
- Adaptable to all eye wear, keep your favorite loupe or eye glasses
- Energy Saving, cost pennies to operate = less HVAC bill = going green
- Cost effective, equipment for your hygienists and assistants too = more efficient team work
- No more overhead light, technical/technological advancement
- Pays for itself, by saving energy and increasing productivity




The First & Only  
Hands Free Operation Headlight

### Pro-Lite LED systems include:

- 1 Lithium ion battery pack, with easy access ON/OFF button, adjustable intensity control and durable metal clip
- 1 LED head lamp with small adapters to mount on any eye wear
- 1 Orange filter with easy clip adapter
- 1 Battery charger with ready light indicator (110/220 volt)

**PRO-Lite**  
THE ULTIMATE HEAD LIGHT



### Pro-Lite The Ultimate Hands Free Headlight System with Carrying Case

ORDER CODE SPL-01 - \$749.00

### Pro-Lite The Ultimate Hands Free Headlight System

ORDER CODE SPL-02 - \$699.00

### Additional Headlight attachment only

ORDER CODE SPL-03 - \$269.00

### Pro Lite Replacement Transformer DC 4.2v--2.0A

ORDER CODE SPL-04 - \$32.99

# Sabra "White-Lite" Transilluminator

## Sabra White-Lite Portable Cordless Transilluminator/Illuminated Mirror

### Features and Benefits:

- Powerful High Intensity Pure White Light at 6000° K
- Easily Detects Tooth Cracks Increasing Case Awareness
- Can Be Used By a Hygienist To Diagnose Possible Caries
- Vividly Displays Root Canals During Endodontic Procedures
- Perfect for Illuminating bone depth in Sinus Procedures
- Cordless Rechargeable Nickel Hydride Battery Operation
- Light Guide is Easily Removable and Autoclavable
- Modular Battery Connection Permits Spare Battery Use
- Optional One Minute Pre-Programmed Auto-Off Cycle
- Single One-Touch Button Operation for Easy ON/OFF
- Unique Battery Monitor Indicates When Batteries Are Low
- Unique Illuminated mirror option lights up the oral cavity for easy visual inspection, and easy clean up with disposable mirrors



### Sabra White-Lite Portable Cordless Transilluminator/Illuminated Mirror Option

ORDER CODE SIU-01 - \$329.00

### Sabra White-Lite Portable Cordless Transilluminator (no mirror option)

ORDER CODE SIU-01a - \$289.00

### Sabra White-Lite Mirror Add on ( probe and 50 mirrors, barrier sleeves )

ORDER CODE SIU-02 - \$99.00

### Refill Sabra White-Lite Disposable Mirrors, 50 per and 50 barrier sleeves


ORDER CODE SIU-03 - \$29.95



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Medical Oximeters



### Oxi-Go FingerTip Oximeter - The Elite

- Dual Color, High resolution, easy to read OLED display
- Display SpO2, pulse rate, and pulse bar
- 4 Display modes for easy reading from all sides

ORDER CODE SOF-01 - \$69.99



### MD300M Handheld Oximeter

- High resolution 2.8" TFT display
- Adjustable audible and visual alarms
- Store Up To 127 Patient IDs
- 72-hour data storage and review
- Spot clock and continuous use
- Suitable for adult and pediatric greater than 3 years of age
- Auto power-off: Can select the optional auto power-off time from 1 minute to 15 minutes.
- Battery charge function with adapter or charge stand (optional)
- Data Analysis Software: Manage patient data easily on your PC.

ORDER CODE SOH-300M - \$349.99



### MD300A Handheld Oximeter

- High resolution LCD display
- Display SpO2, PR, Plethysmogram, Pulse bar
- 24-hours record storage
- Audible & visual alarm; battery-low indicator
- Rubber jacket providing external protection
- 4 AA Alkaline Batteries

Suitable for adults and pediatrics greater than 3 years of age

ORDER CODE SOH-300A - \$314.99



### ChoiceMMED 300K1 Handheld Oximeter

- Suitable for adult and pediatric patients 3 years of age or older
- Spot check and continuous monitoring modes
- Displays SpO2, Pulse Rate, Pulse Bar
- Up to 99 patient information and 72-hour record storage
- Audio and visual alarm
- Battery-low indicator
- Supports 2 AA alkaline
- Low power consumption
- Digital LED displays

ORDER CODE SOH-300K1 - \$329.00



### Monitor Roll Stand

ORDER CODE PM-60-1 - \$499.00



Call Toll Free at 1-800-888-4435

## Defibrillator



### Cardiac Science Powerheart AED G3 Plus Automatic Package

**System Includes:**

- 7 year Device Warranty
  - Intellisense Lithium Medical Grade Battery w/ FULL 4 YEAR operational guarantee
  - 2 sets Adult defibrillation Electrodes ( 1 pre-connected, 1 Spare )
  - 1 AED Ready Kit that includes one way CPR mask, razor, scissors, 2 sets of gauze, antiseptic wipes, 2 sets of gloves
  - 1 Carrying case for G3 w/ shoulder strap
  - 1 Quick Start ready tool kit w/ Rescue Link and MD link Software w/ 1 serial port communications cable to connect Powerheart to PC ( AED software)
  - 1 AED Prescription provided by Cardiac Science medical Director at NO COST, as required by federal law
- ORDER CODE G3Plus - \$1295.00**

## INR Testing & Supplies



See [www.sabradent.com](http://www.sabradent.com) for refills supplies.

### Hemosense Inratio® PT/INR Monitoring System Professional Kit

**System Includes:**

- Monitor, Carrying Case
- AC Power Supply
- 10 UniStick
- 2 Lancets
- 10 MicroSafe Capillary Tubes
- CLIA Waived

The compact and portable system offers rapid 1 minute PT/INR results at the point of care. From a small fingerstick blood sample, INRatio2 offers proven precision and accuracy. Combined with the only available quantitative and reportable on-board QC, it enables effective and efficient management of anticoagulation patients with confidence and ease of the point of care. Physician Sale Only!

**ORDER CODE HEM 0200431 - \$849.50**

**Lancets, Unistick, 21G, 50/bx**  
**ORDER CODE HEM902009 - \$19.25**

**INRatio Test Strips, Bulk Pack, CLIA Waived, 48/bx (PHYSICIAN SALE ONLY)**  
**ORDER CODE HEM0100139 - \$249.75**

## Glucometer & Supplies



### Bayer Contour USB Meter

- Easy-to-use blood glucose meter with No Coding™ technology
- USB port for quick downloads and trend tracking
- Intuitive pre- and post-meal markers to help with food choices
- Knowledge you can share with your healthcare professional

**ORDER CODE BHD7393M - \$19.99**



### Bayer Contour Next Blood Glucose Test Strips, 25 count

**Innovative Strip Technology**

- Ready to test, right out of the box making daily testing easy.
- Next Generation CONTOUR® NEXT test strips deliver exceptional accuracy for close to professional lab results that give you greater confidence when you test.
- No Coding™ technology eliminates errors due to user miscoding.

**ORDER CODE BHD7310 - \$33.29**



### Duracell Coppertop Battery, Alkaline, Size AA, 24/pk, 6 pk/cs

**ORDER CODE MN1500BKD - \$89.99**



### Duracell Coppertop Battery, Alkaline, Size AAA, 24/pk, 6 pk/cs

**ORDER CODE MN2400BKD - \$89.99**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.



# TMJ Therapy Device

## Passive Motion Jaw Mobilization Device



### E-Z Flex II Air Cushioned TMJ Exerciser

E-Z Flex II is an air cushioned TMJ Exerciser that is designed to restore range of motion and function to the jaw. It is designed specifically for the temporomandibular joint and surrounding muscles. The unique use of flexible tubing and an air pump provides therapeutic passive motion to the temporomandibular joint that is gentle, safe and easy to use with total control by the user.

- Air cushioned force is safe, gentle and flexible
- Safety stops enable controlled progress toward improved jaw opening and function
- User can visually measure progress while using the device providing immediate feedback and control
- Modified syringe hand pump allows user to have hands in a relaxed natural position away from the mouth during use
- Sustained stretch can be achieved without pressing on hand pump
- Jaw muscle strengthening with closure against gentle resistance can be performed by removing elastics
- Unique system of elastics permits control of jaw muscle toning. Removing more elastics creates greater resistance to jaw closure, permitting a controlled progression of toning the muscles of the jaw.

### There are numerous conditions in which E-Z Flex II is indicated:

- Disorders of the Temporomandibular Joint (TMJ)
- Jaw muscle disorders
  - Muscle Spasm
  - Clenching
  - Bruxism (grinding of teeth)
  - Muscle stress & strain
- Tension headaches caused by jaw muscle strain
- Rehabilitation of jaw mobility & function following jaw surgery:
  - TMJ surgery-arthroscopy, arthrotomy
  - TMJ arthrocentesis (joint irrigation)
  - Orthognathic surgery (correction of skeletal jaw disorders)
- Rehabilitation of jaw mobility & function after
  - Trauma (following the healing phase)
- Rehabilitation of jaw mobility & function due to trismus from:
  - Dental/facial infection
  - Dental injection
- Maintenance of jaw mobility and function:
  - During & after radiation therapy to head and neck
  - Certain intra-cranial (brain) surgeries
- Surgi-scaler Add-on (optional)

### E-Z Flex II Air Cushioned TMJ Exerciser

**ORDER CODE EZ-01 - \$395.00**

**Includes:** TMJ Excersizer, Head wrap, 2 cold/heat packs for pre and post therapy, calipers to measure progress.

Please call for wholesale pricing for practitioners who want to re-sell this device direct to patients. If you want to prescribe this medical device to your patients - please call for a prescription pad to be supplied at no charge. ( Most major medical insurance companies are now covering this type of therapy device )



# Sabra Mobile-Aire Portable Cart System



\* implant unit pictured for suggestion, not included with purchase



## Mobile-Aire Cart Now with Bio-Free Tubing

Tubing that prevents Bio-film and contaminant growth within your water lines for up to one year. Change the Pure Tube cartridge and the antimicrobial water bottle annually and stay ahead of the CDC & ADA guideline.

- New EZ-Mount Hose System
- New Self Contained Bio-free Water System (optional)
- 2 Wet Handpiece Lines
- 1 Air/Water Syringe
- Quick Connects to any Air or N2 Schrader Connector
- Ritter or Mayo Tray for Instruments
- Height Adjustable from 30" to 38"
- Built in Fiber Optics (optional)
- New Sturdier "H" Style Frame
- Surgi-scaler Add-on (optional)

## ORDER CODE

### SA-4100 New Sabra Mobile-Aire Portable Flip-Top Cart Unit

2 wet lines, air / water syringe, H style frame w/ EZ-Mount hose system - **\$1995.00**

### SA-4100B Sabra Mobile-Aire Portable Flip-Top Cart Unit

2 wet lines, air / water syringe, H style frame w/ EZ-Mount hose system w/ New Bio-Free Tubing (installed) - **\$2095.00**

### BLP-1000 Sabra Brite Lite Plus

L.E.D. Curing Light and Handpiece Light Source Illumination System (includes: 1 variable control power pack, 1 low voltage transformer) - **\$389.00**

### BLHK-EZ Sabra Brite-Lite EZ Mount Hose Auxiliary Kit

Includes: 1 Sabra EZ-Mount optic tubing, 1 BLB-01 optic bulb, 1 Sabra-Quick coupling - **\$369.00**

### BLHK-EZB Sabra Brite-Lite EZ Mount Hose Bio-Free Auxiliary Kit

Includes: 1 Sabra Bio-Free EZ-Mount optic tubing, 1 BLB-01 optic bulb, 1 Sabra-Quick coupling - **\$369.00**

### BLEZ-4 Sabra EZ Mount 4 Hole Replacement Tubing

Light Grey - **\$62.00**

### BLEZ-4B Sabra EZ Mount Bio-Free 4 Hole Replacement Tubing

Light Grey - **\$72.00**

### BLP-01 Sabra Brite-Lite Plus Curing System

Optic & assembly, EZ-Mount L.E.D. tubing and 8mm probe - **\$889.00**

### BLP-02 / BLP-03 Spare Probe

In either 8mm or 11mm - **\$229.00**

### BR360AB 360 day Pure Tube and Bio-Free Sterisil Antimicrobial bottle - \$229.00

### UC1000 Water quick connect for Piezo / Ultra sonic scaler unit

(includes: tubings, adaptors, fittings & labor) - **\$89.00**

### UC3000 131/2" x 93/4" Plastic Ritter tray for Sabra SA-4100 cart unit - \$12.00

### UC4000 191/8" x 125/8" Stainless steel Mayo tray for Sabra SA-4100 cart - \$67.99

### UC2000 Anti Skid Pad / tray holder- n / c with either tray - \$6.00 without

### 126-008 Foot Control Carrier - \$79.00



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Schrader Style Check Units (Female Couplers)

All Pipe Fittings are National Pipe Threads



**Schrader Style "Y" Assembly**

Air  
ORDER CODE 6252 - \$165.94



**1/8" NPT Male**  
Air  
ORDER CODE 6251 - \$73.70



**1/8" NPT Female**  
Air  
ORDER CODE 6250 - \$73.70



**1/4" Hose Barb**  
Air  
ORDER CODE 6254 - \$73.70



**1/4" NPT Male**  
Air  
ORDER CODE 6253 - \$73.70



**DISS Hand Nut/Tight**  
Air  
ORDER CODE 6256 - \$78.70



**DISS Male**  
Air  
ORDER CODE 6255 - \$73.70



**6202S Connection for N2**  
1/8" NPT male Schrader quick connect valve  
ORDER CODE P-1811-1L - \$49.85



**1/4" Poly QD Male quick disconnect**  
with shut off valve  
ORDER CODE P-1070AS - \$16.00



**Replacement Syringe Button Kit**  
ORDER CODE 113-403 - \$22.00



**Replacement Water Bottle**  
1 ltr  
ORDER CODE P-1811-1L - \$32.00



**Foot Control Carrier**  
Attach to your Sabra cart and get your foot pedal up and off the floor  
ORDER CODE 126-008 - \$79.00



**1/4" Umbilical Tubing**  
for Sabra Cart  
ORDER CODE T-1805 - \$.99 per foot



**Air Water Syringe Tips**  
10 PACK  
ORDER CODE 113-237 - \$52.90




**Replacement Foot Control**  
ORDER CODE SCP-FOOT Standard, Sabra Cart - \$79.00



**Foot Pedal Tubing**  
Per Foot  
ORDER CODE 115-720 - \$1.59




# Sabra Quality Carbide Burs




|         |     |   |   |   |   |
|---------|-----|---|---|---|---|
|         | 1/2 | 2 | 4 | 6 | 8 |
| FG      | .   | . | . | . | . |
| RA      | .   | . | . | . | . |
| HP      | .   | . | . | . | . |
| FG-SU   | .   | . | . | . | . |
| FG-SU-L | .   | . | . | . | . |
| FG-SS   | .   | . | . | . | . |
| HP-L    |     | . | . | . | . |
| HP-XL   |     |   | . | . | . |
| HP-XXL  |     | . | . | . | . |
| HP-S    |     | . | . | . | . |

Round




|        |     |     |      |     |
|--------|-----|-----|------|-----|
|        | 701 | 702 | 702L | 703 |
| FG     | .   | .   | .    | .   |
| RA     | .   | .   | .    | .   |
| HP     | .   | .   | .    | .   |
| FG-SU  | .   | .   | .    | .   |
| FG-SS  | .   | .   | .    | .   |
| HP-L   | .   | .   | .    | .   |
| HP-XL  | .   | .   | .    | .   |
| HP-XXL | .   | .   | .    | .   |
| HP-S   | .   | .   | .    | .   |

Crosscut Fissure  
TAPER




|         |      |      |       |
|---------|------|------|-------|
|         | 1702 | 1703 | 1703L |
| FG      | .    | .    | .     |
| RA      |      |      |       |
| HP      | .    | .    | .     |
| FG-SU   | .    | .    | .     |
| FG-SU-L | .    |      |       |
| FG-SS   |      |      |       |
| HP-L    |      |      | .     |
| HP-XL   |      |      |       |
| HP-XXL  |      |      | .     |
| HP-S    |      |      | .     |

Crosscut Fissure  
TAPERED-ROUNDED HEAD




|        |     |     |     |
|--------|-----|-----|-----|
|        | 557 | 558 | 559 |
| FG     | .   | .   | .   |
| RA     | .   | .   | .   |
| HP     | .   | .   | .   |
| FG-SU  | .   | .   | .   |
| FG-SS  | .   | .   | .   |
| HP-L   |     | .   | .   |
| HP-XL  |     |     |     |
| HP-XXL |     |     |     |
| HP-S   |     | .   | .   |

Crosscut Fissure  
STRAIGHT




|        |     |      |
|--------|-----|------|
|        | 958 | 958C |
| FG     | .   | .    |
| RA     | .   | .    |
| HP     |     |      |
| FG-SU  | .   | .    |
| FG-SS  |     |      |
| HP-L   |     |      |
| HP-XL  |     |      |
| HP-XXL |     |      |
| HP-S   |     |      |

End Cutting




|        |     |    |
|--------|-----|----|
|        | 33½ | 35 |
| FG     | .   | .  |
| RA     | .   | .  |
| HP     | .   | .  |
| FG-SU  | .   | .  |
| FG-SS  | .   | .  |
| HP-L   |     |    |
| HP-XL  |     |    |
| HP-XXL |     |    |
| HP-S   |     |    |

Inverted Cone



|        |     |     |
|--------|-----|-----|
|        | 330 | 331 |
| FG     | .   | .   |
| RA     | .   | .   |
| HP     |     |     |
| FG-SU  |     |     |
| FG-SS  | .   | .   |
| HP-L   |     |     |
| HP-XL  |     |     |
| HP-XXL |     |     |
| HP-S   |     |     |

Pear



|        |     |
|--------|-----|
|        | 245 |
| FG     | .   |
| RA     |     |
| HP     |     |
| FG-SU  |     |
| FG-SS  | .   |
| HP-L   |     |
| HP-XL  |     |
| HP-XXL |     |
| HP-S   |     |

Amalgam Prep  
ELONGATED PEAR

| Model # | Description                                 | Price  |
|---------|---|--------|
| FG      | STANDARD HIGH SPEED BURS (19MM)             | \$1.69 |
| RA      | LATCH BURS                                  | \$1.69 |
| FG-SU   | SURGICAL LENGTH (25MM)                      | \$2.79 |
| FG-SU-L | SURGICAL LENGTH LONG (30MM)                 | \$3.29 |
| FG-SS   | SHORT SHANK BURS (16MM)                     | \$1.69 |
| HP      | LOW SPEED BURS (1 3/4 INCH/44.5MM DIAMETER) | \$1.69 |
| HP-L    | (51MM)                                      | \$3.29 |
| HP-XL   | (65MM)                                      | \$3.29 |
| HP-XXL  | (70MM)                                      | \$3.29 |
| HP-S    | STRYKER J-NOTCH BUR (44.5MM)                | \$3.49 |
| HP-SFS  | STRYKER FLAT SHORT BUR (44.5MM)             | \$3.49 |
| HP-SFL  | STRYKER FLAT LONG BUR (58.22MM)             | \$3.49 |

Burs available in 5 pack or bulk pack



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Sabra Israeli Quality Multi-Use Diamonds



"." in column indicates availability

|                     |     |     |     |     |
|---------------------|-----|-----|-----|-----|
|                     |     |     |     |     |
| Catalog Number      | A-1 | A-2 | A-3 | A-4 |
| Extra Fine (white)  |     |     |     |     |
| Fine Grit (Yellow)  |     |     | .   | .   |
| Medium Grit (Blue)  | .   | .   | .   | .   |
| Coarse Grit (Green) |     |     | .   | .   |

**Ball**

|                     |     |     |     |     |
|---------------------|-----|-----|-----|-----|
|                     |     |     |     |     |
| Catalog Number      | B-1 | B-2 | B-3 | B-4 |
| Extra Fine (white)  |     |     |     |     |
| Fine Grit (Yellow)  | .   | .   |     |     |
| Medium Grit (Blue)  | .   | .   | .   | .   |
| Coarse Grit (Green) |     |     | .   |     |

**Inverted Cones**

|                     |     |     |     |
|---------------------|-----|-----|-----|
|                     |     |     |     |
| Catalog Number      | C-1 | C-2 | C-3 |
| Extra Fine (white)  | .   | .   |     |
| Fine Grit (Yellow)  | .   | .   | .   |
| Medium Grit (Blue)  | .   | .   | .   |
| Coarse Grit (Green) | .   | .   | .   |

**Tapers Round End**

|                     |     |     |     |     |
|---------------------|-----|-----|-----|-----|
|                     |     |     |     |     |
| Catalog Number      | D-1 | D-2 | D-3 | D-4 |
| Extra Fine (white)  |     |     |     |     |
| Fine Grit (Yellow)  |     |     | .   |     |
| Medium Grit (Blue)  | .   | .   | .   | .   |
| Coarse Grit (Green) | .   | .   | .   | .   |

**Cylinders**

|                     |     |     |     |     |     |     |     |     |
|---------------------|-----|-----|-----|-----|-----|-----|-----|-----|
|                     |     |     |     |     |     |     |     |     |
| Catalog Number      | E-1 | E-2 | E-3 | E-4 | E-5 | E-6 | E-7 | E-8 |
| Extra Fine (white)  | .   | .   | .   |     |     |     |     |     |
| Fine Grit (Yellow)  | .   | .   | .   | .   | .   | .   | .   |     |
| Medium Grit (Blue)  | .   | .   | .   | .   | .   | .   | .   | .   |
| Coarse Grit (Green) |     | .   | .   | .   | .   | .   | .   |     |

**Flame**

|                     |     |     |
|---------------------|-----|-----|
|                     |     |     |
| Catalog Number      | F-1 | F-2 |
| Extra Fine (white)  |     |     |
| Fine Grit (Yellow)  | .   | .   |
| Medium Grit (Blue)  | .   | .   |
| Coarse Grit (Green) | .   | .   |

**Flat End Tapers**

|                     |      |      |
|---------------------|------|------|
|                     |      |      |
| Catalog Number      | F-1R | F-2R |
| Extra Fine (white)  |      |      |
| Fine Grit (Yellow)  |      |      |
| Medium Grit (Blue)  | .    | .    |
| Coarse Grit (Green) |      | .    |

**Round End Tapers**

|                     |     |     |     |     |     |     |     |     |     |
|---------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|
|                     |     |     |     |     |     |     |     |     |     |
| Catalog Number      | K-1 | K-2 | L-1 | M-1 | M-2 | M-3 | M-4 | M-5 | M-6 |
| Extra Fine (white)  |     | .   |     |     |     |     |     |     |     |
| Fine Grit (Yellow)  | .   | .   |     |     |     | .   | .   |     |     |
| Medium Grit (Blue)  | .   | .   | .   | .   | .   | .   | .   | .   | .   |
| Coarse Grit (Green) | .   | .   |     |     |     | .   | .   |     |     |

**Gingival Curettage Set**

|                     |     |     |     |     |
|---------------------|-----|-----|-----|-----|
|                     |     |     |     |     |
| Catalog Number      | G-1 | G-2 | G-3 | G-4 |
| Extra Fine (white)  |     |     |     |     |
| Fine Grit (Yellow)  |     |     |     |     |
| Medium Grit (Blue)  | .   | .   | .   | .   |
| Coarse Grit (Green) |     |     |     | .   |

**Wheels**

|                                      |    |    |    |    |    |    |    |    |    |    |    |
|--------------------------------------|----|----|----|----|----|----|----|----|----|----|----|
| Available in Medium Grit ONLY (BLUE) |    |    |    |    |    |    |    |    |    |    |    |
| Catalog Number                       | S1 | S2 | S3 | S4 | S5 | T1 | T2 | T3 | T4 | T5 | T6 |

**Gingival Curettage Set**

|                    |     |     |     |     |     |     |     |     |     |     |
|--------------------|-----|-----|-----|-----|-----|-----|-----|-----|-----|-----|
| Super Coarse (RED) |     |     |     |     |     |     |     |     |     |     |
| Catalog Number     | FP2 | CP2 | EP5 | EP6 | SP5 | KP2 | MP3 | BR2 | BR3 | BR4 |

**For Gross Reduction - Super Coarse**

|                        |        |
|------------------------|--------|
| Extra Fine . . . . .   | \$4.95 |
| Fine . . . . .         | \$4.95 |
| Medium . . . . .       | \$4.95 |
| Coarse . . . . .       | \$5.95 |
| Super Coarse . . . . . | \$6.75 |

**VOLUME DISCOUNTS**  
AVAILABLE - *Please ask*

**Satisfaction Guarantee**



# Sabra Sterile Single Use Diamonds

## Sterile Single Use Diamonds

Made with Natural DeBeers diamonds, with 25% more diamond surface exposed for greater efficiency for single use cutting.

All styles \$1.29 per bur

Packaged individually in sterile packs, 25 per box. - \$32.95 per box of 25

1-5 boxes- \$32.25 per box

6-10 boxes- \$29.95 per box

11 or more - \$ 27.95 per box

### ROUND END TAPER

|                  |           |            |           |           |            |           |           |            |           |
|------------------|-----------|------------|-----------|-----------|------------|-----------|-----------|------------|-----------|
|                  |           |            |           |           |            |           |           |            |           |
| ORDER/ISO #      | 856.012   | 850.012    | S855.014  | 856.014   | 850.014    | S855.016  | 856.016   | 850.016    | 856.018   |
| GRITS            | C M       | C M        | SC        | C M F     | SC M F     | SC M      | SC M F XF | SC M F     | SC F      |
| Head/Lgh/Tip(mm) | 1.2/7/0.7 | 1.2/10/0.7 | 1.4/6/0.9 | 1.4/8/0.9 | 1.4/10/0.9 | 1.6/6/1.1 | 1.6/8/1.1 | 1.6/10/1.1 | 1.8/8/1.4 |
| Neo. #           | 1112.7    | 1112.10    | 1114.6    | 1114.8    | 1114.10    | 1116.6    | 1116.8    | 1116.10    | 1118.7    |

### ROUND END TAPER

### FLAT END TAPER

|                  |           |           |           |           |           |           |           |           |           |           |
|------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|                  |           |           |           |           |           |           |           |           |           |           |
| ORDER/ISO #      | 856N.018  | 856L.018  | 856.021   | 856.022   | S856.024  | 856.024   | 847.012   | 847L.014  | 847.014   | 847.016   |
| GRITS            | C         | SC M F    | SC        | SC        | SC        | SC        | C         | C         | SC M F    | SC M F    |
| Head/Lgh/Tip(mm) | 1.8/8/1.2 | 1.8/9/1.4 | 2.1/8/1.6 | 2.2/8/1.8 | 2.4/7/1.8 | 2.4/7/1.8 | 1.2/8/0.8 | 1.4/9/0.7 | 1.4/8/1.0 | 1.6/8/1.2 |
| Neo. #           | 1118.8    | 1118.9    | 1121.8    | 1122.8    | 2424CS    | 2424      | 0914.9    | 0914.8    | 0914.8    | 0916.8    |

### FLAT END TAPER

### MODIFIED FLAT END TAPER

### MODIFIED CHAMFER

|                  |            |           |            |            |           |           |           |           |           |            |
|------------------|------------|-----------|------------|------------|-----------|-----------|-----------|-----------|-----------|------------|
|                  |            |           |            |            |           |           |           |           |           |            |
| ORDER/ISO #      | 848.016    | 847.018   | 848.018    | S846KR.016 | 847KR.016 | 845KR.018 | 847KR.018 | 845KR.025 | 878K.012  | 879K.012   |
| GRITS            | SC M       | SC M      | SC M       | SC M       | SC M F    | SC M      | SC        | M         | C M       | C M        |
| Head/Lgh/Tip(mm) | 1.6/10/1.1 | 1.8/8/1.4 | 1.8/10/1.4 | 1.6/6/1.1  | 1.6/8/1.1 | 1.8/4/1.4 | 1.8/8/1.4 | 2.5/4/2.0 | 1.2/8/0.8 | 1.2/10/0.8 |
| Neo. #           | 0916.10    | 0918.8    | 0918.10    | 0816.6S    | 0816.8    | 0818.4    | 0818.8    | 2525      | 1712.8    | 1712.10    |

### MODIFIED CHAMFER

### MOD. FLAT END CYL.

### FLAT END CYLINDER

|                  |           |            |           |           |           |            |           |           |           |           |
|------------------|-----------|------------|-----------|-----------|-----------|------------|-----------|-----------|-----------|-----------|
|                  |           |            |           |           |           |            |           |           |           |           |
| ORDER/ISO #      | 878K.014  | 879K.014   | S877K.016 | 878K.016  | 878K.018  | 879K.018   | 878K.021  | 837KR.014 | S835.010  | 835.010   |
| GRITS            | C M F     | C M        | C M       | SC M F    | SC F      | SC M       | C         | C F       | C M       | C M       |
| Head/Lgh/Tip(mm) | 1.4/8/1.0 | 1.4/10/1.1 | 1.6/6/1.2 | 1.6/8/1.2 | 1.8/8/1.4 | 1.8/10/1.4 | 2.1/8/1.6 | 1.4/8/1.4 | 1.0/4/1.0 | 1.0/4/1.0 |
| Neo. #           | 1714.8    | 1714.10    | 1716.6S   | 1716.8    | 1718.8    | 1718.10    | 1721.8    | 0614.8    | 0710S     | 0710      |

### FLAT END CYLINDER

### ROUND END CYLINDER

### MOD. BEVELED CYL.

### BEVELED CYLINDER

|                  |           |           |           |           |           |           |           |           |           |           |
|------------------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|-----------|
|                  |           |           |           |           |           |           |           |           |           |           |
| ORDER/ISO #      | 835.012   | 837.012   | 837.016   | 881-009   | 881.012   | 881.014   | 878.012   | 878.014   | 885.010   | 885.012   |
| GRITS            | C         | C         | C         | C F       | SC        | SC        | C F       | C         | M         | C F       |
| Head/Lgh/Tip(mm) | 1.2/4/1.2 | 1.2/7/1.2 | 1.6/8/1.6 | 0.9/7/0.9 | 1.2/7/1.2 | 1.4/8/1.4 | 1.2/8/1.2 | 1.4/8/1.4 | 1.0/8/1.0 | 1.2/8/1.2 |
| Neo. #           | 0712      | 0712.7    | 0716.8    | 1210.7    | 1212.7    | 1214.8    | 1800.8    | 1804.8    | 1812.8    | 1812.8    |

### BEVELED CYLINDER

### ROUND

### INVERTED CONE

|                  |            |         |         |         |         |         |         |         |         |         |
|------------------|------------|---------|---------|---------|---------|---------|---------|---------|---------|---------|
|                  |            |         |         |         |         |         |         |         |         |         |
| ORDER/ISO #      | 886.012    | 801.010 | 801.012 | 801.016 | 801.018 | 801.023 | 805.012 | 805.016 | 807.016 | 807.018 |
| GRITS            | C M        | C       | M       | C M     | C M     | C F     | M       | C       | C       | C       |
| Head/Lgh/Tip(mm) | 1.2/10/1.2 | 1.0/0.8 | 1.2/1.0 | 1.6/1.4 | 1.8/1.6 | 2.3/2.0 | 1.2/0.9 | 1.6/1.4 | 1.6/4.0 | 1.8/5.0 |
| Neo. #           | 1812.10    | 0110    | 0112    | 0116    | 0118    | 0123    | 0312    | 0316    | 0316.4  | 0318.5  |

### DOUBEL INVERTED CONE

### PEAR

### END CUTTER

|                  |         |          |         |         |          |         |          |           |           |           |
|------------------|---------|----------|---------|---------|----------|---------|----------|-----------|-----------|-----------|
|                  |         |          |         |         |          |         |          |           |           |           |
| ORDER/ISO #      | 813.014 | S830.008 | 830.008 | 830.010 | S830.012 | 830.012 | 830L.012 | 10839.012 | 10839.014 | 10839.016 |
| GRITS            | C       | M        | M       | M       | M        | M       | C        | M         | M         | M         |
| Head/Lgh/Tip(mm) | 1.4/2.2 | 0.8/1.6  | 0.8/1.6 | 1.0/2.7 | 1.2/2.7  | 1.2/2.7 | 1.2/4    | 1.2       | 1.4       | 1.6       |
| Neo. #           | 0416    | 0508S    | 0508    | 0512S   | 0512S    | 0512    | 5012     | 5012      | 5014      | 5016      |



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.

# Sabra Sterile Single Use Diamonds



All styles \$1.29 per bur

Packaged individually in sterile packs, 25 per box. - \$32.95 per box of 25

1-5 boxes- \$32.25 per box


6-10 boxes- \$29.95 per box

11 or more - \$ 27.95 per box

|                  | DOUBEL INVERTED CONE |                 |                 | PEAR            |                 |                 |                    | END CUTTER       |                  |                  |                |
|------------------|----------------------|-----------------|-----------------|-----------------|-----------------|-----------------|--------------------|------------------|------------------|------------------|----------------|
|                  |                      |                 |                 |                 |                 |                 |                    |                  |                  |                  |                |
| ORDER/ISO #      | <b>813.014</b>       | <b>S830.008</b> | <b>830.008</b>  | <b>830.010</b>  | <b>S830.012</b> | <b>830.012</b>  | <b>830L.012</b>    | <b>10839.012</b> | <b>10839.014</b> | <b>10839.016</b> |                |
| GRITS            | <b>C</b>             | <b>M</b>        | <b>M</b>        | <b>M</b>        | <b>M</b>        | <b>M</b>        | <b>C</b>           | <b>M</b>         | <b>M</b>         | <b>M</b>         |                |
| Head/Lgh/Tip(mm) | <b>1.4/2.2</b>       | <b>0.8/1.6</b>  | <b>0.8/1.6</b>  | <b>1.0/2.7</b>  | <b>1.2/2.7</b>  | <b>1.2/2.7</b>  | <b>1.2/4</b>       | <b>1.2</b>       | <b>1.4</b>       | <b>1.6</b>       |                |
| Neo. #           | 0416                 | 0508S           | 0508            |                 | 0512S           | 0512            |                    | 5012             | 5014             | 5016             |                |
|                  | INTER-PROXIMAL       |                 | NEEDLES         |                 |                 |                 | WHEEL              |                  | EGG              |                  |                |
|                  |                      |                 |                 |                 |                 |                 |                    |                  |                  |                  |                |
|                  |                      |                 | Pit & Fissure   |                 |                 |                 |                    |                  |                  |                  |                |
| ORDER/ISO #      | <b>392.016</b>       | <b>S858.008</b> | <b>858.008</b>  | <b>859L.012</b> | <b>858.010</b>  | <b>858.012</b>  | <b>859.014</b>     | <b>909.035</b>   | <b>909.042</b>   | <b>379.023</b>   |                |
| GRITS            | <b>M F</b>           | <b>F</b>        | <b>M F XF</b>   | <b>C F</b>      | <b>M</b>        | <b>C M F XF</b> | <b>C M F XF</b>    | <b>C</b>         | <b>C</b>         | <b>SC M F XF</b> |                |
| Head/Lgh/Tip(mm) | <b>1.6/5</b>         | <b>0.8/3</b>    | <b>0.8/3</b>    | <b>1.2/11</b>   | <b>1.0/8</b>    | <b>1.2/8</b>    | <b>1.4/10</b>      | <b>3.5/1.3</b>   | <b>4.2/1.5</b>   | <b>2.3/4.2</b>   |                |
| Neo. #           | 1416                 | 1300FS          | 1310.3          | 1312.11         |                 | 1312.8          | 1314.10            | 2035             | 2042             | 1900             |                |
|                  | EGG                  |                 | FOOTBALL        |                 | DEPTH CUTTER    |                 |                    | FLAME            |                  |                  |                |
|                  |                      |                 |                 |                 |                 |                 |                    |                  |                  |                  |                |
| ORDER/ISO #      | <b>379.018</b>       | <b>368.016</b>  | <b>368.020</b>  | <b>368.023</b>  | <b>828.026</b>  | <b>834.016</b>  | <b>834.021</b>     | <b>889.009</b>   | <b>862.010</b>   | <b>S862.012</b>  | <b>861.012</b> |
| GRITS            | <b>C F</b>           | <b>C F XF</b>   | <b>C F SC</b>   | <b>M F XF</b>   | <b>M</b>        | <b>M</b>        | <b>M</b>           | <b>M F</b>       | <b>C M F</b>     | <b>C M</b>       | <b>C XF</b>    |
| Head/Lgh/Tip(mm) | <b>1.8/3.4</b>       | <b>1.6/3.2</b>  | <b>2.0/5.0</b>  | <b>2.3/5.2</b>  | dep: 0.5        | dep: 0.3        | dep: 0.5           | <b>0.9/3.0</b>   | <b>1.0/8.0</b>   | <b>1.2/8.0</b>   | <b>1.2/7.0</b> |
| Neo. #           | 1908                 | 1916            | 1920            | 1923            | 6005            | 7003            | 7005               | 1509.3           | 1510.8           | 1512.6           | 1512.7         |
|                  | FLAME                |                 |                 |                 |                 |                 | OCCLUSAL REDUCTION |                  | GROSS REDUCTION  |                  |                |
|                  |                      |                 |                 |                 |                 |                 |                    |                  |                  |                  |                |
| ORDER/ISO #      | <b>862.012</b>       | <b>863.012</b>  | <b>862N.014</b> | <b>862.014</b>  | <b>862.016</b>  | <b>862L.018</b> | <b>811.033</b>     | <b>811.037</b>   | <b>6051.17</b>   | <b>6055.18</b>   |                |
| GRITS            | <b>C M F XF</b>      | <b>C M F XF</b> | <b>C M</b>      | <b>C</b>        | <b>C M F XF</b> | <b>C F</b>      | <b>SC</b>          | <b>SC</b>        | <b>C</b>         | <b>C</b>         |                |
| Head/Lgh/Tip(mm) | <b>1,2/8.0</b>       | <b>1.2/10</b>   | <b>1.4/7.5</b>  | <b>1.4/8</b>    | <b>1.6/8</b>    | <b>1.8/9.0</b>  | <b>3.3/5.0</b>     | <b>3.7/7.0</b>   | <b>1.7/8/1.7</b> | <b>1.8/8/1.3</b> |                |
| Neo. #           | 1512.8               | 1512.10         | 1514.7          | 1514.8          | 1516.8          | 1518.9          | 2133               | 2137             | 2217             | 2218             |                |



# Sabra Cerami-Polishers



- Diamond Impregnated Ceramic Polishers
- Our line has assorted polisher shapes, grits and sizes, allowing you to logically and quickly gross trim, reduce and contour (blue) and pre-polish (pink). Then, use our fine grit (gray) polishers to achieve that final, high gloss, "wet look".
- Optimal speed range: 5,00 - 10,000 rpm  
All \$8.99 Per Polisher



**SPF - Coarse Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Flame Shape, Latch Bur, Reduction and Contouring - (Blue)**

|                         |      |
|-------------------------|------|
| Size                    | 040  |
| Length (mm)             | 10.0 |
| Polish                  | Blue |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPF - Medium Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Flame Shape, Latch Bur, Pre-Polish - (Pink)**

|                         |      |
|-------------------------|------|
| Size                    | 040  |
| Length (mm)             | 10.0 |
| Polish                  | Pink |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPF - Fine Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Flame Shape Latch Bur, High Gloss Polish - (Gray)**

|                         |      |
|-------------------------|------|
| Size                    | 040  |
| Length (mm)             | 10.0 |
| Polish                  | Gray |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPI - Coarse Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Inverted Cone Shape, Latch Bur, Reduction and Contouring - (Blue)**

|                         |      |
|-------------------------|------|
| Size                    | 060  |
| Length (mm)             | 7.5  |
| Polish                  | Blue |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPI - Medium Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Inverted Cone Shape, Latch Bur, Pre-Polish - (Pink)**

|                         |      |
|-------------------------|------|
| Size                    | 060  |
| Length (mm)             | 7.5  |
| Polish                  | Pink |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPI - Fine Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Inverted Cone Shape Latch Bur, High Gloss Polish - (Gray)**

|                         |      |
|-------------------------|------|
| Size                    | 060  |
| Length (mm)             | 7.5  |
| Polish                  | Gray |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPD - Coarse Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Disc, Latch Bur, Reduction and Contouring - (Blue)**

|                         |      |
|-------------------------|------|
| Size                    | 100  |
| Length (mm)             | 2.5  |
| Polish                  | Blue |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPD - Medium Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Disc, Latch Bur, Pre-Polish - (Pink)**

|                         |      |
|-------------------------|------|
| Size                    | 100  |
| Length (mm)             | 2.5  |
| Polish                  | Pink |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |



**SPD - Fine Sabra Intra-Oral Ceramic Polishers w/ Diamond Grit-Disc, Latch Bur, High Gloss Polish - (Gray)**

|                         |      |
|-------------------------|------|
| Size                    | 100  |
| Length (mm)             | 2.5  |
| Polish                  | Gray |
| Right-Angle Latch Shank |      |
| 21 Standard             |      |
| 5 Pack                  |      |





# Diamond Separating Strips



High Quality Diamond Separating Strips



|           |      | 5 pcs | SINGLE SIDED |         |         |         |
|-----------|------|-------|--------------|---------|---------|---------|
| WIDTH     | 2 mm |       | ○ U 342      | ● F 342 | ● C 342 | 342     |
|           | 3 mm |       | ○ U 343      | ● F 343 | ● C 343 | 343     |
|           | 4 mm |       | ○ U 344      | ● F 344 | ● C 344 | 344     |
|           | 6 mm |       | ○ U 346      | ● F 346 | ● C 346 | 346     |
| THICKNESS |      |       | 0,07 mm      | 0,10 mm | 0,14 mm | 0,17 mm |

2mm=\$9.99 per 5 pack, 3mm=\$11.99 per 5 pack, 4mm=\$12.99 per 5 pack, 6mm=\$16.99 per 5 pack

## Diamond Separating Strips



**Fully Coated**


|           |      | 10 pcs | SINGLE SIDED |         |         |
|-----------|------|--------|--------------|---------|---------|
| WIDTH     | 2 mm |        | ● F 352      | ● C 352 | 352     |
|           | 4 mm |        | ● F 354      | ● C 354 | 354     |
|           | 6 mm |        | ● F 356      | ● C 356 | 356     |
| THICKNESS |      |        | 0,10 mm      | 0,14 mm | 0,17 mm |

**With Gap**

|           |      | 10 pcs | SINGLE SIDED |          |         |
|-----------|------|--------|--------------|----------|---------|
| WIDTH     | 2 mm |        | ● F 352L     | ● C 352L | 352L    |
|           | 4 mm |        | ● F 354L     | ● C 354L | 354L    |
|           | 6 mm |        | ● F 356L     | ● C 356L | 356L    |
| THICKNESS |      |        | 0,10 mm      | 0,14 mm  | 0,17 mm |

2mm=\$39.95 per 10 pack, 4mm=\$42.95 per 10 pack, 6mm=\$49.95 per 10 pack 2mm=\$39.95 per 10 pack, 4mm=\$42.95 per 10 pack, 6mm=\$49.95 per 10 pack

## Perforated Diamond Separating Strips



**Fully Coated**

|           |      | 10 pcs | SINGLE SIDED |         | DOUBLE SIDED |
|-----------|------|--------|--------------|---------|--------------|
| WIDTH     | 3 mm |        | ● F 363      | ● C 363 | 363          |
|           | 4 mm |        | ● F 364      | ● C 364 | 364          |
|           | 6 mm |        | ● F 366      | ● C 366 | 366          |
| THICKNESS |      |        | 0,10 mm      | 0,14 mm | 0,17 mm      |
|           |      |        |              |         | 0,20 mm      |

**With Gap**

|           |      | 10 pcs | SINGLE SIDED |         |         |
|-----------|------|--------|--------------|---------|---------|
| WIDTH     | 3 mm |        | ● F 373      | ● C 373 | 373     |
|           | 4 mm |        | ● F 374      | ● C 374 | 374     |
|           | 6 mm |        | ● F 376      | ● C 376 | 376     |
| THICKNESS |      |        | 0,10 mm      | 0,14 mm | 0,17 mm |

Single Sided: 2mm=\$39.95 per 10 pack, 4mm=\$42.95 per 10 pack, 6mm=\$49.95 per 10 pack 3mm=\$46.95 per 10 pack, 4mm=\$49.95 per 10 pack, 6mm=\$54.95 per 10 pack  
 Double Sided: 4mm=\$66.95 per 10 pack, 6mm=\$79.95 per 10 pack

## Abrasive Stainless Steel Strips




|           |      | 10 pcs | SINGLE SIDED |     |         |
|-----------|------|--------|--------------|-----|---------|
| WIDTH     | 4 mm |        |              | 444 | \$14.99 |
|           |      |        |              |     |         |
| THICKNESS |      |        |              |     | 0,05 mm |



Call Toll Free at 1-800-888-4435

# Lindeman Type Bone Cutters Bone Cutting, Contouring & Seratted Cutters

“•” in column indicates availability




**ORDER CODE 161**  
FG161-016 Sabra Bone Cutting Bur, 19mm  
RA161-016 Sabra Bone Cutting Bur, 26mm  
HP161-016 Sabra Bone Cutting Bur, 44.5 MM  
**\$25.89 each**  
**\$23.89 each, 2 or more**



**ORDER CODE 267**  
FG267-016 Sabra Bone Cutting bur, 19mm  
HP267-016 Sabra Bone Cutting Bur, 44.5mm  
**\$25.89 each**  
**\$23.89 each, 2 or more**



**ORDER CODE 269**  
FG269-016 Sabra Bone Cutting Bur, 19mm  
**\$25.89 each**  
**\$23.89 each, 2 or more**




**ORDER CODE 162**  
FG162-016 Sabra Bone Cutting Bur, 19mm  
HP162-016 Sabra Bone Cutting Bur, 44.5mm  
RA162-016 Sabra Bone Cutting Bur, 22mm  
RA162-016SU Sabra Bone Cutting Bur, 26mm  
**\$25.89 each**  
**\$23.89 each, 2 or more**




**ORDER CODE 162A**  
FG162A-016 Sabra Bone Cutting Bur, 19mm  
HP162A-016 Sabra Bone Cutting Bur, 44.5mm  
**\$25.89 each**  
**\$23.89 each, 2 or more**

## All Port Bone Cutter Burs



**ORDER CODE H141**  
Head Size 1/10mm - 031,035,040,050  
**Sold per bur at \$19.99 each. 5 or more- \$18.99 each**  
\* 035S is serrated


## Sabra Implant Burs



**Implant Long Bur**  
**ORDER CODE FG 700XXL - \$6.49**

**Implant Pilot Bur**  
**ORDER CODE RA6SS - \$5.99**  
**ORDER CODE RA8SS - \$5.99**

## Surgical Length Diamonds




**ORDER CODE**  
Medium 801-head size-SU - **\$8.99**  
Coarse 6801-head size-SU - **\$8.99**  
Super 5801-head size-SU - **\$10.89**  
ie 5801-029SU  
859-010SU - **\$9.89**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.


## 30mm Surgical Long Burs



ORDER CODE FGSU-L  
Round Bur 30mm shank- **\$3.29**



## Finishing & Trimming Burs



ORDER CODE FGSU-25M  
Sabra surgical length finishing and trimming bur- **\$6.99**

## Bur Kits



### Sabra PFM Crown Removal Kit

The Sabra PFM removal kit is a cost effective way to remove old crowns. Cut a trough with the Sabra super coarse diamond, just wide enough to allow access for your Super metal cutter. Used in conjunction these 2 burs will be used 5-6 uses each, saving you money, time and your sanity.

Includes: 6 CP-2 super coarse diamonds, 6 MC-L super metal cutting burs, 1 bur block

ORDER CODE SPFM - Sabra Dental Products PFM Crown Removal Kit - **\$79.99**



### Dr El Chaar Crown Lengthening Kit

A kit containing all the burs necessary for this treatment to simplify the procedure so that every practitioner can do it, affordably with time saving efficiency as all the tools are here, in order, in one simple kit.

Includes: 3 each of the following: FG 7004 SU, FG 274 SU, FG 4 SU, FG 958 SU in Sabra Bur Block

ORDER CODE EC-CL Kit - **\$89.99**

### Dr. El Chaar Peri-Implantitis Kit

In the ever increasing incidence of peri-implantitis, the need to have tools to simplify the treatment of peri-implantitis is increasing. This kit is made to simplify the treatment and make it available to any practitioner involved in implant dentistry. The burs will help decontaminate the implant surface completely and efficiently reducing chair time and make it very affordable by having everything needed present at time of the surgery in one simple kit.

Includes: 3 each of the following burs 379.018UF, FG7008, 41.018UF, 801-023XF

ORDER CODE EC-PI Kit - **\$139.99**

## Profiler Implant Bur Tool



### Profiler Multifunctional Implant Bur Tool

The Profiler is a universal, multifunctional, time saving, rotary instrument which can be used in your contra-angle with many implant systems to create a soft or hard tissue emergence profile.

- Crestal Reduction: As a rotary instrument, it cuts uneven crestal bone into a leveled implant osseous platform by moving the drill in a buccal-lingual motion. For gross reduction, the countersink area should be used. A black coating/diamond-like hardened finish maintains the cutting edge.
- Osteomy Locator: Pin-points an osteotomy implant location of choice as the preferred site without slipping
- Osteomy Lateral Redirection: Observe bucco-lingual cortical bone plates and redirect osteotomy as desired, using 2mm of the cutting edge
- Countersink: Following osteotomy drilling, it countersinks 2mm deeper, in a tapered mode, for the abutment/implant emergence profiler. Countersink area is also used for gross reduction of ridge
- Osteocompression: The upper portion functions as a rotary osteotome to compact the bone in a controlled fashion

\* Hold your contra angle with two hands; (two-hand technique) one finger concentrically over the Profiler.

ORDER CODE PL-PRO8 - **\$135.00**



Call Toll Free at 1-800-888-4435

# Trephines & Mucosa Membrane Punches

## Trephines

Manufactured from highest grade stainless steel, SABRA Trephines are available in a variety of diameters, with and without laser depth markings and feature windows for enhanced view and for ejection of bone fragments.



| ORDER CODE | Interior Diameter | Exterior Diameter | Teeth | Body Length | Grading          |
|------------|-------------------|-------------------|-------|-------------|------------------|
| 5015003U0  | 1.7mm             | 2.3mm             | 7     | 22mm        | 7 / 10 / 13 / 16 |
| 5015004U0  | 2.3mm             | 2.8mm             | 7     | 22mm        | 7 / 10 / 13 / 16 |
| 5015005U0  | 2.8mm             | 3.3mm             | 9     | 22mm        | 7 / 10 / 13 / 16 |
| 5015006U0  | 3.3mm             | 3.8mm             | 9     | 22mm        | 7 / 10 / 13 / 16 |
| 5015007U0  | 4.0mm             | 4.5mm             | 11    | 22mm        | 7 / 10 / 13 / 16 |
| 5015008U0  | 4.3mm             | 4.8mm             | 11    | 22mm        | 7 / 10 / 13 / 16 |
| 5015009U0  | 4.8mm             | 5.8mm             | 9     | 22mm        | 7 / 10 / 13 / 16 |
| 5015015U0  | 5.0mm             | 5.5mm             | 11    | 22mm        | 7 / 10 / 13 / 16 |
| 5015010U0  | 6.0mm             | 7.0mm             | 12    | 22mm        | 7 / 10 / 13 / 16 |
| 5015011U0  | 7.0mm             | 8.0mm             | 18    | 22mm        | 7 / 10 / 13 / 16 |
| 5015012U0  | 8.0mm             | 9.0mm             | 18    | 22mm        | 7 / 10 / 13 / 16 |
| 5015013U0  | 9.0mm             | 10.0mm            | 18    | 22mm        | 7 / 10 / 13 / 16 |
| 5015014U0  | 10.0mm            | 11.0mm            | 19    | 22mm        | 7 / 10 / 13 / 16 |

ALL  
Trephines  
\$99.99



| ORDER CODE | Interior Diameter | Shaft       |
|------------|-------------------|-------------|
| 5015017U0  | 2.0mm             | 2.35 x 15mm |
| 5015018U0  | 5.0mm             | 2.35 x 15mm |
| 5015019U0  | 7.5mm             | 2.35 x 15mm |
| 5015020U0  | 10.0mm            | 2.35 x 15mm |

## SABRA Mucosa Membrane Punches

SABRA stainless steel Mucosa Punches feature sharp punch edges and are available in multiple diameters.



### Mucosa Membrane Punch

ORDER CODE 5015022U0  
3mm, interior diameter  
\$99.99

### Mucosa Membrane Punch

ORDER CODE 5015023U0  
4mm, interior diameter  
\$99.99

### Mucosa Membrane Punch

ORDER CODE 5015024U0  
5mm, interior diameter  
\$99.99

### Mucosa Membrane Punch

ORDER CODE 5015025U0  
6mm, interior diameter  
\$99.99




Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Low Temp Implant Burs



Thermal injury when drilling bone is well documented and known, but was not addressed until now! Introducing the Low Temp Osseo Drill™: Its unique design enables a less penetrating force, which means less time for heat generation inside the drilled bone. More importantly, it means no damage to the bone as may be evident with other standard bone drills that heat above 40° C. Studies have shown that the Low Temp's unique drill design resulted in a temperature "safe zone" for bone, under 40° C (see graph). Less Heat & Less Force guarantees remarkably enhanced implant integration superior to any other bone drills. Use the Low Temp Osseo Drill™ to decrease the drilling temperature and to increase your rate of success. **\$79.95 Per Drill. Buy 3 or More \$69.95**



### D = 2.0 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 200-85     |
| 10.0mm | 200-100    |
| 11.5mm | 200-115    |
| 13.0mm | 200-130    |
| 15.0mm | 200-150    |

### D = 2.3 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 230-85     |
| 10.0mm | 230-100    |
| 11.5mm | 230-115    |
| 13.0mm | 230-130    |
| 15.0mm | 230-150    |

### D = 2.75 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 275-85     |
| 10.0mm | 275-100    |
| 11.5mm | 275-115    |
| 13.0mm | 275-130    |
| 15.0mm | 275-150    |

### D = 3.0 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 300-85     |
| 10.0mm | 300-100    |
| 11.5mm | 300-115    |
| 13.0mm | 300-130    |
| 15.0mm | 300-150    |

### D = 3.15 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 315-85     |
| 10.0mm | 315-100    |
| 11.5mm | 315-115    |
| 13.0mm | 315-130    |
| 15.0mm | 315-150    |

### D = 3.25 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 325-85     |
| 10.0mm | 325-100    |
| 11.5mm | 325-115    |
| 13.0mm | 325-130    |
| 15.0mm | 325-150    |

### D = 4.25 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 425-85     |
| 10.0mm | 425-100    |
| 11.5mm | 425-115    |
| 13.0mm | 425-130    |
| 15.0mm | 425-150    |


### D = 5.25 mm

| L      | ORDER CODE |
|--------|------------|
| 8.5mm  | 525-85     |
| 10.0mm | 525-100    |
| 11.5mm | 525-115    |
| 13.0mm | 525-130    |
| 15.0mm | 525-150    |

## Osseo™ Tip Bone Grafter

Osseo-Tip™ collects all the bone chips into a removable cup that would otherwise have been wasted during the drilling procedure. After drilling you simply empty the cup of autogenous bone and it is ready for immediate use.

The Osseo-Tip™ Bone Grafter proved to be an effective way to harvest sufficient quantities of autogenous bone without the need to employ a bone mill or morselizer. It allowed for bone to be harvested in a safe and rapid manner with no hand fatigue and without creating a secondary surgical site. Osseo-tip™ is designed to fit your standard handpiece and is fully Autoclavable. **\$129.95 Per Drill. Buy 3 or More \$116.10**



# Implant Drivers and Tools

## Torque Control

Tightening with the correct torque minimizes the risk of the screw working loose and improves the precision during bridge fabrication. Torque control can also reduce the risk for stress being built into the restoration and prevents component fractures.

## Easy to use

Insert the screw or abutment into the implant with the appropriate driver and thumb knob manually and tighten with finger pressure. Final seat with the Sabra adjustable torque wrench with the torque setting recommended by the implant system. Pre-load the screw or abutment by using the torque wrench in the same way as a ratchet wrench. The wrench will slip when the selected torque is applied.

## Customized

The Sabra adjustable torque wrench may be purchased individually or as a customized kit. The kit can be organized to be used for several implant systems.

### SABRA Adjustable Torque Wrench

- 9 Torque Settings: 15, 20, 25, 30, 35, 40, 45, 50, and 60 Ncm
- Autoclavable
- No Lubrication
- Ratchet Mechanism
- CE Mark
- Torque settings lock into select Location for accuracy and repeatability
- Ratchet clutch disengages (slips) when torque is applied to avoid over torquing
- Torque settings are provided for use with all implant designs
- Wrench head accepts the standard 4mm square connection

ORDER CODE SIA-01 - \$489.00



### Sabra Adjustable Torque Wrench

9 settings w/ 9 Tools/ Drivers

**Including:** Thumb knob, ITI/Strauman Torx Long Driver, ITI/Strauman Torx Short Driver, Torx Driver Short, Torx Driver long, .050 Hex Tool Standard, .050 Hex Tool Tapered Long, .048 Hex Tool Tapered Short, .048 Hex Tool Tapered Long.

ORDER CODE SIA-01A - \$769.89



## 4mm Square Driver Handle


One of Sabra's newest products is the 4mm Square Screw Driver Handle. This stainless steel driver handle is approximately 4.5 inches in length and accommodates the full selection of 4mm square drivers. This instrument assists in dental procedures involving the placement of implants and other related components when working in the maxilla or sinus. Four grooves in the driver handle ensure positive-grip action when either tightening or loosening various implant components. This single unit instrument is fully autoclavable.

ORDER CODE SIA-02 - \$129.00



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.

## Sabra 1:1 Multi Driver



### Sabra 1:1 Multi Driver

Multi Use hand wrench for better access

#### Can be used to:

Install implants, abutment screws, root canals, drilling, reaming, countersinking, hole enlargement, and many more typical dental office procedures  
Shape and size allow optimal intra oral access  
Magenetic drive handle securely connects via hex driver to shaft on Doriet style control head and can be removed quickly for cleaning and sterilization  
Completely autoclavable, chemclavable, and dry heat sterilizable.

**ORDER CODE SIA-05** Sabra Manual Multi Wrench - **\$149.89**

**ORDER CODE SIA-05A** Sabra Multi Wrench with light w/50 Plastic Barriers - **\$169.89**

## Anthogyr Driver



### TORQ CONTROL Adjustable Torque Wrench 10-35NCM

- Recommended with any type of implants
- High Precision
- Perfect control of torque
- 7 torques values: 10/1/20/25/30/32/35NCM
- Optimum Security
- Automatic declutching
- Prevents the screw from loosening and fracture
- Ergonomic design
- Lightweight for a better freedom of movement: only 135g!
- Optimized access in mouth thanks to the micro-head (posterior areas)
- Reliability
- High precision design and manufacturing
- Easy maintenance
- It can withstand sterilization (in an autoclave at 135° C)
- One piece design with smooth surface to limit bacterial retention

**ORDER CODE 15501** - **\$850.00**

### SABRA-TORQ Torque Wrench

- Affordable
- Easy to Use
- Achieves consistent pre-loading of implant abutments and screws
- Easy Maintenance-Stainless Steel-Sterilize by Autoclave or chemclave

Tightening with the correct torque minimizes the risk of the screw working loose and improves the precision during bridge fabrication. Torque control can also reduce the risk of stress being built into the restoration and prevent component fractures. The Sabra-Torq Wrench is a pre-calibrated torque wrench system that prevents over and under tightening of implant screws and abutments. It is available in 10Ncm, 20Ncm, 30Ncm, 32Ncm and 35Ncm settings to meet the needs of a variety of implant systems. An assortment of drivers fit all the hex or slot head screws regardless of implant type. Easy to use Insert the screw or abutment into the implant with the driver and tighten with finger pressure. Remove the driver and place into the Sabra-Torq Wrench. Pre-load the screw by using the wrench in the same way as a ratchet wrench. The wrench clutch will disengage when the pre-calibrated torque is applied.

**ORDER CODE TW10** 10NCM Sabra-Torq Torque Wrench - **\$159.00**

**ORDER CODE TW20** 20NCM Sabra-Torq Torque Wrench - **\$159.00**

**ORDER CODE TW30** 30NCM Sabra-Torq Torque Wrench - **\$159.00**

**ORDER CODE TW35** 35NCM Sabra-Torq Torque Wrench - **\$159.00**

### Direction Indicator

**ORDER CODE DI1620** Direction Indicator, 1.6/2.0mm ( Paralell Pins ) - **\$14.99**

**ORDER CODE DI2025** Direction Indicator, 2.0/2.5mm ( Paralell Pins ) - **\$14.99**

### Reversable Ratchet Wrench

**ORDER CODE RWA** Reversable Ratchet Wrench - **\$109.00**






# Implant Drivers, Extensions & Adaptors

## SHAFT LENGTHS:

SHORT: 6.4mm  
LONG: 12.7mm



| TIPS | HEX DRIVERS | SPECS   |
|------|-------------|---|
|      |             | .035 Hex Tapered<br>Short: HCT<br>Long: HCLT<br>\$34.95       |
|      |             | .035 Hex Short<br>Short: HC<br>Long: HCL<br>\$34.95           |
|      |             | .040 Hex Short<br>Short: HSM<br>Long: HSML<br>\$34.95         |
|      |             | .048 Hex Tapered<br>Short: HNT<br>Long: HNLT<br>\$34.95       |
|      |             | .048 Hex Short<br>Short: HN<br>Long: HNL<br>\$34.95           |
|      |             | .050 Hex Tapered<br>Short: HTT<br>Long: HTLT<br>\$34.95       |
|      |             | .052 Hex Tapered<br>Short: HT2SMT<br>Long: HT2SMLT<br>\$34.95 |
|      |             | .062 Hex Short<br>Short: HS<br>Long: HSL<br>\$34.95           |
|      |             | .070 Hex Long<br>Long: HTSI<br>\$34.95                        |
|      |             | .080 Hex Short<br>Short: HU<br>Long: HUL<br>\$34.95           |

| TIPS | SPECIALTY DRIVERS | SPECS   |
|------|-------------------|---|
|      |                   | Torx Short: PROC-S<br>Long: PROC-L<br>NobelBiocare<br>\$45.99 |
|      |                   | Torx Short: ITIS<br>Long: ITIL<br>ITI/Straumann<br>\$45.99    |
|      |                   | .050 Square Short: CER<br>Long: CERL<br>\$34.95               |
|      |                   | .017 Slot Short: SN<br>Long: SNL<br>\$34.95                   |
|      |                   | Female Driver .080 Hex<br>Short: HFD<br>\$34.95               |
|      |                   | TG Post for 3i<br>ITI/Straumann<br>TGD<br>\$62.00             |

| TIPS | LATCH STYLE DRIVERS | SPECS  |
|------|---------------------|--|
|      |                     | .050 latch Style, Square Driver<br>T17/24<br>\$39.95     |
|      |                     | .040 Latch Style, Tapered Hex Driver<br>H40RA<br>\$39.95 |
|      |                     | .048 Latch Style, Tapered Hex Driver<br>HTRA<br>\$39.95  |
|      |                     | .050 Latch Style, Tapered Hex Driver<br>HTRA<br>\$39.95  |

| TIPS | EXTENSIONS | SPECS                     |
|------|------------|---------------------------|
|      |            | 3/8" Extension<br>\$34.95 |
|      |            | 5/8" Extension<br>\$34.95 |

| TIPS | ADAPTORS | SPECS   |
|------|----------|---|
|      |          | NA NobelBiocare Torque Wrench Adaptor<br>\$62.99    |
|      |          | ITIA ITI/Straumann Torque Wrench Adaptor<br>\$62.99 |
|      |          | LS Latch to Square Adaptor<br>\$39.95               |

| MANUAL KNOB | SPECS                 |
|-------------|-----------------------|
|             | Thumb Knob<br>\$39.99 |

To use Driver tools in latch type handpiece attachment



**ORDER CODE IMDSN**  
Latch Type Driver Mount, Short, 1.39CM Shank - \$59.99  
**ORDER CODE IMDLN**  
Latch Type Driver Mount, Long, 2.29CM Shank - \$59.99

Product photos not to scale



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Universal Drivers



Long and short handled drivers work with all major implant systems. Finally, universal drivers for abutments and cover screws that will work with 95% of all implant systems sold in the United States. No more wasted time or rescheduling appointments while you try and figure out what type of implant system was placed in a patient. With these drivers you have the instrument you need. Our long handled drivers are designed for anterior mandible cases. They are long enough to access the chimney in screw-retained prosthesis.

All Drivers are color coded for easy identification. They are autoclavable and either sold individually or in sets of 7. Short handled drivers are 1.75 inches long. Long handled drivers are 4.25 inches long, and are great for lab use as well.



**ORDER CODE OTA540** Short Handled Driver, Set of 7 - \$300.00  
Individual Drivers, Short - \$75.00

|                         | Tips                                 | Short  | Long   |
|-------------------------|--------------------------------------|--------|--------|
| Straight Screwdriver    | Straight Systems                     | OTA535 | OTA530 |
| SCS Screwdriver         | Straumann                            | OTA536 | OTA533 |
| UNIgrip Driver          | Replace Select                       | OTA534 | OTA527 |
| Hex 1.3mm/.050 inch     | Lifecore, Biohorizons, Astra, Zimmer | OTA524 | OTA531 |
| Hex 1.2mm/.048 inch     | 3i System                            | OTA537 | OTA532 |
| Hex 0.9mm/.035 inch     | Friadent, Ankylos, 3i cover Screw    | OTA529 | OTA528 |
| Square 1.3mm/.0505 inch | Gold Screw/Nobel Biocare             | OTA525 | OTA526 |



**ORDER CODE OTA539** Long Handled Driver, Set of 7 - \$300.00  
Individual Drivers, Long - \$75.00



## RA Driver Bit Kit

### RA Driver Bit Kit

Compatible with all major implant systems  
 For use with any right angled driver or standard latch type contra angle system  
 Stop wasting time rescheduling patients because you have to identify what system was placed by another dentist. This system covers 95% of the systems sold in the US. Laser etched for easy ID.

- BIT1 - Straight Screwdriver
- BIT2 - SCS Screwdriver, Strauman
- BIT3 - UniGrip Driver, Replace Select
- BIT4 - Hex 1.3mm/.050", Lifecore, Bio-Horizons, Astra, Zimmer
- BIT5 - Hex 1.2mm/.048", 3i System
- BIT6 - Hex 0.9mm/.035", Friadent, Ankylos, 3i Cover Screw
- BIT7 - Square 1.3mm/.050", Gold Screw/ Nobel Biocare

**ORDER CODE OTABITS Set of 7 - \$300.00**  
**Individual bits - \$75.00 per**



## Osseous Glider "Autogenous Bone Collection"

The disposable bone scraper that always has a fresh blade. Tired of harvesting bone shavings with a scraper that won't stay sharp? Then it's time you tried the Osseous Glider. The 360 degree blade is marked into four sections. As the blade starts to dull, you can expose a new cutting surface by simply twisting the blade with an .050 hex driver. There's no easier way to precisely harvest bone in an atraumatic manner. Because of the restrictions in the posterior of the mouth, the straight Osseous Glider is recommended for the ramus, symphysis and tuberosity areas. The angled Osseous Glider is optimized for the external oblique ridge. If you've never used an angled glider in this area of the mouth, you won't believe how much easier it is to use than a straight Glider. The secret is to hold the angled Glider at a 35 degree angle on the ridge.

**ORDER CODE OTA200** Osseous Glider, Straight - \$50.00  
**ORDER CODE OTA205** Osseous Glider, Angled - \$50.00



## Impla Parallel Implantology Drill Guide Set

Re-establishing harmony between function and aesthetics should be the main objective of each implant supported restoration. Precise pre-operative planning and surgical implant insertion achieve this ideal. Exact positioning of multiple implant axes result in accurate fit, superior function and excellent aesthetics.

- Can be used for any implant system
- Parallel, reciprocal implant alignment
- Wide choice of horizontal implant positions according to surgical aspects
- Operation time is reduced
- Prosthetic result is optimized
- Autoclavable

**ORDER CODE 536100 - \$129.95**  
 Cover for Impla Paralell  
**ORDER CODE 536101 - \$19.60**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
 ©2013 Sabra Dental Products, All Rights Reserved.

# Bone Mill Instruments & Kit



## SABRA Bone Mill Instruments

SABRA Bone Mill reduces harvested bone to a granular consistency that is ideal for adaptation to the grafting site. The TAVA™ Bone Mill is able to be easily disassembled for cleaning and sterilization, and the kit includes a sinus lift plugger/spoon instrument as well as a cleaning/collection device for use with the unit.



**Bone Mill\***  
48mm, (2 pcs.)\*  
ORDER CODE 5013582U0 - \$1363.50



**Sinus-Lift Plugger and Spoon**  
double-ended\*  
ORDER CODE 5013554U0 - \$160.60



**Cleaning Instrument to Collect Autologous Bones\***  
ORDER CODE 5015027U0 - \$151.50

## SABRA Bone Mill Forceps

SABRA Bone Mill Forceps are ideal for reducing smaller pieces of harvested bone to a granular level. The unique design of these forceps, with parallel orientation of the chamber to the milling surface, ensures that the entire unit of bone is milled. The Bone Mill Forceps is easily disassembled for cleaning and sterilizing.



**Bone Mill Forceps**  
ORDER CODE 5015028U0 - \$484.80



**Bone Mill Kit**  
ORDER CODE 5015026U0 - \$2222.00

\*Included in the complete kit

## SABRA Augmentation Applicators

SABRA Augmentation Applicator, constructed of stainless steel and available in various tip diameters, is ideal for placing autogenous bone or bone replacement materials in the defect site.



**Augmentation Applicator**  
ORDER CODE 5015029U0 2.5 mm diameter - \$131.30  
ORDER CODE 5015030U0 3.5 mm diameter - \$131.30  
ORDER CODE 5015031U0 4.5 mm diameter - \$131.30  
ORDER CODE 5015032U0 6.0 mm diameter - \$131.30





# Sapian Root Remover System

## Sapian Root Remover System - Atraumatic root extraction

### Advantages of the Sapian root remover system:

- Fast simple effective
- The most versatile atraumatic extraction device
- Minimally invasive, optimizes extraction site for immediate implant placement
- Less traumatic than traditional extraction methods, resulting in faster healing times and greatly reduced post op pain
- Easy to use excellent printed color and video instructions with kit
- Currently used in dental schools nationwide
- Tactile feel connects you to the procedure
- Reusable extraction screws can be straightened if necessary with included straightening tool. Screws can endure many procedures if cared for
- Preferred atraumatic extraction device for many dentists and the world

### The Sapian root remover kit includes:

- Anterior driver
- Collet Nut for anterior driver
- Manual latch head driver (universal use, and for installing extracts screws)
- Pry bar
- Gear oil vial and spout for manual latch head driver
- Extraction screw block with built in screw straightener
- Aluminum screw block containing:
  - Three thin extraction screws, 1 ball, 2 ball and 3 ball lengths
  - Three thick extracts screws, 1 ball, 2 ball and 3 ball lengths
  - Finger driver for removing roots from extractions screws
  - FG557SU high speed cutting bur used as a pilot drill
  - Short tapered pilot drill
  - Long tapered pilot drill
- Fulcrum tray with removable silicone elastomer cushion pad
- Soft silicone pin mat cushion for screw lock and fulcrum tray
- Autoclavable stainless steel instrument cassette

ORDER CODE RR0001 Sapian Root Remover System - \$1895.00

**Sapian Root Remover System Components**

1. Anterior driver
2. Collet Nut for Anterior Driver
3. Manual Latch Head Driver
4. Pry Bar
5. Gear oil vial and spout for Manual Latch Head Driver
6. Extraction Screw block with built-in screw straightener
7. Aluminum Screw Block containing:
  - a. Three Thin Extraction Screws: one 1-ball, one 2-ball and one 3-ball.
  - b. Three Thick Extraction Screws: one 1-ball, one 2-ball and one 3-ball.
  - c. Finger Driver
  - d. # 557 Surgical high speed cutting bur
  - e. 35.5 mm total length RA tapered pilot drill
  - f. 41mm total length RA tapered pilot drill
8. Fulcrum Tray with removable Silicone Elastomer Cushion pad
9. Soft Silicone Pin Mat Cushion for Screw Block and Fulcrum Tray
10. Autoclavable Stainless Steel Sterilization Cassette

For the atraumatic/minimally invasive price of: **\$1895.00**

# Bone Expander Kit

**Includes:** 5 expanders 2.6mm, 3.0mm, 3.4mm, 3.8mm, 4.3mm with markings at 8mm, 10mm, 11.5mm, and 13mm, ratchet wrench, pilot drill, ratchet wrench extender, manual knob, and latch style handpiece adaptor, in an autoclavable case.

**An Alternative to Osteotomes:** Bone Expander drills are an alternative to osteotomes for the expansion and condensing of the atrophic mandible and maxilla in preparation for dental implant insertions. Expanders are also an alternative to the maxillary sinus elevation technique.

**Surgical Trauma Decreased:** Expanders are driven into the bone with a ratchet wrench or low speed handpiece. This decreases the surgical trauma of osteotomes. Bone expanders improve clinical outcome by improving stability, maintaining bone density and increasing fixation.

ORDER CODE SBE-01 - \$495.00

**Bone Expander Kit**

5 bone expander diameters

4 score lines per expander

13mm  
11.5mm  
10mm  
8mm

2.6 3.0 3.4 3.8 4.3mm

Ratchet Wrench

Pilot Drill

Ratchet Wrench Extender

Manual Knob

Latch Style Hand Piece Adaptor

(Kit Components Above- Not to scale)



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Mini Plate Kit Titanium Mesh, Tack and Screw Kit



## TriStar™ (MiniPlate™) Bone Graft Fixation and Titanium Mesh Full Kit

- \* **Kit includes:** graft, screws and accessories needed for onlay bone grafting procedures
- \* **Full Kit includes:** compact 5 in 1 sterilization cassette and 50 screws / tenting screws
- \* Kit also includes your choice of two different titanium mesh sheets

The MiniPlate™ Bone Graft Fixation and Titanium Mesh Full Kit includes the following components:

Sterilization Cassette (box only), Replacement Screwdriver handle with spring-loaded design, Replacement Quick release blade, Twist Drill for use with screwdriver 1.3 x 5.0mm length

50 total screw: *Self Drilling Tack Screws* : 5 pack of 1.5 x 3mm, 5 pack of 1.5 x 4mm, 5 pack of 1.5 x 5mm. *Self Tapping Screws*: 5 pack of 2.0x 6mm, 5 pack of 2.0 x 8mm, 5 pack of 2.0 x 10mm, 5 pack of 2.0 x 12mm, 5 pack of 2.0 x 14mm. *Tenting Screw (Self Drilling)*: 5 pack of 2.0x 4mm, 5 pack of 2.0 x 6mm.

*Your Choice of titanium mesh:* 40 x 60 x .02mm Titanium Mesh Sheet (non-sterile) or 18 x 25 x .02mm Titanium Mesh Sheet (non-sterile - 3 mesh per box)

**ORDER CODE MP-FKIT - \$2399.00**



## TriStar™ (MiniPlate™) Bone Graft Fixation and Titanium Mesh Starter Kit

### Kit Includes:

All Handle and driver components as Full Kit

Your Choice of 15 screws:

5/pk: Self Drilling Tack Screw 1.5 x 3mm, 4mm, 5mm

5/pk: Self Tapping Screw 2.0 x 6mm, 8mm, 10mm, 12mm, 14mm

5/pk: Tenting Screw (Self Drilling) 2.0 x 4mm, 6mm

Your Choice of titanium mesh: 40 x 60 x .02mm Titanium Mesh Sheet (non-sterile), 18 x 25 x .02mm Titanium Mesh Sheet (non-sterile - 3 mesh per box)

Call and we'll assemble the Kit that's right for you. 1-800-888-4435

**ORDER CODE: MP-SKIT - \$1099.00**

**ORDER CODE: MP-BOX** Fixation and Titanium mesh kit sterilization kit ( empty, cassette and lid only ) - **\$250.00**

**ORDER CODE MESH1825** Titanium Mesh Sheet- non sterile 18 x 25 x .02mm 3 pack - **\$144.00**

**ORDER CODE MESH4060** Titanium Mesh Sheet- non sterile 40 x 60 x.02mm - **\$144.00**

**ORDER CODE MP-BLD** Quick Release Blade for Fixation screw/tack kit ( Blade Only ) - **\$59.00**

**ORDER CODE MP-HNDL** ScrewDriver Handle with Press-fit tip design - **\$186.00**

**ORDER CODE MP-BIT** Twist Drill for 2.0mm self tapping screws ,1.3 x 5.0mm length (Driver handle required) - **\$59.00**



**ORDER CODE MPSTS-20.6** Self Drilling Screw 2.0 x 6mm, 5 pack - **\$150.00**

**ORDER CODE MPSTS-20.8** Self Drilling Screw 2.0 x 8mm, 5 pack - **\$150.00**

**ORDER CODE MPSTS-20.10** Self Drilling Screw 2.0 x 10mm, 5 pack - **\$150.00**

**ORDER CODE MPSTS-20.12** Self Drilling Screw 2.0 x 12mm, 5 pack - **\$150.00**

**ORDER CODE MPSTS-20.14** Self Drilling Screw 2.0 x 14mm, 5 pack - **\$150.00**



**ORDER CODE MPTAC-15.3** Self Drilling Tack 1.5 x 3mm, 5 pack - **\$125.00**

**ORDER CODE MPTAC-15.4** Self Drilling Tack 1.5 x 4mm, 5 pack - **\$125.00**

**ORDER CODE MPTAC-15.5** Self Drilling Tack 1.5 x 5mm, 5 pack - **\$125.00**



**ORDER CODE MPTNS-15.4** Tenting Screw, 2.0 x 4mm, self drilling 5 pack - **\$170.00**

**ORDER CODE MPTNS-15.6** Tenting Screw, 2.0 x 6mm, self drilling 5 pack - **\$170.00**



# Sinus Lift Tools

## Sinus Lift Balloons

Looking for the secret to perfect sinus lifts? Try the Sinus Lift Balloon. These unique balloons make it easy to perform an atraumatic membrane elevation. All you have to do is remember this simple formula while inflating the balloon with the syringe: 1cc of saline equals 6mm of membrane elevation. It couldn't be any easier! This innovative product gives predictable results, decreases surgical time and minimizes the possibility of membrane tears. A typical single-tooth or multi-tooth procedure that generally takes 30 to 45 minutes will probably take 10 to 15 minutes with the Sinus Lift Balloon. Graft material required is in direct proportion to the amount of fluid used to inflate the balloon: 1cc of saline will require 1cc of graft material. ( 1cc of saline will elevate the schneiderian membrane by 6mm. )

Three different configurations are available. The straight model features a 3.1mm diameter shaft and has a 4cc capacity. It is typically used for the osteotome socket or closed technique through an extraction site. The micro-mini model has a 1.9mm shaft and has a 1.5cc capacity. It is typically used for entry through the implant osteotomy site when placing a very small implant. The angled model also has a 3.1mm shaft and 4cc capacity. It is designed for an antral approach using the Caldwell-Luc window or open technique.

- ORDER CODE SME14003** Balloon for Sinus Lifts, Angled, 1 Per - **\$105.99**
- ORDER CODE SME14004** Balloon for Sinus Lifts, Micro-Mini, 1 Per - **\$105.99**
- ORDER CODE SME14002** Balloon for Sinus Lifts, Straight 1 Per - **\$105.99**

Please log on and see video of procedure at:  
<http://www.osseoustech.com/media/video-balloon.html>



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.

# AmnioMTM



## AmnioMTM By AmnioGenix

### Injectable multipotential tissue matrix

Engineered by nature to potentiate wound healing and bone regeneration.

AmnioMTM is a novel multipotential extra cellular matrix (ECM) which acts as a natural scaffold. While providing structural support to cells AmnioMTM assists in regulating cell adhesion, cell division, and cell migration.

#### Potential applications:

- Wound covering
- Ridge augmentation
- Soft tissue repair
- Periodontal defects
- Boney defects
- TMJ

#### Bio-active design:

AmnioMTM is an injectable, liquid human amnion allograft that provides a micro scaffold for tissue repair and remodeling. AmnioMTM retains the natural composition of extra cellular matrix molecules including collagen, structural proteins, growth factors, glycoproteins and glycosaminoglycans.

AmnioMTM has anti-inflammatory, anti-microbial, anti-bacterial, and anti-adhesive as well as immunoprivileged properties.

Presence of multipotent cells may have bi-potential osteogenic and adipogenic differentiation capabilities.

Amniotic tissue matrix may suppress potent pro-inflammatory cytokines, IL-1a & IL-1b, leading to accelerated healing response times.

Multipotent cells identified to be present within amniotic membrane tissue: Amniotic epithelial cells, Amniotic mesenchymal cells.

#### Product characteristics:

FDA registered product regulated under 21 CFR Part 1271 and Section 361 of the Public Health Service act.

Cryopreserved liquid allograft, derived from live birth tissue and cells, which is not dependent on the metabolic activity of living cells. Extensive donor selection and screening, according to AATB guidelines, is performed to ensure clinical safety.

Growth factors identified to be present within amniotic membrane tissue: EGF, TGF-B, FGF, PDGF A&B, VEGF, IGF 1&2



| Component   | Feature  | Clinical Benefit   |
|---|--|--|
| Amniotic Morselized Tissue  | Has been shown to be immune privileged (Class II antigens rarely expressed) anti-microbial, anti-bacterial.  | Donor source from live healthy patients resulting in less complications and donor rejection versus traditional allografts, xenografts. And compared to autografts, there is no need for a second surgical procedure to harvest bone resulting in less site morbidity, pain and associated costs. |
| Three Dimensional Extracellular Matrix Scaffold   | Provides mechanical properties such as collagen I-VI and elastin creating osteoconductive environment for cell attachment and hosting.   | Provides structural matrix for use in filling soft tissue defects and bone voids.  |
| Fibronectin, Laminins, Vitronectin and Cytokines.   | Proteins from the ECM that signal cells to attach and migrate. Serve as adhesion ligands and allows for intercellular communication.   | Rarely evokes an immune response with patients experiencing less pain, swelling and inflammation.  |
| Epidermal Growth Factor (EGF), Transforming Growth Factor (TGF-β), Fibroblast Growth Factor (FGF), Platelet Derived Growth Factor A&B ( PDGF A&B) | EGF proliferates epithelial cells. TGF-β acts as an anti-inflammatory, reduces scar tissue, speeding up healing. FGF proliferates cells for collagen matrix. PDGF A&B proliferates connective tissue and enhances soft tissue. | Has been shown to accelerate wound healing rates, bone formation and soft tissue attachment. Also healthcare economics are positive when comparing AmnioMTM to other osteoinductive and osteogenic materials or synthetic products.  |
| Hyaluronic Acid   | It's been reported to increase osteoblastic bone formation in vitro through increased mesenchymal cell differentiation and cell migration.   | Has the potential to reduce scar formation and improve vascularity, resulting in better tissue formation.  |

ORDER CODE AG-010025 - 0.25 ml. - \$545.00

ORDER CODE AG-010050 - 0.50 ml. - \$1090.00

ORDER CODE AG-010125 - 1.25 ml. - \$2725.00

ORDER CODE Cryofreezer - \$3199.00

\*If Cryo freezer is purchased, AmnioMTM prices decrease significantly - please ask



Call Toll Free at 1-800-888-4435

# Grafting / Allograft Materials

We at Sabra are proud to have secured a relationship with Bone Bank Allografts, and their partner tissue processing facility Texas Human Biologics. Together they are a leading supplier of high quality sterile bone and soft tissue allograft to medical professionals for the past 20 years and over 600,000 successful transplantations.



With accreditation by the American Association of Tissue Banks (AATB), registration with the U.S. Food and Drug Administration (FDA) and in compliance with the FDA's Good Tissue Practices, Bone Bank Allografts continues to pursue their mission: To enhance the quality of patient care by providing safe, high quality allograft tissues to health care professionals for transplant and research. They have done so by exceeding the expectations set by the industry and providing exceptional sterile allograft products and providing top level service to the medical community

Sabra now supplies the full line of SteriGraft™ allografts including the SteriShield™ amnion patch, SteriFuse™ DBM, and the complete line of SteriGraft™ dental allograft bone products.

Bone and soft tissue allografts in the SteriGraft™ line are validated to a sterility assurance level (SAL) of 10<sup>-6</sup> in accordance with ISO 11137-2 Method 1.

#### Donor Screening & Processing:

A meticulous donor screening system has been developed to assure that only uniform, high quality tissues will be recovered and processed.

**Donor Screening and Review:** An extensive donor screening process is conducted to review medical / social history. Only those donors that meet our standards and screening requirements are qualified to be accepted for donation. This rigorous screening process is performed in an effort to eliminate potential donors with risk factors which may contribute to the transmission of infectious disease.

**Donor Recovery and Processing:** Our donor recovery partners and their staff are held to the high standards implemented by the American Association of Tissue Banks (AATB). Specific policies and procedures are followed to ensure that donor tissue is recovered appropriately and aseptically.

#### Laboratory Testing:

In accordance with FDA regulations and AATB Standards, donor tissues are tested by laboratories certified under the Clinical Laboratory Improvement Amendments of 1988 (CLIA) and carefully analyzed using the latest technology in order to confirm that the donor is free from the following infectious diseases:

Human Immunodeficiency Virus Type I Ab  
HIV I Antibody  
HIV I Nucleic Acid Test  
Human Immunodeficiency Virus Type II Ab  
HIV II Antibody  
Human T-Lymphotropic Virus Type I  
HTLV I Antibody

Human T-Lymphotropic Virus Type II  
HTLV II Antibody  
Hepatitis B Virus  
HBc Antibody (total)  
HBs Antigen  
Hepatitis C Virus  
HCV Antibody  
HCV Nucleic Acid Test  
Syphilis  
Rapid Plasma Reagin Test  
FTA Confirmation

All donors are evaluated for hemodilution which could potentially affect test results. Only a licensed physician determines final donor suitability.



#### Tissue Processing:

The GraftShield™ process utilizes Improved and efficient techniques during tissue processing. Our processing partner, Texas Human Biologics, operates from a state of the art processing facility which makes the following steps possible.

**Tissue Cleansing:** To prepare for processing, tissues are sanitized using GraftCleanse™. This proprietary blend of cleansing compounds significantly reduces bio-burden and allows for an aesthetic, natural appearance.

**Sterilization:** GraftShield™ utilizes a low dose radiation sterilization method that is validated in accordance with ISO 11137-2 Method 1. Terminal sterilization is achieved with allografts in their final packaging.

**Custom Processing Machines:** We have designed and built custom equipment that enables individual equipment sterilization between cases, thus eliminating the risk of cross contamination - a significant advance in the safety of processing human tissue allografts.

**Allograft Performance:** Each allograft is evaluated at the time of processing, packaging, labeling, and again just prior to release to ensure that high quality tissue is being made available for transplantation.

#### Sterile Packaging:

SteriGraft™ Allograft Tissues from Bone Bank Allografts are sterilized in their final packaging to ensure the highest patient safety. Our packaging processes ensure that allografts arrive to the surgery in the best possible condition to facilitate a successful transplantation.

#### Quality Assurance:

To comply with the standards set out by the American Association of Tissue Banks (AATB) and the U.S. Food and Drug Administration's (FDA), Good Tissue Practices, Bone Bank Allografts, its processing partner (THB), and the various donor processing organizations we work with all are subjected to an extensive quality process.



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# SteriGraft Allograft

## Cortical Cancellous Allograft

### Mineralized Cortical/Cancellous Bone Powder 250-1000 mic

**Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Bone grafting applications in dental procedures.

|  |          |
|--|----------|
| ORDER CODE 5M25025 SteriGraft - Mineralized Cortical - Cancellous Bone Powder 250-1000 mic, Vol 0.25cc - | \$35.00  |
| ORDER CODE 5M25050 ----- Vol 0.5cc -   | \$53.00  |
| ORDER CODE 5M25100 ----- Vol 1.0cc -   | \$70.00  |
| ORDER CODE 5M25250 ----- Vol 2.5cc -   | \$131.00 |
| ORDER CODE 5M25500 ----- Vol 5.0cc -   | \$197.00 |



### Deminerized Cortical/Cancellous Bone Powder 250-1000 mic

**Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Bone grafting applications in dental procedures.

|  |          |
|--|----------|
| ORDER CODE 7M25025 SteriGraft - Deminerized Cortical - Cancellous Bone Powder 250-1000 mic, Vol 0.25cc - | \$44.00  |
| ORDER CODE 7M25050 ----- Vol 0.5cc -   | \$67.00  |
| ORDER CODE 7M25100 ----- Vol 1.0cc -   | \$90.00  |
| ORDER CODE 7M25200 ----- Vol 2.0cc -   | \$137.00 |
| ORDER CODE 7M25250 ----- Vol 2.5cc -   | \$166.00 |
| ORDER CODE 7M25500 ----- Vol 5.0cc -   | \$215.00 |

### Mineralized Cortical/Cancellous Crushed 1-4mm

**Source:** Proximal and distal ends of long bones, ilium, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Filling of bone void defects, augments prosthetic or other implants used for spinal fusion, oral or maxillofacial reconstruction, or ankle and foot reconstruction.

|  |          |
|--|----------|
| ORDER CODE 449 SteriGraft - Cortical - Cancellous Crushed 1-4mm, Vol 5cc - | \$125.00 |
| ORDER CODE 464 ----- Vol 15cc -  | \$325.00 |
| ORDER CODE 465 ----- Vol 30cc -  | \$425.00 |



### Mineralized Cortical/Cancellous Crushed 4-10mm

**Source:** Proximal and distal ends of long bones, ilium, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Filling of bone void defects, augments prosthetic or other implants used for spinal fusion, oral or maxillofacial reconstruction, or ankle and foot reconstruction.

|  |          |
|--|----------|
| ORDER CODE 459 SteriGraft - Cortical - Cancellous Crushed 4-10mm, Vol 10cc - | \$200.00 |
| ORDER CODE 467 ----- Vol 15cc -  | \$300.00 |
| ORDER CODE 468 ----- Vol 30cc -  | \$400.00 |

## Cortical Allograft

### Mineralized Cortical Bone Powder 250-1000 mic

**Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Bone grafting applications in dental procedures.

|   |          |
|---|----------|
| ORDER CODE 4M25025 SteriGraft - Mineralized Cortical Bone Powder 250-1000 mic, Vol 0.25cc - | \$34.00  |
| ORDER CODE 4M25050 ----- Vol 0.5cc -  | \$51.00  |
| ORDER CODE 4M25100 ----- Vol 1.0cc -  | \$68.00  |
| ORDER CODE 4M25250 ----- Vol 2.5cc -  | \$128.00 |



### Deminerized Cortical Bone Powder 250-1000 mic

|   |          |
|---|----------|
| ORDER CODE 6M25025 SteriGraft - Deminerized Cortical Bone Powder 250-1000 mic, Vol 0.25cc - | \$39.00  |
| ORDER CODE 6M25050 ----- Vol 0.5cc -  | \$58.00  |
| ORDER CODE 6M25100 ----- Vol 1.0cc -  | \$78.00  |
| ORDER CODE 6M25200 ----- Vol 2.0cc -  | \$122.00 |
| ORDER CODE 6M25250 ----- Vol 2.5cc -  | \$147.00 |
| ORDER CODE 6M25500 ----- Vol 5.0cc -  | \$186.00 |

## Cancellous Allograft

### Mineralized Cancellous Crushed 1-4mm

**Source:** Proximal and distal ends of long bones, ilium, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Filling of bone defects, prosthetic or bone graft augmentation for spinal fusion, oral and maxillofacial, ankle and foot, or general orthopaedic reconstruction.

|   |          |
|---|----------|
| ORDER CODE 316 SteriGraft - Cancellous Crushed 1-4mm, Vol 5cc - | \$175.00 |
| ORDER CODE 309 ----- Vol 10cc -                                 | \$275.00 |
| ORDER CODE 308 ----- Vol 15cc -                                 | \$375.00 |
| ORDER CODE 306 ----- Vol 30cc -                                 | \$475.00 |
| ORDER CODE 307 ----- Vol 60cc -                                 | \$875.00 |



### Mineralized Cancellous Crushed 4-10mm

**Source:** Proximal and distal ends of long bones, ilium, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Filling of bone defects, prosthetic or bone graft augmentation for spinal fusion, oral and maxillofacial, ankle and foot, or general orthopaedic reconstruction.

|   |          |
|---|----------|
| ORDER CODE 317 SteriGraft - Cancellous Crushed 4-10mm, Vol 10cc - | \$250.00 |
| ORDER CODE 312 ----- Vol 15cc -                                   | \$350.00 |
| ORDER CODE 313 ----- Vol 30cc -                                   | \$450.00 |
| ORDER CODE 310 ----- Vol 60cc -                                   | \$850.00 |

### Mineralized Cancellous Cubes 5mm

**Source:** Proximal and distal ends of long bones, ilium, femur, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Filling of bone defects, prosthetic or bone graft augmentation for spinal fusion or general orthopaedic reconstruction, oral and maxillofacial reconstruction, or ankle and foot reconstruction.

|  |          |
|--|----------|
| ORDER CODE 318 SteriGraft - Cancellous Cubes 5mm, Vol 10cc - | \$275.00 |
| ORDER CODE 315 ----- Vol 15cc -                              | \$375.00 |
| ORDER CODE 304 ----- Vol 30cc -                              | \$475.00 |
| ORDER CODE 305 ----- Vol 60cc -                              | \$875.00 |




Call Toll Free at 1-800-888-4435

# SteriGraft Allograft

## Dense Cancellous Block 12 x 14 Paralell

**Source:** Proximal / distal tibia, distal femur, femoral head or humeral head. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Spinal fusion or general orthopaedic reconstruction.



|  |            |  |            |
|--|------------|--|------------|
| ORDER CODE 72205 SteriGraft - Dense Cancellous Block 12 x 14 Paralell, 5mm | - \$850.00 | ORDER CODE 72206 SteriGraft - Dense Cancellous Block 12 x 14 Paralell, 6mm | - \$850.00 |
| ORDER CODE 72207 ----- 7mm   | - \$850.00 | ORDER CODE 72208 ----- 8mm   | - \$850.00 |
| ORDER CODE 72209 ----- 9mm   | - \$850.00 | ORDER CODE 72210 ----- 10mm  | - \$850.00 |
| ORDER CODE 72211 ----- 11mm  | - \$850.00 | ORDER CODE 303 SteriGraft - Cancellous Block 30 x 15 x 15mm                | - \$600.00 |

## Soft Tissue Allografts

### Pericardium

**Source:** Pericardial sac. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature.

**Common Surgical Applications:** Pericardial defect, dura matter repair, and periodontal reconstruction.

|   |            |
|---|------------|
| ORDER CODE 720 SteriGraft - Pericardium 15 x 20mm (+/- 0.5mm) | - \$100.00 |
| ORDER CODE 721 ----- 20 x 30mm (+/- 0.5mm)                    | - \$150.00 |
| ORDER CODE 722 ----- 30 x 40mm (+/- 0.5mm)                    | - \$250.00 |
| ORDER CODE 723 ----- 60 x 60mm (+/- 5mm)                      | - \$450.00 |
| ORDER CODE 724 ----- 80 x 80mm (+/- 5mm)                      | - \$650.00 |



### Fascia Lata

**Source:** Lateral surface of the thigh. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Description:** Fascia may or may not include the iliotibial bands. The width (W) will be perpendicular to the grain, and the length (L) will be taken along the grain. **Common Surgical Applications:** Ligament reconstruction, repair of craniotomy and dural defects, eye lid repair, or bladder suspension.

ORDER CODE 701 SteriGraft - Fascia Lata W = 40mm, L = 60mm - \$300



### Amnion Patch

**Source:** Human amniotic membrane. **Preservation Method:** Desiccated (Air Dried). Containers must be stored at room temperature. **Description:** Amnion is a versatile tissue with application in a variety of surgical procedures to provide physical barrier protection to wounds. Literature has documented amnion can reduce inflammation, contains collagen, has a variety of growth factors, has been shown to reduce pain, is biocompatible, and resorbs within the body. Amnion has also been used to reduce the post operative risks associated with scar tissues. SteriShield™ Protective Amnion Patch is an easy to handle, dry graft that rehydrates within the wound setting. Features an orientation notch in the upper right hand corner to ensure proper placement on the wound bed. **Common Surgical Applications:** A wide variety of procedures can benefit from the SteriShield™ Protective Amnion Patch.

|   |            |
|---|------------|
| ORDER CODE 730 SteriShield™ Protective Amnion Patch 1cm x 1cm | - \$199.00 |
| ORDER CODE 731 ----- 2cm x 2cm                                | - \$339.00 |



## Specialty Fashioned Allografts

### Uni-Cortical Cancellous Blocks 12 x 14 Lordosed 5 degrees

**Source:** Femoral head, patella, distal tibia, talus or calcaneus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Cervical spinal fusion.


|   |            |
|---|------------|
| ORDER CODE 70305 SteriGraft - Unicortical Dense Cancellous Block 12 x 14 Lordosed 5 degree, 5mm | - \$850.00 |
| ORDER CODE 70306 ----- 6mm  | - \$850.00 |
| ORDER CODE 70307 ----- 7mm  | - \$850.00 |
| ORDER CODE 70308 ----- 8mm  | - \$850.00 |
| ORDER CODE 70309 ----- 9mm  | - \$850.00 |
| ORDER CODE 70310 ----- 10mm   | - \$850.00 |
| ORDER CODE 70311 ----- 11mm   | - \$850.00 |



### Uni-Cortical Cancellous Blocks 12 x 14 Paralell

**Source:** Femoral head, patella, distal tibia, talus or calcaneus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Common Surgical Applications:** Cervical spinal fusion.

|  |            |
|--|------------|
| ORDER CODE 71305 SteriGraft - Unicortical Dense Cancellous Block 12 x 14 Paralell, 5mm | - \$850.00 |
| ORDER CODE 71306 ----- 6mm   | - \$850.00 |
| ORDER CODE 71307 ----- 7mm   | - \$850.00 |
| ORDER CODE 71308 ----- 8mm   | - \$850.00 |
| ORDER CODE 71309 ----- 9mm   | - \$850.00 |
| ORDER CODE 71310 ----- 10mm  | - \$850.00 |



### Steriflex™ Wrappable Cortical Bone

**Components:** Femur, tibia, fibula, ilium or humerus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Description:** 100% cortical bone that can be bent, contoured, rolled, trimmed, molded or sewn, making this flexible bone material incredibly versatile in surgical settings. Maintains the structure, stability, and osteoconductivity of human bone and stimulates new bone growth. Rehydration of bone implants for such a desired use has been discussed in the literature for years (Mulliken, Ann. Surg., 1981), but only recently has technology allowed us to commercially produce such an innovative solution. SteriFlex™ is sterilized with the GraftShield™ process for improved allograft safety.

**Common Surgical Applications:** Orthopaedic, neurologic, dental / maxillofacial, trauma and bone augmentation.

|  |            |
|--|------------|
| ORDER CODE 72825 SteriFlex™ Wrappable Cortical Bone, 15x25mm | - \$245.00 |
| ORDER CODE 72835 ----- 15x35mm                               | - \$340.00 |
| ORDER CODE 72850 ----- 15x50mm                               | - \$485.00 |
| ORDER CODE 72875 ----- 15x75mm                               | - \$740.00 |
| ORDER CODE 72810 ----- 15x100mm                              | - \$925.00 |
| ORDER CODE 73825 ----- 20x25mm                               | - \$325.00 |
| ORDER CODE 73835 ----- 20x35mm                               | - \$450.00 |
| ORDER CODE 73850 ----- 20x50mm                               | - \$650.00 |
| ORDER CODE 73875 ----- 20x75mm                               | - \$975.00 |



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.


# Sterigraft Allograft



## Sterisorb™ Compressible Cancellous Allograft Strip

**Components:** Proximal and distal ends of long bones, ilium, femoral head, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Description:** 100% cancellous bone that is compressible, has wicking characteristics, and maintains the biocompatibility of human bone. A new solution for "tight space" bone voids, which assumes characteristics of a sponge by absorbing saline, blood, or bone marrow aspirate. The graft compresses to 1/2 its original size (or less) then expands back to its former size to secure a better fit. This gives the surgeon more flexibility, decreases surgery time, reduces the potential for migration, and eliminates the leaching commonly seen with putty and gels. Encourages excellent osteoconductivity and promotes rapid vascularization of the allograft to incorporate the recipient's natural bone growth factors and stimulate new bone growth. Such rehydration of bone implants has been discussed in the literature for years (Mulliken, Ann. Surg., 1981), but only recently has technology allowed us to commercially produce such an innovative solution. SteriSorb™ is sterilized with the GraftShield™ process for improved allograft safety. **Common Surgical Applications:** Spinal fusion and general orthopaedic surgeries; dental reconstruction surgeries - any application where a flexible bone void filler will aid the surgical process.

|                         |  |             |
|-------------------------|--|-------------|
| <b>ORDER CODE 74801</b> | SteriSorb™ Allograft Strip 10 x 50 x 5mm | - \$1050.00 |
| <b>ORDER CODE 74818</b> | ----- 15 x 20 x 5mm                      | - \$1300.00 |
| <b>ORDER CODE 74815</b> | ----- 15 x 40 x 5mm                      | - \$1500.00 |
| <b>ORDER CODE 74830</b> | ----- 20 x 30 x 5mm                      | - \$1700.00 |
| <b>ORDER CODE 74850</b> | ----- 20 x 50 x 5mm                      | - \$2600.00 |
| <b>ORDER CODE 74840</b> | ----- 25 x 40 x 5mm                      | - \$1800.00 |
| <b>ORDER CODE 74805</b> | ----- 5 x 50 x 5mm                       | - \$950.00  |



## Sterisorb Compressible Cancellous Bone Chips 1-4mm

**Components:** Proximal and distal ends of long bones, ilium, femoral head, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Description:** 100% cancellous bone that is compressible, has wicking characteristics, and maintains the biocompatibility of human bone. A new solution for "tight space" bone voids, which assumes characteristics of a sponge by absorbing saline, blood, or bone marrow aspirate. The graft compresses to 1/2 its original size (or less) then expands back to its former size to secure a better fit. This gives the surgeon more flexibility, decreases surgery time, reduces the potential for migration, and eliminates the leaching commonly seen with putty and gels. Encourages excellent osteoconductivity and promotes rapid vascularization of the allograft to incorporate the recipient's natural bone growth factors and stimulate new bone growth. Such rehydration of bone implants has been discussed in the literature for years (Mulliken, Ann. Surg., 1981), but only recently has technology allowed us to commercially produce such an innovative solution. SteriSorb™ is sterilized with the GraftShield™ process for improved allograft safety. **Common Surgical Applications:** Spinal fusion and general orthopaedic surgeries; dental reconstruction surgeries - any application where a flexible bone void filler will aid the surgical process.



|                         |   |             |
|-------------------------|---|-------------|
| <b>ORDER CODE 71801</b> | SteriSorb™ Compressible Bone Chips 1-4mm particles, Vol 1cc | - \$225.00  |
| <b>ORDER CODE 71802</b> | ----- Vol 2cc   | - \$425.00  |
| <b>ORDER CODE 71805</b> | ----- Vol 5cc   | - \$850.00  |
| <b>ORDER CODE 71810</b> | ----- Vol 10cc  | - \$1400.00 |
| <b>ORDER CODE 71815</b> | ----- Vol 15cc  | - \$1850.00 |
| <b>ORDER CODE 71830</b> | ----- Vol 30cc  | - \$3500.00 |



## Sterisorb Compressible Cancellous Cubes

**Components:** Proximal and distal ends of long bones, ilium, femoral head, calcaneus or talus. **Preservation Method:** Lyophilized. Containers must be stored at ambient temperature. **Description:** 100% cancellous bone that is compressible, has wicking characteristics, and maintains the biocompatibility of human bone. A new solution for "tight space" bone voids, which assumes characteristics of a sponge by absorbing saline, blood, or bone marrow aspirate. The graft compresses to 1/2 its original size (or less) then expands back to its former size to secure a better fit. This gives the surgeon more flexibility, decreases surgery time, reduces the potential for migration, and eliminates the leaching commonly seen with putty and gels. Encourages excellent osteoconductivity and promotes rapid vascularization of the allograft to incorporate the recipient's natural bone growth factors and stimulate new bone growth. Such rehydration of bone implants has been discussed in the literature for years (Mulliken, Ann. Surg., 1981), but only recently has technology allowed us to commercially produce such an innovative solution. SteriSorb™ is sterilized with the GraftShield™ process for improved allograft safety. **Common Surgical Applications:** Spinal fusion and general orthopaedic surgeries; dental reconstruction surgeries - any application where a flexible bone void filler will aid the surgical process.

|                         |   |             |
|-------------------------|---|-------------|
| <b>ORDER CODE 70808</b> | SteriSorb™ Cube Compressible Bone Allografts, 8 x 8 x 8mm | - \$850.00  |
| <b>ORDER CODE 70810</b> | ----- 10 x 10 x 10mm                                      | - \$900.00  |
| <b>ORDER CODE 70812</b> | ----- 12 x 12 x 12mm                                      | - \$950.00  |
| <b>ORDER CODE 70814</b> | ----- 14 x 14 x 14mm                                      | - \$1000.00 |



## Demineralized Bone Matrix

### DBM / Demineralized Bone Matrix

Now available from Bone Bank Allografts, a DBM Putty and Crunch with the following features: • 100% Human DBM bone void filler • No filler / carrier necessary due to a unique processing method • Both putty and crunch are flowable and formable • Does not wash away easily • Produced under strict c-GMP and ISO 13485 practices • Recovery and screening performed according to AATB and FDA guidelines • Terminally sterilized using validated gamma irradiation • 2 year shelf life

#### DBM PUTTY / DBM CRUNCH

|                            |                   |             |
|----------------------------|-------------------|-------------|
| <b>ORDER CODE JHG-007P</b> | DBM Putty 0.7cc   | - \$150.00  |
| <b>ORDER CODE JHG-007P</b> | DBM Putty 0.7cc   | - \$160.00  |
| <b>ORDER CODE JHG-02P</b>  | DBM Putty 2.5cc   | - \$370.00  |
| <b>ORDER CODE JHG-05P</b>  | DBM Putty 5.0 cc  | - \$875.00  |
| <b>ORDER CODE JHG-10P</b>  | DBM Putty 10.0cc  | - \$1300.00 |
| <b>ORDER CODE JHG-01C</b>  | DBM Crunch 1.0cc  | - \$160.00  |
| <b>ORDER CODE JHG-02C</b>  | DBM Crunch 2.5cc  | - \$370.00  |
| <b>ORDER CODE JHG-05C</b>  | DBM Crunch 5.0cc  | - \$875.00  |
| <b>ORDER CODE JHG-010C</b> | DBM Crunch 10.0cc | - \$1300.00 |

DBM Crunch is similar to the DBM Putty, with the addition of a larger size particle ( up to 4mm) mixed in for a more voluminous fill.





# Synthetic Bone Grafting Materials

## NanoGen

- The First Nanotechnology based bone graft material
- Particles composed of nanocrystalline calcium sulfate
- Biodegradable, osteoconductive guided bone regeneration material
- Simple to use, Extremely cost effective
- Comes with Setting solution ( sterile saline )



Most bone grafts on the market are microcrystalline. NanoGen is the first nanotechnology-based bone graft. Particles are composed of nanocrystalline calcium sulfate grains manufactured using proprietary technology. As it undergoes controlled degradation, NanoGen encourages mineralization of newly regenerated bone. It is a biodegradable, osteoconductive, guided bone regeneration material that is simple to use, extremely cost effective.

Add saline to NanoGen to form a putty, graft it into the defect and close it - Simple.

NanoGen is a proprietary product developed by Orthogen, LLC. Microcrystalline calcium sulfate is converted to grains of calcium sulfate (size range 200 to 900 nanometers), which are tightly compressed together to form granules in sizes ranging from 400 to 1000 microns. This unique microscopic structure imparts unique properties to the product. NanoGen undergoes controlled degradation over a period of 12 weeks as compared to 4 to 6 weeks for traditional calcium sulfate. As it degrades, it leads to the formation of calcium phosphate which is the main inorganic component of bone and which also serves as a stimulus for further bone regeneration. Osteoblasts attach to the calcium phosphate and deposit more bone in the defect. Since NanoGen degrades over a controlled period as compared to traditional calcium sulfate, it provides the stimulation for bone regeneration over a much longer period. Hence robust bone regeneration is observed in defects grafted with NanoGen. Because it is completely degradable, bone defects grafted with NanoGen are eventually completely replaced with regenerated vital bone. Research shows that availability of vital bone is important for the placement and eventual success of dental implants.

Each box of NanoGen has 2 patient kits. Each kit carries 1 gram of NanoGen material and regular setting solution.

**ORDER CODE Nano-Gen** Nanocrystalline Calcium Sulfate, 2 Gram Kit - \$139.00 or Buy 3 or more - \$129.00



## DentoGen

Calcium Sulfate Hemihydrate Bone Graft Material, Bioresorbable and Osteoconductive

Biocompatible & osteoconductive, Easy to handle and place, Acts as a soft tissue barrier, Cost Effective, 4-6 week Resorbion time.

PROMOTES BONE REGENERATION - IN A WIDE RANGE OF CLINICAL INDICATIONS

Post-extraction • Infra-osseous defects • Apicoectomy • Sinus lifts • Root perforations • Open apices • Dehiscences and fenestrations

Dentogen can be used in combination with other grafts. The handling characteristics and technical performance of other grafts, such as allograft, are improved when used in combination with Dento-Gen.

Each box of DentoGen has 2 patient kits. Each kit has 1 gram of DentoGen material, regular setting solution ( saline ) and fast setting solution ( Potassium Sulfate. )

**ORDER CODE DENT01** Dentogen Calcium Sulfate Hemihydrate, 2 Gram Kit - \$83.00

**ORDER CODE DENT01NFS** Dentogen Calcium Sulfate Hemihydrate, 2 Gram Kit, no fast set solution - \$78.00

## UNIGRAFT®

A sterile, osteoconductive bone mineral for use in periodontal and maxillofacial surgery. Unigraft® is a synthetic bone mineral for use in the repair of osseous defects. It is a bioceramic granule (commonly referred to as bioactive glass) that is easy to use and facilitates new bone formation in osseous defects. After implantation, Unigraft® granules are transformed into a biological hydroxyapatite scaffold that steadily releases bone mineral ingredients, including calcium, phosphorus and silicon. The development and ingrowth of new bone at the implantation site of Unigraft® is favored, due to the affinity of osteoblasts to the biological hydroxyapatite scaffold and the steady supply of bone minerals. Generation of new osseous tissue progresses with time, within and around the scaffold until the biological hydroxyapatite scaffold is replaced and filled with the newly formed bone.

### Unigraft® is recommended for:

- Augmentation or reconstructive treatment of the alveolar ridge
- Filling of infrabony periodontal defects
- Filling of defects after root resection, apicoectomy, and cystectomy
- Filling of extraction sockets to enhance preservation of the alveolar ridge
- Elevation of the maxillary sinus floor
- Filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration (GTR) and Guided Bone Regeneration (GBR)
- Filling of per-implant defects in conjunction with products intended for Guided Bone Regeneration (GBR)
- Filling of large oral and maxillofacial osseous cavities

ADVERSE REACTIONS: No adverse reactions have been reported



**UNIGRAFT, Osteo Conductive Synthetic Bone Mineral**

**ORDER CODE V03-0105** .37 Gram - \$165.00 5 Per Pack

**ORDER CODE V03-0205** 1 Gram - \$265.00



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.



# Synthetic Bone Grafting Materials



## SynOss™ Putty - Synthetic Mineral Collagen Composite Bone Graft

SynOss™ Putty Synthetic Mineral Collagen Composite is a bone graft matrix with an additional characteristic that enables it to become moldable putty upon hydration. It is indicated for use in oral surgical applications involving bone repair such as augmentation or reconstructive treatment of the alveolar ridge and for the filling of periodontal defects in conjunction with products intended for Guided Tissue Regeneration and Guided Bone Regeneration. SynOss™ putty is calcium phosphate based mineral with a carbonate apatite structure similar to natural bone combined with Type I collagen derived from bovine Achilles tendon. The mineral particles are dispersed within collagen fibers forming a 3-dimensional matrix. It is supplied dry and forms a moldable putty upon hydration. It is fully resorbed during the natural process of bone formation and remodeling.


**Product Features:**

- Carbonite apatite plus Type I collagen
- Moldable putty upon hydration
- Osteoconductive, 3 dimensional structure for bone forming cells
- Mineral particle are contained within the defect site

- ORDER CODE SMB050** 9.0mm x 8.0mm 0.5cc - **\$106.00**
- ORDER CODE SMB100** 11.0mm x 10.5mm 1.0cc - **\$170.00**
- ORDER CODE SMB200** 11.0mm x 21.0mm 2.0cc - **\$240.00**



SynOss™ Putty Synthetic Mineral Collagen Composite Bone Graft, hydrated



## SynOss™ Granules - Synthetic Mineral Bone Graft

SynOss™ Granules Synthetic Mineral Bone Graft is an osteoconductive calcium phosphate based bone graft material intended for use in oral surgical applications involving bone repair such as augmentation or reconstructive treatment of the alveolar ridge, filling of periodontal defects, filling of defects after root resection or extraction sockets, and elevation of maxillary sinus floor.

**Product Features:**

- Synthetic alternative to allografts and xenografts
- Carbonate apatite structure similar to natural bone
- Osteoconductive 0.35 - 1.0mm particle size granules

SynOss™ Granules provide clinicians and their patients with an ideal alternative to human allograft and animal origin bone graft material. The similarity in structure between SynOss™ Granules and natural bone mineral allows the SynOss™ Granules in vivo resorption and remodeling profile to mimic that of natural bone.

- ORDER CODE SM2510** 0.25g 0.5cc - **\$45.00**
- ORDER CODE SM5010** 0.5g 1.0cc - **\$63.00**
- ORDER CODE SM10010** 1.0g 2.0cc - **\$117.00**




# Osteo Tape

## OsteoTape®

OsteoTape® is OsteoGen® with a Collagen Resorbable Bone Graft Matrix in porous preformed bone grafting shapes comprised of highly purified Type I bovine Achilles tendon collagen, combined with crystals of the product OsteoGen®, a synthetic bioactive resorbable graft of the non-ceramic hydroxylapatite category. OsteoTape® simulates the natural collagen matrix and mineral porous structure of human bone. OsteoTape® eliminates the need for secondary surgery that is required to remove non-resorbable (PTFE) membrane. OsteoTape® bone graft in preformed shapes is designed to prevent the problematic migration of particulate bone graft products to major anatomical landmarks (e.g., mandibular canal, sinus perforations, trigeminal nerves and mental foramen).

OsteoTape®, a synthetic bioactive resorbable graft, is an osteoconductive non-ceramic form of hydroxylapatite, the major mineral component of bone and dental enamel. OsteoGen® [Ca<sub>5</sub>(PO<sub>4</sub>)<sub>3</sub>(OH)] has been used for the augmentation and repair of bone defects since 1984. OsteoGen® is a highly crystalline material, with higher purity and crystallinity than ceramic HA as specified by ASTM standards. OsteoGen® has no tricalcium phosphates, amorphous phases, or bone inhibiting pyrophosphates, as found in ceramic HA products which are devoid of the hydroxyl group, due to their sintering process. However, the hydroxyl group is found in human bone and OsteoGen®. By x-ray diffraction, OsteoGen® has been shown to meet and exceed ASTM standards and specifications. OsteoGen® bone graft by dental x-ray will show as radiolucent the day of placement in defect sites. Four to six months later, the defect site will be radiopaque depending on vascularity, patient's age and quantity of product delivered. In many animal and clinical studies, OsteoGen® has been shown to control migration of connective tissue.



### ORDER CODE

**Ot5215-2** OsteoTape Preformed Porous Resorbable Bone Graft 5x21x5mm, 2 per- **\$130.00**

**Ot2060-2** OsteoTape Preformed Porous Resorbable Bone Graft 20x60x2.5mm, 2 per- **\$178.00**

**Ot1621-4** OsteoTape Preformed Porous Resorbable Bone Graft 16x21x1.3mm, 4 per box- **\$272.00**

**Ot2133-4** OsteoTape Preformed Porous Resorbable Bone Graft 21x33x1.3, 4 per- **\$300.00**

Also available in OsteoRidge - Preformed Ridge Shape:

**Otr625-2** OsteoTape OsteoRidge Preformed Porous Resorbable Bone Graft, 2per- **\$130.00**.

# Algisorb™



## Algisorb™

Algisorb, Formerly C-Graft - Your perfect alternative to autogenous bone.

An attractive alternative to autogenous bone, allografts, and xenografts for preserving and grafting alveolar bone, algisorb™ (formerly known as C-Graft™ or Algipore™) is a completely natural bone-grafting material derived from the skeleton of marine algae, a renewable plant harvested off the coasts of South Africa and England. When dried and processed according to strict Food and Drug Administration (FDA) guidelines and Good Manufacturing Procedures, the harvested algae is transformed into calcium phosphate, which can provide osteoconductive scaffolding for bone regeneration in grafting procedures.

According to Dr. Boyd Tomasetti, DMD, Assistant Clinical Professor at the University of Colorado School of Dentistry, a Board Certified Oral/Maxillofacial Surgeon, and past President of the American Association of Oral and Maxillofacial Surgeons, and the Colorado Society of Oral and Maxillofacial surgeons, the processing of the algae eliminates any type of contamination or proteins present in the plant form, making the finished product completely safe for human implantation. "The clinical applications of C-Graft™ relate to dental implant surgery specifically", Tomasetti explains. "I have begun using C-Graft™ almost exclusively in my practice. It certainly has been shown to be superior to bovine derived material and, as an oral surgeon, I have found it to be an excellent alternative to autogenous bone."

### ORDER CODE

**OTA501** Algisorb - .5ml vial - **\$80.00**

**OTA502** Algisorb - 1ml vial - **\$125.00**

**OTA503** Algisorb - 2ml vial - **\$195.00**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Collagen Matrix™




## MatrixDerm™

### Product Features:

- Highly purified intact collagen
- High mechanical strength and structural integrity for easy handling
- Highly biocompatible
- Permeable to macromolecules and nutrients while creating an effective barrier to epithelial cells
- Conformable and repositionable for precise adjustment and placement
- Resorbable in 6-9 months, eliminating the need for removal surgery
- Terminally sterilized and ready for use following brief hydration

MatrixDerm™ Regenerative Collagen Dental Membrane is a nonfriable, conformable, resorbable, membrane matrix consisting of highly purified type I and III collagen derived from intact porcine dermis.



MatrixDerm™ has been designed to simultaneously achieve drapability, structural integrity, and the ability to reposition the membrane for ideal placement prior to closure. MatrixDerm™ conforms to the defect site, but also has sufficient memory to maintain space for bone and tissue ingrowth.

The programmed resorption time of MatrixDerm™ supports bone or tissue regeneration for 6 to 9 months to allow for ingrowth and remodeling of the defect site. The semipermeable membrane allows for nutrient exchange while providing a cell barrier to prevent epithelial down growth.

**ORDER CODE PDM1520** MatrixDerm™ Regenerative Collagen Dental Membrane, 15x20mm - **\$115.00**


**ORDER CODE PDM2030** MatrixDerm™ Regenerative Collagen Dental Membrane, 20x30mm - **\$170.00**

**ORDER CODE PDM3040** MatrixDerm™ Regenerative Collagen Dental Membrane, 30x40mm - **\$250.00**

## MatrixDerm Cap: Regenerative Collagen Dental Membrane

### Product Features:

- Highly purified intact collagen
- Resorbable in 9-12 months
- Highly biocompatible
- Non-pyrogenic
- Curves when hydrated to facilitate placement
- Conformable and repositionable



MatrixDerm™ Cap Regenerative Collagen Dental Membrane Cap is a white, nonfriable membrane matrix engineered from highly purified porcine dermis. It is indicated for oral surgical procedures as a resorbable material for use in extraction sockets and small bone defects.

MatrixDerm™ Cap was specifically designed to meet the needs for socket ridge preservation procedures, and it has a slight curve when hydrated to facilitate placement. It is resorbable in 9-12 months for extended stability. The size, shape, strength, in vivo stability and conformability are particularly useful for guided bone regeneration in socket preservation.


## MatrixMem: Regenerative Collagen Dental Membrane

Resorbable in 4-6 months, Shorter Resorption time

**ORDER CODE CMPT1520** Regenerative Collagen Dental Membrane, 15mm x 20mm - **\$106.00**

**ORDER CODE CMPT2030** Regenerative Collagen Dental Membrane, 20mm x 30mm - **\$153.00**

**ORDER CODE CMPT3040** Regenerative Collagen Dental Membrane, 30mm x 40mm - **\$225.00**



## MatrixDerm Ext: Regenerative Collagen Dental Membrane

Resorbable in 9-12 months, Extended Resorption time

**ORDER CODE PDMX1520** Regenerative Collagen Dental Membrane, 15mm x 20mm - **\$137.00**

**ORDER CODE PDMX2030** Regenerative Collagen Dental Membrane, 20mm x 30mm - **\$187.00**


**ORDER CODE PDMX3040** Regenerative Collagen Dental Membrane, 30mm x 40mm - **\$275.00**





# OTA Resorbable Collagen Biomaterial

- Made from highly purified Type 1 Bovine Collagen
- Designed to contain graft materials to prevent migration
- Repair tissue tears
- Resorb and remodel naturally
- Adjustable stiffness to match function
- Shapeable
- Approximate 6 week resorption time



**CASE STUDY**  
Small capsule is ideal for buccal plate defects.

## Socket Capsule: Site Preservation via socket approach

- Can be easily trimmed during surgery for a custom fit
- Can be wet with saline or blood from the local site
- Resorbs as connective tissue grows into and through the collagen
- Provides containment of graft material in its intended position
- Immediately after tooth extraction for site preservation
- Regenerate deficient or fractured labial/buccal plate after tooth extraction
- Approximately 7mm diameter

**ORDER CODE OTASOCKET - \$105.00**



- 1 Tooth with thin labial/buccal plate
- 2 Empty socket
- 3 Fractured and incomplete labial/buccal plate
- 4 Surgeon cut and shaped Socket Capsule, delivered by Sinus Delivery Instruments
- 5 Bone graft and sutures applied

## Sinus Capsule: Maxillary Sinus Elevation via Crestal Entry


- Various sizes for socket or sinus approach
- Can be wet with saline or blood from the local site
- Can be easily trimmed during surgery for a custom fit
- Repairs tears in sinus membrane
- Wound healing
- Reinforce soft tissues
- Retains bone graft materials in its intended position

**Available in 2 sizes:**

OTA Medium- Approx 11mm in diameter

Large- Approx 16mm in diameter

**ORDER CODE OTASINUS - \$105.00**




- 1 Maxillary Sinus
- 2 Schneiderian Membrane
- 3 Alveolar Ridge
- 4 Crestal incision
- 5 Balloon for Sinus Lifts
- 6 Surgeon cut and shaped Sinus Capsule, delivered by Balloon for Sinus Lifts
- 7 Bone graft and sutures applied

## Wedge- Caldwell-Luc Sinus Lift

- Thicker reinforced edge for attachment to lateral wall bone with tacks
- Can be wet with saline or blood from the local site
- Can be easily trimmed during surgery for a custom fit
- Repairs tears in sinus membrane
- Wound healing
- Creates new sinus floor during major tears
- Approx 23mm x 38mm

**ORDER CODE OTAWEDGE - \$105.00**



- 1 Maxillary Sinus
- 2 Schneiderian Membrane
- 3 Alveolar Ridge
- 4 Lateral Osteotomy
- 5 Tear in Schneiderian Membrane
- 6 Surgeon cut and shaped Wedge
- 7 Bone tack
- 8 Bone graft and sutures applied



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Resorbable & Non Resorbable Barrier Membranes & Titanium Mesh



## Cytoflex Resorb Resorbable Synthetic Barrier Membrane

For reconstruction of the alveolar ridge and periodontal defects

- Triple layered membrane
- High regenerative capacity
- Excellent gingival compatibility
- Ideal dimensional stability
- Microporous surface texture enhances host tissue adhesion
- Flexible barrier membrane adapts easily to tissue contours

Cytoflex resorb exhibits robust functions for GTR and GBR procedures. While the barrier membrane prevents tissue penetration, it was also designed to allow rapid nutrient permeation without hindrance. Human gingival fibroblasts attach to and proliferate naturally on the Cytoflex Resorb surface without losing cellular viability. The barrier membrane is highly adaptable and exhibits excellent dimensional stability in vivo.

**C03-0101** CYTOFLEX RESORB, Resorbable Synthetic Barrier Membrane, 30 x 40mm, 1 per pack - **\$110.00**

**C03-0102** CYTOFLEX RESORB, Resorbable Synthetic Barrier Membrane, 30 x 40mm, 2 per pack - **\$200.00**

**C03-0201** CYTOFLEX RESORB, Resorbable Synthetic Barrier Membrane, 20 x 25mm, 1 per pack - **\$75.00**

**C03-0202** CYTOFLEX RESORB, Resorbable Synthetic Barrier Membrane, 20 x 25mm, 2 per pack - **\$140.00**



## Cytoflex Mesh Precision Titanium Micro Mesh

- Pure Phase, golden color titanium micro mesh
- Easy to trim and conform to 3 dimensional contours
- Excellent compatibility with soft and hard tissue
- Fixable with bone tacks or bone screws

### Applications:

- Reconstruction of the alveolar ridge
- Reconstruction of extraction sockets to preserve the alveolar ridge
- GBR or GTR treatment of periodontal or peri-implant defects

Cytoflex Mesh is designed to facilitate the three dimensional reconstruction of the alveolar ridge. The pliable titanium mesh conforms easily to tissue contours, yet offers sufficient stiffness to retain the graft material in desirable shape when osseous tissue regenerates. Made of pure titanium foil, Cytoflex Mesh exhibits excellent compatibility with both hard and soft tissue. The thin and smooth profile of Cytoflex Mesh reduces the irritation to the surrounding sensitive tissue and minimizes potentially early exposure. Each mesh is comprised of numerous micro pores that facilitate the flap attachment and the regeneration of osseous tissue. The golden color of Cytoflex Mesh is aesthetically compatible when placed underneath the gingival tissue flap. Cytoflex Mesh restores the alveolar bone ridge by retaining the augmentation material and providing a favorable environment for complete, healthy new bone growth with desired contour.

**ORDER CODE C02-0101** CYTOFLEX MESH, Precision titanium mesh, 25mm x 30mm, .05mm thickness, .25mm pore size, 1 per pack - **\$50.00**

**ORDER CODE C02-0105** CYTOFLEX MESH, Precision titanium mesh, 25mm x 30mm, .05mm thickness, .25mm pore size, 5 per pack - **\$200.00**



## Cytoflex Micro Porous ePTFE Barrier Membrane

For reconstruction of the alveolar ridge and periodontal defects

- Excellent handling and rigidity
- Easily adaptable
- Easily retrieved as one piece
- Resists fibroblasts and bacteria penetration

### Advantages of a Micro designed Barrier Surface

- Better host tissue attachment with fewer flap dehincens
- Facilitates tissue regeneration despite flap recessions or incomplete primary closure
- Easily retrievable

Cytoflex non-resorbable barrier membrane is made of micro porous expanded PTFE. The micro porosity of Cytoflex is designed to resist the penetration of fibroblasts and bacteria, while simultaneously allowing the exchange of interstitial fluids through the membrane. The unique micro porous design enhances the gingival tissue attachment and provides a favorable protected environment for neo-vascularization and repopulation of osteoblasts in the bony defects. The flexible material easily conforms to the bony defects of grafted site, and can be easily retrieved as one piece after completion of bone growth. The bacteria-impenetrable membrane protects the tissue regeneration site despite flap recession or if primary closure is not obtained.

**ORDER CODE C01-0101** CYTOFLEX TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 25mm x 30mm, 1 pack - **\$45.00**

**ORDER CODE C01-0105** CYTOFLEX TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 25mm x 30mm, 5 pack - **\$200.00**

**ORDER CODE C01-0301** CYTOFLEX TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 12mm x 24mm, 1 pack - **\$35.00**

**ORDER CODE C01-0305** CYTOFLEX TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 12mm x 24mm, 5 pack - **\$160.00**

### Textured

**ORDER CODE C04-0101** CYTOFLEX TEXTURED TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 12mm x 24mm, 1 pack- **\$45.00**

**ORDER CODE C04-0105** CYTOFLEX TEXTURED TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 12mm x 24mm, 5 pack- **\$200.00**

**ORDER CODE C04-0201** CYTOFLEX TEXTURED TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 25mm x 30mm, 1pack- **\$55.00**

**ORDER CODE C04-0205** CYTOFLEX TEXTURED TEF-GUARD, Non-Resorbable Synthetic Barrier Membrane, 25mm x 30mm, 5 pack- **\$250.00**



Call Toll Free at 1-800-888-4435

# Regenerated Cellulose Membrane & Wound Care



## Benacel Hemostatic Gauze

### Advantages

- Completely resorbable within days
- Prevents dry socket formation
- Non-allergenic and non-immunogenic
- Offered in several sizes for various wound configurations

Benacel dental dressing is made from biocompatible regenerated cellulose and contains no chemical additives. Benacel dental dressing adheres to moist oral mucosa and forms a barrier, protecting the wound from further irritation and pain. The material dissolves in a few days and is safe if swallowed. When placed in the extraction socket, Benacel dental dressing forms a gelatinous scaffold in the blood extrudate, facilitating the development of a stable blood clot to prevent dry socket formation. Benacel dental dressing is suitable for use as a wound dressing in extraction sites and the management of alveolar osteitis (dry socket), and may be used as a wound dressing for the temporary management of oral surgical wound, such as operative, postoperative, donor sites, and traumatic injuries. It may also be used as a wound dressing for the management and protection of oral lesions, including sores, ulcers, and injuries, such as cuts, lacerations, and abrasions or the oral mucosa.

### Sizes and Applications:

#### ORDER CODE

- C-001 10 Per Box, 5mm x 7mm plug** - for deep wounds, biopsy sites and post extraction sockets - **\$32.00**
- C-002 12 Per Box, 15mm x 15mm** - for small oral lesions, including sores, ulcers, cuts, lacerations, and abrasions of the oral mucosa - **\$44.00**
- C-003 12 Per Box, 6mm x 8mm** - for small oral lesions, including sores, ulcers, cuts, lacerations, and abrasions of the oral mucosa - **\$34.00**
- C-005 10 Per Box, 50mm x 50mm** - For oral surgical wounds, such as operative, postoperative, donor sites and traumatic injuries - **\$60.00**

## EMED Hemostatic Gauze

EMED Soluble Hemostatic Gauze is indicated for the use anytime blood is present due to cuts, lacerations or tooth extractions. EMED Soluble Hemostatic Gauze Pads provides fast-acting control of bleeding and minimizes blood loss. Composed natural cellulose, EMED Soluble Hemostatic Gauze Pads rapidly absorbs blood and transforms into a clear gel to seal the wound with a protective transparent layer, bind and activate platelets, provide support for clot formation, and create an environment for wound healing.

- Made Of 100% Natural Cellulose
- Cut or Folded To Desired Size & Shape - Acts In Seconds
- Reduces Healing Time And Recurring Bleeding
- Usually Painless Application & Removal
- Becomes Transparent For Non-Disruptive Wound Monitoring

Packaged 40 dressings per box. 20 per foil package – 2 foil packages per box. Dressings measure ½" x 1 ½" each.

**ORDER CODE EHEM-42-1EMED Soluble Hemostatic Gauze, 40 per box, ½" x 1 ½" each - \$49.99**

## Act Cel Hemostatic Gauze

Advanced Bleeding Control. For Extraction sockets, Dry socket prevention, Orthognathic surgery, Trauma, Biopsies, Periodontal surgery, Endodontic surgery, Other oral surgery procedures.

**ORDER CODE HG020022, 2 x 2 Size, 20 Per Box - \$114.00**

## Hemostatic Solution

- Topical Astringent for control of gingival bleeding
- Ephinephrine FREE - No Cardiac Reaction
- Contains Aluminum Chloride-25%

**ORDER CODE H221-1 Hemostatic solution 15cc- \$5.78**

**ORDER CODE H221-2 Hemostatic solution 30cc- \$8.74**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Collagen, Wound Care & Dressings



## Collaform Singles

- Absorption time: 4-6 weeks on average
- Absorbable collagen for tissue preservation
- Maintains graft in extraction site for ridge preservation
- Socket grafting without primary closure allowing keratinized tissue preservation

**Size:** 12mm x 20mm x 3mm each

**ORDER CODE CFS1220**-5 Collaform singles bovine collagen wound dressing, 5 per box -**\$185.00**



## Collaform Plugs

- Absorption time: 4-6 weeks on average
- Absorbable collagen wound dressing
- Economical, easy to handle

**Size:** 10mm x 20mm each

**ORDER CODE CFP1020** - 5 Collaform plug bovine collagen, 5 per box - **\$70.00**



## Reso-Pac

### Hydrophilic surgical dressing, 2 gram cups

Adhesive cellulose based hydrophilic periodontal dressing. An essential aid in implantology, prosthetics, and periodontology.

**Protects:** Seals sulcus and approximate areas Ideal for sealing of pockets, prosthetics and periodontology

**Lasting adhesion:** Dressing remains in place for up to 3 hours. During that time the treated area is completely separated from the oral environment. The adhesive paste safely secures the marginal periodontium.

**Uncomplicated:** Hydrophilic, ready for use dressing, does not have to be mixed, for quick and easy dispensation.

**Patient friendly:** Reso-PAC remains soft at all times. For this reason pressure spots and insufficient supply of the periodontium or the papillae, as well as pressure, ulcers, do not occur. Neutral in smell and taste and the patient does not feel traumatized as is the case when hard dressings are used.

**Dissolves:** The dressing slowly dissolves by itself and does not need to be removed.

**Good value:** 5 year shelf life. Very economical

**ORDER CODE 155019** UNI-DOSE 2gram cups, 50 per jar - **\$43.99**

**ORDER CODE 155018** UNI-DOSE 2 gram cups, 10 per jar - **\$9.95**



Good applicability of Reso-Pac in hard-to-reach areas.



The dressing after 24 hours.



Area covered with Reso-Pac for additional protection.



Reso-Pac dissolves away without leaving

## Dressol-X, Multi-Use Medicated Surgical Dressing

### Indicated uses:

- Dry Socket Treatment
- Surgical Exposure
- Incision and Drainage
- Post Operative Bleeding
- Prevention of Poor Healing Sockets

### MEDICATED CUSHIONING FOR IMMEDIATE DENTURES!!

It is radiopaque, and non-absorbable. (Active ingredients: acetylsalicylic acid, eugenol) **40" strip**

**ORDER CODE D502** Dressol-X Multi-Use Medicated Surgical Dressing - **\$73.99**



## Lukens™ Bone Wax

Lukens™ BoneWax has been a tradition in the medical industry for nearly 100 years. The peel-open foil packet is specifically designed for easy dispensing and provides optimum safety and convenience for use in the operating room.

**ORDER CODE 901** Yellow Bone Wax, 2.5 grams each - **\$66.65**

**ORDER CODE 902** Beige Bone Wax, 2.5 grams each - **\$66.65**

**ORDER CODE 903** White Bone Wax, 2.5 grams each - **\$66.65**




Call Toll Free at 1-800-888-4435




# Implant Models / Case Presentation

## Combination Models



**S63-RRR**  
Upper arch in crystal clear with natural shaded teeth. Central tooth on implant and a 3 unit bridge supported by 3 implants.  
**ORDER CODE S63-RRR - \$209.99**

**S63-RR**  
Upper arch in crystal clear with natural shaded teeth. Central tooth on implant and a 3 unit bridge supported by 2 implants.  
**ORDER CODE S63-RR - \$198.99**



**S75-BR**  
Upper arch in pink with natural shaded teeth. Central tooth on implant. Window shows 3 unit bridge supported by 3 implants. Opposite side shows 3 unit bridge spanning natural abutments.  
**ORDER CODE S75-BR - \$209.99**

**S75-WRR**  
Upper arch in pink with natural shaded teeth. Central tooth on implant + 3 unit bridge, implant supported, shown in window.  
**ORDER CODE S75-WRR - \$209.99**



**SOB-30**  
Lower combination model in crystal clear. One side has ball attachment in implant and other side has ball attachment in prepared tooth.  
**ORDER CODE SOB-30 - \$129.99**

**SCIM-53BB**  
Lower combination model in crystal clear with pink overdenture. On side has 2 "O" rings which clip onto ball anchors. The other side has "Hader" bar and 1 embedded clip.  
**ORDER CODE SCIM-53BB - \$269.99**



**S85-WRR**  
Upper arch on pink base with natural shaded teeth. Central tooth on implant + 3 unit bridge, implant supported. Buccal side is natural pink. Implants are visible from palatal window only.  
**ORDER CODE S853-WRR - \$247.99**



**SCB-10**  
A variety of restorations all on one model-chrome and acrylic removable partial denture-key and keyway precision attachment-ceramic onlay-laminate veneer-post & core-ceramic crown-porcelain bonded to metal 3 unit bridge-single metal crown-implant in window.  
**ORDER CODE SCB-10 - \$304.99**

**S75-STI**  
Upper arch in pink with natural shaded teeth. Central tooth on implant. Window shows 2nd premolar on implant. Opposite side shows 3 unit bridge spanning natural abutments.  
**ORDER CODE S75-STI - \$302.99**



**SCIM-7**  
Lower combination model crystal clear. Half overdenture with 2 "O" rings, clips onto 2 ball anchors supported by implants and 3 unit bridge is implant supported.  
**ORDER CODE SCIM-7 - \$189.99**

**SCIM-5**  
Lower combination model in crystal clear. Half overdenture clips onto bar supported by 2 implants. 3 unit bridge is supported by 2 implants. Single tooth is implant supported.  
**ORDER CODE SCIM-5 - \$329.99**



# Implant Models / Case Presentation



## Assorted Models



**S33-WRR**

Lower arch in pink with natural shaded teeth + 3 unit bridge, implant supported, shown in clear window.

**ORDER CODE S33-WRR - \$199.99**

**S33-CR**

Lower arch in crystal clear with natural shaded teeth, with 3 unit bridge supported by 2 implants.

**ORDER CODE S33-CR - \$197.99**

**SSTR-1**

Upper quadrant in crystal clear with natural shaded teeth, with 1 implant in Central site.

**ORDER CODE SSTR-1 - \$82.99**

**SSTR-5**

Upper quadrant in crystal clear with natural shaded teeth, and visible anatomical roots and 1 implant in Central site.

**ORDER CODE SSTR-5 - \$94.99**

**SQC-35**

Quadrant Cosmetic before and after model, with 1 side showing discolored & imperfect teeth and the other side restored with 1 removable veneer.

**ORDER CODE SQC-35 - \$65.99**

## Perio And Common Disorder Models



**SPMH-1**

Upper hinged gingiva revealing gum recession and bone loss.

**ORDER CODE SPMH-1 - \$26.99**

**SPMH-2**

Upper hinged gingiva revealing calculus - plaque accumulation and increased bone loss.

**ORDER CODE SPMH-2 - \$26.99**

**SPMH-3**

Upper hinged gingiva revealing gross deep calculus and plaque with hard and soft tissue loss.

**ORDER CODE SPMH-3 - \$26.99**



**SMTD-7**

Multiple tooth disease model in crystal clear with natural shaded teeth and anatomical roots, demonstrating a variety of disorders: missing - shifting teeth - maxillary extrusion - fractured root tips - impacted molars - calculus buildup - gum erosion - abscess formation.

**ORDER CODE SMTD-7 - \$129.99**



**SPTL-A**


Hinged progressive tooth loss model demonstrates progressive deterioration and ultimate tooth loss from lack of treatment. In crystal clear with natural shaded teeth and anatomical roots.

**ORDER CODE SPTL-A - \$104.99**




# Implant Models / Case Presentation

## Overdenture & Implant Models




**S85-4B**  
Upper quadrant in crystal clear with natural shaded teeth, and visible anatomical roots and 1 implant in Central site.  
**ORDER CODE S85-4B - \$199.99**

**S85-FB**  
Upper arch in crystal clear and pink denture with 5 embedded clips that snap onto bar supported by 6 implants.  
**ORDER CODE S85-FB - \$399.99**




**S85-6M**  
Upper arch in crystal clear with pink denture supported by 6 "O" rings, which clip onto ball anchors supported by implants.  
**ORDER CODE S85-6M - \$209.99**

**S85-3M**  
Upper arch in crystal clear with pink denture supported by 3 "O" rings, which clip onto ball anchors supported by implants.  
**ORDER CODE S85-3M - \$169.99**




**SOCBI-4C**  
Lower arch in crystal clear and pink denture with 4 embedded clips that snap onto bar supported by 4 implants.  
**ORDER CODE SOCBI-4C - \$367.99**

**SOCBI-5C**  
Lower arch in crystal clear and pink denture with 5 embedded clips that snap onto bar supported by 4 implants.  
**ORDER CODE SOCBI-5C - \$387.99**



**SOCBI-30**  
Upper arch in crystal clear with pink denture supported by 3 "O" rings, which clip onto ball anchors supported by implants.  
**ORDER CODE SOCBI-30 - \$367.99**

**SOCB-2C**  
Lower arch in crystal clear and pink denture with 1 embedded clip that snaps onto bar supported by 2 implants.  
**ORDER CODE SOCB-2C - \$209.99**



**SOVIC-2B**  
Lower arch in crystal clear. Pink overdenture with 2 "O" rings clips onto ball anchors supported by 2 implants.  
**ORDER CODE SOVIC-2B - \$159.99**

**SOVIC-4B**  
Lower arch in crystal clear. Pink overdenture with 4 "O" ring clips onto ball anchors supported by 4 implants.  
**ORDER CODE SOVIC-4B - \$199.99**



**SFRL-6N**  
Lower arch in crystal clear with 12 unit cement retained bridge supported on 6 implants.  
**ORDER CODE SFRL-6N - \$369.99**

**SFRL-6SR**  
Lower arch in crystal clear with 12 unit ridgetop bridge supported on 6 implants and occlusal screws.  
**ORDER CODE SFRL-6SR - \$369.99**



# Implant Models / Case Presentation



## Overdenture & Implant Models *cont.*



**SOVIC-3B**  
Lower arch in crystal clear. Pink overdenture with 3 "O" rings clips onto ball anchors supported by 3 implants.  
**ORDER CODE SOVIC-3B - \$169.99**




**SHWB-5N**  
Lower arch in crystal clear, high water bridge model with a cemented 24 unit denture bridge supported on 5 implants.  
**ORDER CODE SHWB-5N - \$219.99**



**SFRU-16U**  
Upper arch in crystal clear with fixed removable 16 unit ridgetop bridge supported by 6 implants and occlusal screws.  
**ORDER CODE SFRU-16U - \$799.99**



**SHWB-5SR**  
Lower arch in crystal clear. Pink overdenture with 4 "O" ring clips onto ball anchors supported by 4 implants.  
**ORDER CODE SHWB-5SR - \$219.99**



**SHWB-4SR**  
Lower arch in crystal clear, high water bridge model with a 12 unit denture bridge supported on 4 implants and occlusal screws. **ORDER CODE SHWB-4SR - \$199.99**

**SHWB-6SR**  
Lower arch in crystal clear, high water bridge model with a 14 unit denture bridge supported on 6 implants and occlusal screws. **ORDER CODE SHWB-6SR - \$229.99**

**SHWB-6N**  
Lower arch in crystal clear, high water bridge model with a cemented 14 unit denture bridge supported on 6 implants.  
**ORDER CODE SHWB-6N - \$229.99**

## Locator Models



**SULOC-4**  
Upper arch in crystal clear with pink overdenture featuring 4 implants with LOCATOR abutments and LOCATOR attachments. The female abutment is threaded on and may be removed using the genuine tool supplied by ZEST ANCHORS. Comes in crystal clear case.  
**ORDER CODE SULOC-4 - \$356.99**



**SLOC-2**  
Lower arch in crystal clear with pink overdenture featuring two implants with LOCATOR abutments and LOCATOR attachments. It is an excellent tool for demonstrating the retentive qualities of the LOCATOR. The female abutment is threaded on and may be removed using the genuine tool supplied by ZEST ANCHORS. Comes in crystal clear case.  
**ORDER CODE SLOC-2 - \$261.99**



**SLOC-4**  
Lower arch in crystal clear with pink overdenture featuring four implants with LOCATOR abutments and LOCATOR attachments. It is an excellent tool for demonstrating the retentive qualities of the LOCATOR. The female abutment is threaded on and may be removed using the genuine tool supplied by ZEST ANCHORS. Comes in a crystal clear case.  
**ORDER CODE SLOC-4 - \$356.99**




**SLOC-4Bar**  
Lower arch in crystal clear featuring 4 LOCATOR attachments supported by a Hader Bar. This model may be ordered with the posterior implants straight or angled towards the median line. Comes in a crystal clear display case.  
**ORDER CODE SLOC-4Bar - \$409.99**





# Implant Models / Case Presentation

## Mandibles And Maxillas



### SMM4-P

Set of 4 mounted mandibles from full dentition to extreme atrophied ridge and deteriorating condyles - demonstrating the problem with lower dentures and possible need for implants.

ORDER CODE SMM4-P - \$69.99

### SMM4-X

Set of 4 mounted maxillas from full dentition to extreme atrophied ridge demonstrating the problem with upper dentures and possible need for implants.

ORDER CODE SMM4-X - \$98.99



### SL-4 Sinus Lift

Sinus Lift model shows membrane and is useful in demonstrating sinus lift procedures.

ORDER CODE SL-4 - \$59.99



### SFDM-RUB

Free drilling Mandible with rubber coating simulating tissues suitable for hands on procedures such as suturing.

ORDER CODE SFDM-RUB - \$24.99

12ct. - \$22.99

24ct. - \$20.99

### SFDM

Free drilling Mandible suitable for hands on procedures such as drilling and implant placing.

ORDER CODE SFDM - \$12.99

### SPML

Mandible with Perio disorders featuring various types of interosseous defects around the dentition as well as Lingual Tori and furcation.

ORDER CODE SPML - \$29.99

### SPLM-11

Perio lesion mandible featuring vertical & horizontal Perio defects - furcation - and atrophied ridge.

ORDER CODE SPLM-11 - \$29.99



### SBG-5 Bone Graft

Maxillary model showing 3 implants secured by bone grafting material.

ORDER CODE SBG-5 - \$49.99

## Study Models



### SCV-32H

Hinged upper and lower model in crystal clear with 32 natural shaded teeth and anatomical roots.

ORDER CODE SCV-32H - \$89.99



### STW-32H

Hinged upper and lower model, with 32 teeth in white.

ORDER CODE STW-32H - \$49.99



### STP-32H

Hinged upper and lower model, in pink, with 32 natural shaded teeth.

ORDER CODE STP-32H - \$79.99

## Dentures



### SED-200

Pink upper arch on black base.

ORDER CODE SED-200 - \$19.99

### SED-400

Pink lower arch on black base.

ORDER CODE SED-400 - \$19.99

### SLD

Full lower denture in pink acrylic.

ORDER CODE SLD - \$42.99

### SUD

Full upper denture in pink acrylic.

ORDER CODE SUD - \$42.99



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.



# Sutures



**Absorbable:** Plain Gut, Chromic Gut, Polysin, Polysin FA, Monoderm  
**Non-Absorbable:** Polyviolene, Nylon, Silk Polypropylene

- Similar to Vicryl™ but greater tensile strength
- Unique polyglide coating for smooth, NO DRAG Suture passage through tissue

## PGA Polyglycolic Acid Synthetic Absorbable - Braided

Resorbs in 21-28 Days

|  | Needle      | Size | Length     | Color  | Vicryl™ Comparison |               | Per Box         | 10 Or More      |
|--|-------------|------|------------|--------|--------------------|---------------|-----------------|-----------------|
|  | C6 (FS-2)   | 3-0  | 27" (68cm) | Undyed | J423G              | <b>M-423B</b> | \$ 52.50 Box 12 | \$ 49.90 Box 12 |
|  | C6 (FS-2)   | 4-0  | 27" (68cm) | Undyed | J422G              | <b>M-422B</b> | \$ 52.50 Box 12 | \$ 49.90 Box 12 |
|  | C3 (P-3)    | 5-0  | 18" (45cm) | Undyed | J493G              | <b>M-493B</b> | \$ 94.50 Box 12 | \$ 89.90 Box 12 |
|  | C3 (P-3)    | 6-0  | 18" (45cm) | Undyed | J492G              | <b>M-492B</b> | \$ 94.50 Box 12 | \$ 89.90 Box 12 |
|  | PC31 (PS-2) | 3-0  | 18" (45cm) | Undyed | J497G              | <b>M-497B</b> | \$ 94.50 Box 12 | \$ 89.90 Box 12 |
|  | PC31 (PS-2) | 4-0  | 18" (45cm) | Undyed | J496G              | <b>M-496B</b> | \$ 94.50 Box 12 | \$ 89.90 Box 12 |
|  | PC31 (PS-2) | 5-0  | 18" (45cm) | Undyed | J495G              | <b>M-495B</b> | \$ 94.50 Box 12 | \$ 89.90 Box 12 |
|  | T14 (SH-1)  | 3-0  | 27" (68cm) | Green  | J311G              | <b>M-311B</b> | \$ 46.90 Box 12 | \$ 44.50 Box 12 |
|  | T14 (SH-1)  | 4-0  | 27" (68cm) | Green  | J310G              | <b>M-310B</b> | \$ 46.90 Box 12 | \$ 44.50 Box 12 |

## Gut Plain

Resorbs in 5-10 Days

|  | Needle     | Size | Length     | Ethicon™ Comparison |                | Per Box          | 10 Or More       |
|--|------------|------|------------|---------------------|----------------|------------------|------------------|
|  | C6 (FS-2)  | 3-0  | 27" (68cm) | H822G               | <b>554B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C6 (FS-2)  | 4-0  | 27" (68cm) | H822G               | <b>551B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C6 (FS-2)  | 4-0  | 18" (45cm) |                     | <b>591B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C6 (FS-2)  | 3-0  | 18" (45cm) |                     | <b>592B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C3 (P-3)   | 5-0  | 18" (45cm) | 686G                | <b>596B</b>    | \$ 65.90 Box 12  | \$ 62.50 Box 12  |
|  | C31 (X-1)  | 3-0  | 18" (45cm) | 612G                | <b>594B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | CP1 (PC-1) | 6-0  | 18" (45cm) | 1916G               | <b>M-1415B</b> | \$ 114.50 Box 12 | \$ 107.90 Box 12 |

## Chromic Gut

Resorbs in 7-14 Days

|  | Needle     | Size | Length     | Ethicon™ Comparison |                | Per Box          | 10 Or More       |
|--|------------|------|------------|---------------------|----------------|------------------|------------------|
|  | C6 (FS-2)  | 3-0  | 27" (68cm) | 636G                | <b>560B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C6 (FS-2)  | 4-0  | 27" (68cm) | 635G                | <b>559B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C6 (FS-2)  | 3-0  | 18" (45cm) |                     | <b>525B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C6 (FS-2)  | 4-0  | 18" (45cm) |                     | <b>558B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |
|  | C3 (P-3)   | 4-0  | 18" (45cm) | 1654G               | <b>595B</b>    | \$ 59.90 Box 12  | \$ 56.90 Box 12  |
|  | C3 (P-3)   | 5-0  | 18" (45cm) | 687G                | <b>1248B</b>   | \$ 59.90 Box 12  | \$ 56.90 Box 12  |
|  | C1 (P1)    | 6-0  | 12" (30cm) |                     | <b>M-1242B</b> | \$ 114.50 Box 12 | \$ 107.90 Box 12 |
|  | T4 (SH-1)  | 4-0  | 14" (35cm) |                     | <b>589B</b>    | \$ 41.90 Box 12  | \$ 39.90 Box 12  |
|  | T28 (RB-1) | 5-0  | 14" (35cm) |                     | <b>590B</b>    | \$ 32.90 Box 12  | \$ 31.25 Box 12  |





# Sutures

## Suture Needle Guide

  
Reverse  
Cutting

  
Precision  
Reverse  
Cutting

  
Taper  
Point



  
Taper  
Cutting




## Silk

|   | Needle                              | Size              | Length                                 | Ethicon™ Comparison  |   | Per Box   | 10 Or More  |
|---|-------------------------------------|-------------------|--|----------------------|---|---|---|
|    | C8 (FS)                             | 2-0               | 18" (45cm)                             | 685-G                | <b>770B</b>                               | \$ 32.90 Box 12                                       | \$ 31.25 Box 12                                       |
|    | C6 (FS-2)<br>C6 (FS-2)<br>C6 (FS-2) | 3-0<br>4-0<br>5-0 | 18" (45cm)<br>18" (45cm)<br>18" (45cm) | 7772<br>683G<br>682G | <b>784B</b><br><b>781B</b><br><b>754B</b> | \$ 32.90 Box 12<br>\$ 32.90 Box 12<br>\$ 32.90 Box 12 | \$ 31.25 Box 12<br>\$ 31.25 Box 12<br>\$ 31.25 Box 12 |
|    | C7 (FS-1)                           | 3-0               | 18" (45cm)                             | 684G                 | <b>786B</b>                               | \$ 32.90 Box 12                                       | \$ 31.25 Box 12                                       |
|    | C3 (P-3)                            | 5-0               | 18" (45cm)                             | 640G                 | <b>M-1265B</b>                            | \$ 89.50 Box 12                                       | \$ 84.50 Box 12                                       |
|   | C17 (C-3)<br>C17 (C-3)              | 4-0<br>6-0        | 18" (45cm)<br>18" (45cm)               | 736G<br>711G         | <b>787B</b><br><b>752B</b>                | \$ 32.90 Box 12<br>\$ 32.90 Box 12                    | \$ 31.25 Box 12<br>\$ 31.25 Box 12                    |
|  | C31 (X-1)                           | 3-0               | 18" (45cm)                             | 632G                 | <b>785B</b>                               | \$ 32.90 Box 12                                       | \$ 31.25 Box 12                                       |
|  | C1 (P1)                             | 6-0               | 18" (45cm)                             | 639G                 | <b>M-1260B</b>                            | \$ 89.50 Box 12                                       | \$ 84.50 Box 12                                       |

## Polypropylene

|   | Needle               | Size       | Length                   | Ethicon™ Comparison |                                  | Per Box                            | 10 Or More                         |
|---|----------------------|------------|--------------------------|---------------------|----------------------------------|------------------------------------|------------------------------------|
|  | C3 (P-3)<br>C3 (P-3) | 5-0<br>6-0 | 18" (45cm)<br>18" (45cm) | 8698G<br>8695G      | <b>M-1085B</b><br><b>M-1082B</b> | \$ 99.50 Box 12<br>\$ 99.50 Box 12 | \$ 94.90 Box 12<br>\$ 94.90 Box 12 |
|  | C1(P-1)<br>C1(P-1)   | 6-0<br>7-0 | 18" (45cm)<br>18" (45cm) | 8697G<br>8696G      | <b>M-1079B</b><br><b>2468N</b>   | \$ 99.50 Box 12<br>\$ 89.50 Box 12 | \$ 94.90 Box 12<br>\$ 84.50 Box 12 |

## Nylon Black Monofilament

|   | Needle                 | Size       | Length                   | Ethicon™ Comparison |                                | Per Box                            | 10 Or More                         |
|---|------------------------|------------|--------------------------|---------------------|--------------------------------|------------------------------------|------------------------------------|
|  | C6 (FS-2)<br>C6 (FS-2) | 4-0<br>5-0 | 18" (45cm)<br>18" (45cm) | 662G<br>661G        | <b>M-922B</b><br><b>M-918B</b> | \$ 41.90 Box 12<br>\$ 41.90 Box 12 | \$ 39.90 Box 12<br>\$ 39.90 Box 12 |

## Polyviolene-Green Braided Polyester for Minimized Tissue Trauma

|   | Needle    | Size | Length     | Ethicon™ Comparison |              | Per Box         | 10 Or More      |
|---|-----------|------|------------|---------------------|--------------|-----------------|-----------------|
|  | C6 (FS-2) | 4-0  | 18" (45cm) | R683G               | <b>1180B</b> | \$ 34.50 Box 12 | \$ 32.50 Box 12 |
|  | Y-5 (V-5) | 5-0  | 18" (45cm) |                     | <b>797B</b>  | \$ 65.90 Box 12 | \$ 62.50 Box 12 |



# Scalpel Blades

## Glass Van Surgical Blades



Blade No. 11, 12, 12B, 15 & 15C fit handle numbers 3, 3L, 5, 7 & 9

- Supplied sterile, in boxes of 100 pcs
- Individually V.C.I paper wrapped & foil pack
- Available in both stainless steel & carbon steel

Blade No. 20, 21, 22, 23 & 25 fit handle numbers 4, 4L & 6

- Supplied sterile, in boxes of 100 pcs
- Individually V.C.I paper wrapped & foil pack
- Available in both stainless steel & carbon steel

**ALL BLADES**

1-5 boxes - **\$24.99**

6-10 boxes - **\$21.99**

11 or more boxes - **\$19.99**

**Also available in all the above sizes:**

**Sterile, disposable, plastic handled scalpels**

1-5 boxes - **\$8.79**

6 or more boxes - **\$7.79**



Blade No. 11



Blade No. 12



Blade No. 12B



Blade No. 15



Blade No. 15C



Blade No. 20



Blade No. 21



Blade No. 22



Blade No. 23



Blade No. 25



# Instruments

## Scalers and Curettes - \$25.88

We carry a full line of all the listed scalars and curettes. We will also customize to fit your exact needs.



### Universal Curettes

- Barnhart
- Columbia
- Goldman Fox
- Langer
- Gracey Universal
- McCall
- Younger Good
- Gracey

### Sickle Scalers


- H Series
- Goldman Sickles
- McCall Sickle
- Morse Sickle
- MTJZ Sickle
- N&Z Style Sickles
- 204 Series Sickles
- H5/Jacquette Sickles
- Jacquette Sickles
- H5 Sickle/Hoe Combinations

Sabra carries a full line of hirschfeld files, perio surgical curettes, and perio surgical chisels, single and double ended explorers, probes, dressing pliers, mouth mirrors, cavity prep instruments, hatchets, excavators, amalgam carriers, crown removers, calipers, gauges, crown scissors, burnishers, pluggers, condensers, plastic filling and composite instruments, carving instruments, cord packers, waxing instruments, root and silver point forceps, rubber dam instruments, scalpel handles and tappers. Please call or check at [www.sabradent.com](http://www.sabradent.com) for details.

## Composite Instruments (Titanium Coated) - \$37.00



## Elevators - \$62.00



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
 ©2013 Sabra Dental Products, All Rights Reserved.



# Instruments



## Crown Removers



SE-1 - \$93.13  
Crown splitter

DECR-1 - \$49.38  
Crown Remover

DECR-2 - \$51.25  
Crown Remover

## Forceps



K-1192 - \$135.00  
150 Cryer

K1198- \$135.00  
151 Cryer

K1167- \$135.00  
23 Cowhorn

K-2162- \$144.50  
Root Forcep

## Ronguers



K-45967 - \$199.50  
Micro Friedman

K-4111 - \$199.50  
Friedman Ronguer

K4166- \$199.50  
Blue Menthol Ronguer

## Scissors



k-4019 - \$51.65  
Iris

k-4017 - \$66.25  
GF

k-4001 - \$64.25  
Kelly

## Periotome



1 - \$54.75

2 - \$54.75

3 - \$54.75

4 - \$54.75

5 - \$54.75

FRP3 - \$54.75

FRP4 - \$54.75

Kurtz FRP1- \$55.55

KurtzFRP2 - \$55.55

## Periosteals



9 Molt - \$50.25

9L - \$50.25

24Glick - \$50.25

Wood WDS 1/h2 - \$ 50.25

PR- 3 - \$50.25



# Surgical Drape Packs / Small Set Up Packs

## Set Up Pack 1



**ORDER CODE 40566 SSI - \$149.25**

**(packed 15 per case)**

- 1 - Gown, large non-RF
- 2 - Hand Towels, ABS White
- 1 - CSR wrap, 24" x 24"
- 1 - Isolation Gown
- 1 - Bib 3 ply, blue

## Set Up Pack II



**ORDER CODE 40567SSII - \$209.25**

**(packed 15 per case)**

- 1 - CSR wrap 30" x 30"
- 1 - Bib 3 ply, Blue
- 1 - Light handle covers 4" x 6"
- 20 - Sponges 2" x 2"
- 2 - Sabra Sleeves, hose covers
- 1 - White saliva ejector
- 1 - HVE tip
- 2 - Utility drape 14" x 25" w/ tape
- 2 - Cotton tip applicators 6"
- 2 - Hand towels, ABS White
- 1 - CSR wrap, 24" x 24"
- 1 - Gown, large, non-RF
- 1 - Isolation Gown (overwrap)

## Set Up Pack III



**ORDER CODE 40568SSIII - \$224.25**

**(packed 15 per case)**

- 1 - Mayo stand cover, 23" x 55"
- 1 - Suture bag
- 2 - Utility drape 14" x 25" w/ tape
- 1 - Half Drape 60" x 44"
- 1 - Table Cover 44" x 90"
- 1 - CSR wrap 24" x 24"
- 2 - Hand towels, ABS white
- 1 - Gown, Large, non-RF
- 1 - Isolation Gown (overwrap)

## Essential Pack



**ORDER CODE 40563EP - \$159.80**

**(packed 4 per case)**

- 1 - Table cover 44" x 90"
- 1 - Drape, Split sheet 55" x 75"
- 3 - Sabra Sleeve hose covers
- 3 - Blue cloth O.R. towels
- 2 - Gowns, Sterile, back large
- 2 - Hand towels, ABS white
- 1 - Mayo stand cover 23" x 55"

## Split Drape Basic Implant Pack-2 packs in one

**ORDER CODE 40809SIP - \$195.80**

**(packed 4 per case)**

### Surgery Pack

- 1 - CSR wrap 45" x 45"
- 4 - Light handle covers 4" x 6"
- 1 - Mayo stand cover 23" x 55"
- 1 - Blue cloth towel
- 1 - Suture bag
- 4 - Utility drape 14" x 25" w/ tape
- 25 - Sponges 2" x 2"
- 8 - Sponges 4" x 4"
- 1 - Bouffant cap
- 2 - Cotton tip applicators 6"
- 4 - Sabra Sleeves hose cover
- 1 - Drape, split sheet 55" x 75"
- 1 - Surg-OVac surgical tip
- 1 - Saliva ejector
- 1 - HVE tip
- 1 - Headrest cover 10" x 13"
- 1 - Machine cover 21" x 36"

**ORDER CODE 40809SIP - \$195.80**

**(packed 4 per case)**

### 2 packs in one

- 1 - CSR wrap 24" x 24"
- 2 - Face masks with ear loops
- 2 - Gown, large non-RF
- 2 - Hand towels, ABS white
- 2 - Bouffant caps



**Convenient and Controls Cost • Streamlines Inventory – Turns 10 Inventory Items into 1  
Saves Set Up and Ordering Time • Lowest Prices - Highest Quality – Smallest Case Quantities**

**Call To Make Your Own Custom Pack**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Surgical Drape Packs / Small Set Up Packs



## Head and Neck Pack



**ORDER CODE 40562HNP - \$154.00**  
(packed 4 per case)

- 1 - Table cover 44" x 90"
- 1 - Drape, split sheet 55" x 75"
- 1 - Suture bag
- 1 - Drape, head w/ tape
- 2 - Gowns, large non-RF
- 2 - Hand towels, ABS white
- 1 - Mayo stand cover 23" x 55"

## Standard Drape Pack



**ORDER CODE 40560SDP - \$159.80**  
(packed 4 per case)

- 1 - Tray, large foam 14" x 18"
- 1 - CSR wrap 45" x 45"
- 1 - Drape, split sheet 55" x 75"
- 2 - Light handle covers 4" x 6"
- 2 - Sheet foil 9" x 10.75"
- 2 - Sabra sleeves hose covers
- 1 - Cap, bouffant 24" blue
- 2 - Gowns, large non-RF
- 2 - Hand towels, ABS white

# Sterile Separates (single wrapped)



### Sterile Surgical Gown / Towel

**ORDER CODE SS79878 – Large, case of 20 - \$185.80**

**ORDER CODE SS79888 – X Large, case of 20 - \$187.50**

Blue O.R. Fluid resistant/ White cuff, Velcro Closure on back collar, and Tie on waist closure



### Isolation Gown / Towel

**Sterile**

**ORDER CODE SS6010B – Uni-size, case of 20 - \$195.80**

Blue O.R. Fluid resistant/ White cuff, Velcro Closure on back collar, and Tie on waist closure

**Non-Sterile (not single wrapped)**

**ORDER CODE 6010B – Uni-size , case of 50 - \$63.99**

Blue O.R. Fluid resistant/ White cuff, Velcro Closure on back collar, and Tie on waist closure

**ORDER CODE 6010P – Uni-size, case of 50 - \$63.99**

Blue O.R. Fluid resistant/ White cuff, Velcro Closure on back collar, and Tie on waist closure



### Sterile Patient Drapes

**ORDER CODE SS49128, case of 20 - \$198.50**

Fenestrated drape 38" x 55"

**ORDER CODE SS170946, case of 20 - \$199.60**

EENT Split sheet 55" x 75"


2.5" x 27" adhesive split Absorbent pad w/ tube holder



### Sterile Table Cover

**ORDER CODE SS170864, case of 20 - \$69.90**

Table cover 44" x 90"



### Sterile Mayo Stand Cover

**ORDER CODE SS170915, case of 20 - \$87.80**

Mayo Stand cover 23" x 55"



### Sterile Machine Cover

**ORDER CODE SS107787, case of 20 - \$54.50**

Machine cover 21" x 36"



### Sterile Sabra Sleeves

**ORDER CODE SS9335, case of 25 - \$69.00**

Sabra Sleeve hose covers 2.5" x 66"

**ORDER CODE 9935, case of 25- Non Sterile - \$39.99**

Sabra Sleeves hose covers 2.5" x 66"

**Call To Make Your Own Custom Pack**



Call Toll Free at 1-800-888-4435

# Face Masks & Shields



## **DEFEND+PLUS Pleated Procedural Ear-Loop Mask**

The latest addition to our mask lineup, the DEFEND+PLUS mask meets all criteria for ASTM standards for Moderate Barrier protection. Perfect for dust, airborne particles, moisture, spray and splatter. The same great features of our Diffuser Pleated Mask with increased protection (BFE and PFE). Provides moderate barrier filtration.

**ORDER CODE MK-4046 Blue**  
50 per box - \$36.60



## **"Diffuser" (Anti-Fog) Pleated Ear-Loop & Tie-On Mask**

Our most popular mask contains an antifog foam nosepiece that absorbs moisture to prevent eyewear from fogging. Constructed of a soft, 3-ply, pleated material for ultimate comfort. The fluid-resistant masks have a BFE and PFE of greater than 99%.

**ORDER CODE MK-1046 Ear - Loop Blue**  
**MK-1056 Ear - Loop Pink**  
**MK-1076 Tie - On Blue**  
50 per box - \$10.45



## **Economy Procedure Pleated Ear-Loop Mask**

The same mask as our popular "Diffuser" earloop mask without the antifog nosepiece. This economy mask is made of a lightweight, 3-ply material that is fluid resistant with a BFE and PFE of greater than 99%. Ultrasonically sealed earloops offer great comfort and durability.

**ORDER CODE MK-2007 Blue**  
**MK-2008 Pink**  
**(10% smaller)**  
**MK-2010 White**  
50 per box - \$8.25

**ORDER CODE MK-2000 Blue**  
**MK-2001 Pink**  
100 per box - \$14.40



## **Breathe E-Z Pleated Ear-Loop Mask**

These fluid-resistant masks contain a soft, 3-ply, pleated material for ultimate "cool-breathing" comfort. This mask is made with a fiberglass-free material and a latex-free elastic ear-loop to provide maximum comfort. Our easiest mask to breathe through yet they have a BFE and PFE of greater than 99%. Provides primary barrier filtration.

**ORDER CODE MK-1246 Blue**  
50 per box - \$4.95



## **Premium Molded Mask**

The fluid resistant masks have a BFE and PFE of greater than 99%. These masks are available in either regular size or a smaller size in pink - both with a unique "non-glare", adjustable nosepiece for a comfortable fit.

**ORDER CODE MK-1006 Blue**  
**MK-1016 Pink**  
**(10% smaller)**  
50 per box - \$9.80



## **Premium Ear-Loop Cone Masks**

Innovative earloop design cone mask allows the user to wear the mask with loops around the ears instead of the back of the head. These masks have a non-latex elastic band and are designed for easy breathing. They are highly pliable, fluid resistant, and feature a non-glare nosepiece. 99% + PFE & 99% + BFE.

**ORDER CODE MK-1036 Blue**  
**MK-1066 Pink**  
**(10% smaller)**  
50 per box - \$10.25



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Face Masks & Shields



**Valuline Molded Cone Mask**  
 These comfortable molded masks are the easiest to breathe through. Generally used for examinations and any non-splatter, dustless, or fluidless procedures. Offers superior protection and a BFE/PFE of greater than 99%.  
**ORDER CODE MK-1106 Blue**  
 50 per box - \$6.15



**Pleated Ear-Loop Mask w/Shield**  
 Practical, comfortable pleated ear-loop mask that has its own face shield attached. The clear, optically correct face shield provides excellent visibility while protecting you against splashing. It is both anti-fog treated and fluid resistant. It's 99% + PFE and BFE rating make this mask reliable and dependable, while the fiberglass and latex free material make this mask hyperallergenic. Provides primary barrier filtration.  
**ORDER CODE MK-1010 Blue**  
 100 per box - \$143.15  
**ORDER CODE MK-1011 Blue**  
 50 per box (economy) - \$51.45



**Face Shields**  
 This Face Shield is made to fit over glasses and face masks. The adjustable velcro headband allows for a comfortable fit. The shield protects against splashing. It is anti-fog treated and latex and fiberglass free.

**ORDER CODE MK-3000 Full Shield**  
 25 per box - \$69.45  
**ORDER CODE MK-3001 Half Shield**  
 25 per box - \$65.85



**Economy Molded Mask**  
 This quality mask offers a pliable nose piece with a single band for a cool and comfortable fit, while offering superior protection. This lightweight material is also fiberglass-free and the band is latex-free.

**ORDER CODE MK-1007 Blue, MK-1008 Pink**  
 (10% smaller)  
 50 per box - \$8.80  
**ORDER CODE MK-1000 Blue, MK-1001 Pink**  
 (10% smaller)  
 100 per box - \$15.40



**N95 Particulate Respirator**  
 N95-level respirator is at least 95% efficient when tested against sodium chloride (salt). Used for solid and non-oil based particles. Applications include grinding, sanding, bagging and general processing of various minerals and other substances that do not contain oil. NIOSH approved. Provides N95 maximum barrier filtration.  
**ORDER CODE MK-N95**  
 20 per box - \$18.80



**Form-Fit Mask**  
 This unique contour design without a metal nose piece prevents fogging of eyewear. These masks are ergonomically designed, fluid resistant and are latex and fiberglass free. BFE of greater than 99%.  
**ORDER CODE MK-1206 Blue**  
 50 per box - \$11.60  
**ORDER CODE MK-1216 Pink**  
 50 per box - \$11.60



# Disinfectants / Ultrasonic Cleansers / Soaps



## DEFEND + PLUS Disinfecting/Deodorizing/Cleaning Wipes

These disinfecting wipes, with new 2 minute kill time\* contain a stable, low pH formulated disinfectant and deodorant for use on hard, non-porous surfaces, in hospitals, emergency medical settings, laboratories, medical offices, dental offices – where control of cross contamination is required. Containing a dual chain quaternary ammonium/alcohol formula, the product is proven effective against the following pathogens in 2 minutes\* at 20 degrees Celsius/60 degrees Fahrenheit, on hard, non-porous surfaces when used as directed.

### \*2 minute kill time:

- Escherichia coli O157:H7
- Human Hepatitis B virus
- Human Hepatitis C virus
- Herpes Simplex virus 2 (genital herpes virus)
- HIV-1 (AIDS virus) – 1 minute kill.
- Influenza A (A2 Japan)
- Methicillin Resistant Staphylococcus aureus (MRSA)

- Mycobacterium tuberculosis
- Pseudomonas aeruginosa
- Salmonella enterica
- Staphylococcus aureus
- Vancomycin Resistant Enterococcus faecalis (VRE)

**See product label for complete information.**

### 5 minute kill time:

- Trichophyton mentagrophytes (Athlete's foot fungus)

**ORDER CODE SO-9000**  
**160 Per Canister - \$9.25**  
**6" x 6" Sheets**



## Omni-Care7™ Medicated Lotion Soap

The gentle and non-drying formulation found in this soap with Biostructure™ contains 7 natural ingredients including aloe and vitamin E. The non-toxic and non-irritating soap also contains the broad-spectrum antimicrobial agent PCMX and has a pleasant baby powder scent.

**ORDER CODE SO-9300**  
**1 Gallon Bottle - \$20.55**  
**ORDER CODE SO-9310**  
**16 oz. Bottle with Pump - \$9.25**



## Tartar & Stain Remover #4 Ultrasonic Solution

Use this potent stain remover on dentures, bridges, crowns and orthodontic appliances. The ready to use formula quickly breaks up and removes calculus and tartar while eliminating stains caused by tobacco and food. Use in a beaker placed inside the ultrasonic unit for best results.

**ORDER CODE SO-9600**  
**1 Gallon Bottle - \$20.55**



## General Purpose Cleaner #1 Ultrasonic Solution

This nonammoniated, concentrated formula, effectively cleans dried blood, tissue, pumice, food, rouge, Tripoli & buffing compounds from the surface of any appliance while providing a rust inhibitor.

**ORDER CODE SO-9400**  
**1 Gallon Bottle - \$20.55**



## Temporary Cement Remover #6 Ultrasonic Solution

This ready to use solution requires no mixing. The strong, fast-acting formula removes cement and stone quickly and completely from instruments, trays and other appliances.

**ORDER CODE SO-9800**  
**1 Gallon Bottle - \$20.55**



## PS-10 Pre-Soak

Used as a pre-soak cleans instruments prior to sterilization.

**ORDER CODE SO-9200**  
**1 Gallon Bottle - \$30.85**



## SRG Evacuation System Cleaner

Used on a daily basis this superior evacuation system cleaner will keep lines running free and clear at maximum suction. Clinically proven to clean lines and discourage build-up. The formula is non-foaming, biodegradable, includes a germicide and has a neutral pH while maintaining the ability to dissolve organic matter and eliminate odors. Has a fresh wintergreen fragrance. 128 treatments per gallon.

- Approved for use with amalgam separators
- Daily use keeps lines clean and at maximum suction
- Eliminates odors
- Neutral pH

**ORDER CODE SO-9100**  
**1 Gallon Bottle with Pump - \$40.70**



## Ultrasonic Enzymatic Tablets

These easy-to-use enzymatic tablets dissolve in a gallon of water in under a minute to create a powerful and effective ultrasonic cleaning solution. Made from a biodegradable, non-chlorine, non-corrosive formulation that contains fast acting protease enzymes. The solution will greatly assist in the removal of blood, bone, tissue, food particles and other soils that may be resistant to chemical detergents. Virtually eliminates handling and storage of heavy liquid concentrates. Contains a pleasant fragrance and is a neutral pH. May also be used as an instrument pre-soak or as an evacuation system cleaner.

- Biodegradable, neutral pH
- Non-chlorine
- Non-corrosive

**ORDER CODE UT-1000**  
**64 Tablets per Box - \$39.95**

\*Universal Diameter.



# Sterilization



### DEFEND+PLUS Pouches – with Built-in Internal Indicator

DEFEND+PLUS pouches are constructed with a built-in, lead-free Internal Indicator, giving you peace of mind when it comes to the sterilization of your instruments. The Internal Indicator ensures the correct sterilization temperature is reached where you really need it...inside the instrument compartment. This eliminates the need for an internal indicator strip inside the pouch with your instruments.

Awarded “Best Product 2010” by Dental Product Shopper and “Editors’ Choice” by the Dental Advisor.

- ORDER CODE SP-0550 (2.25" x 4") - \$6.95
- SP-0550 (2.75" x 10") - \$8.80
- SP-1500 (3.5" x 5.25") - \$8.15
- SP-2500 (3.5" x 10") - \$8.80
- SP-3500 (5.25" x 10") - \$15.90
- SP-4500 (7.5" x 13") - \$33.80
- SP-5500 (9" x 15") - \$39.90
- SP-8000 (13" x 20") - \$63.95
- 200 Per Box



### Sterilization Wrap (CSR)

DEFEND Sterilization Wrap (CSR) is perfect for general purpose sterilization, as it provides a barrier against air and water-borne bacteria. This CSR wrap conforms easily around various shapes and sizes of cassettes, allowing excellent penetration of various methods of sterilization. Compatible with steam, ETO and gamma. Petroleum- and memory-free material with a cloth-like feel.

- ORDER CODE SW-1200 (12" x 12") - \$83.25
- SW-1500 (15" x 15") - \$117.40
- 1000 Per Box
- ORDER CODE SW-2000 (20" x 20") - \$99.90
- SW-2400 (24" x 24") - \$133.85
- 500 Per Box
- ORDER CODE SW-3000 (30" x 30") - \$108.10
- 250 Per Box



### Autoclave Tape – Sterilization Indicator Tape

This tape is for use in autoclave, dry heat or chemical sterilizers and will seal autoclave bags, tubing, pouches and CSR wraps. The tape will indicate when sterile conditions have been attained. Extra tacky tape to fix on all types of paper and plastic.

- ORDER CODE # AT-2001 (1/2") - \$3.80
- AT-2002 (3/4") - \$4.95
- AT-2003 (1") - \$6.15
- 60 Yards Per Roll

**Pro-Sterile**

Reusable & See through container

Can be steam sterilized to 135 C

Environmental Friendly

User friendly

last forever

Cost Effective

### Pro-Sterile

The Pro-Sterile is the newest in Green efforts. This is a re-usable autoclave tube. Simply insert a handpiece or an instrument, insert an autoclave marker close and autoclave. Let cool like normal then re-use again and again - save money - save the planet

- ORDER CODE # SPT-01
- \$7.95 each
- \$6.95 each for 10 or more
- \$5.95 each for 30 or more
- \$4.95 each for 60 or more



# Barrier Products




## Digital X-Ray Sensor Sleeves

All the popular sensor sleeves are now available. Protect your investment in digital imaging with Defend sleeves which boast leak proof construction with extra soft, tear resistant plastic.

### ORDER CODE

- BF-8000 (1 5/8" x 8")
- Fits Trophy, Gendex & Diagnodent - \$23.85
- BF-8100 (2 1/2" x 10")
- Fits Various Systems - \$9.25
- BF-8500 (1 3/8" x 8")
- Fits Trophy & Gendex Systems - \$7.20
- 500 Box



(Barrier Film Not Included)

## Barrier Film Dispenser

Reduce & Reuse for a Better Work Environment! Made of high-impact material, the DEFEND Barrier Film Dispenser helps organize infection control supplies. Great for DEFEND Barrier Film (BF-1600, BF-2600), now available in environmentally-friendly condensed packaging Reusable, reduces waste of cardboard dispenser boxes. Minimizes office clutter. Usable with all DEFEND Barrier Film (BF-1600, #BF-2600, BF-1500 or BF-2500).

### ORDER CODE

- BF-1800 - \$23.85



## Handpiece Barrier Sleeves

These Handpiece Barrier Sleeves are made to fit over most handpiece instruments. Made of high quality clear plastic, the Handpiece Barrier Sleeves are simple to use and easy to dispose.

### ORDER CODE

- IN-7006 (1.5" x 8") - \$10.25
- 500 Per Box



## DEFEND Chair Sleeves

DEFEND Chair Sleeves are made of a clear plastic and are designed to easily fit over most chairs. These disposable chair sleeves are highly durable and easy to use.

### ORDER CODE

- BF-4000 (27.5" X 24") - 1/2 Chair - \$26.30
- BF-4500 (29" x 79") - Full Chair - \$35.55
- 250 Per Box



## Headrest Covers

These Headrest Covers are designed to fit over the top of headrests. They are made of clear plastic and are easy to apply and dispose of.

### ORDER CODE

- BF-9000 (11" x 9.5" x 2") - \$7.20
- BF-9500 (14" x 9.5" x 2") - \$8.20
- 250 Box



## Syringe Sleeve Covers

Fits most 3-way air/water syringes and have a pre-punched opening making application easy. They can also be used to cover handpieces, saliva ejectors, ultrasonic scalers and more.

### ORDER CODE

- BF-3000 (2.5" x 10") - \$8.80
- 500 Per Box



## X-Ray Head Sleeves

These X-Ray Head Sleeves are designed to fit over most X-Ray heads. The barrier is made of clear plastic and is easy to apply.

### ORDER CODE

- BF-6000 (15" x 26") - \$22.15
- 250 Per Box



## Light Handle Sleeves

These Light Handle Sleeves are made of clear plastic and are easy to apply. Made to fit most light handles.

### ORDER CODE

- BF-5000 (4" x 5.75") - \$8.70
- 500 Per Box



## Tray Sleeves

These Tray Sleeves cover the entire instrument tray and are made of clear plastic which is easy to apply and easier to dispose of.

### ORDER CODE

- BF-7000 (10.5" x 14") - \$8.80
- BF-7500 (11 5/8" x 16") - \$23.20
- 500 Per Box





# Gloves



### Blackjack Powder-Free Textured Latex Gloves

Our Blackjack Powder-Free Textured Latex Gloves are polymer coated for an easy on / off. The black color and textured finish for enhanced grip make the Blackjack gloves great for the tattoo professional.

**ORDER CODE**

- LG-8002 - X-Small (6) - \$9.50
- LG-8003 - Small (6 1/2) - \$9.50
- LG-8004 - Medium (7 1/2) - \$9.50
- LG-8005 - Large (8 1/2) - \$9.50
- 100 Per Box
- LG-8006 - X-Large (9) - \$9.50
- 90 Per Box



### Powdered Vinyl Exam Gloves

Our Powdered Vinyl Exam Gloves are lightly powdered and are latex free. They are made of high-quality stretch vinyl which are tear resistant and have an extra soft feel. Vinyl gloves are great alternative for individuals who may be sensitive to latex.

**ORDER CODE**

- VG-3003 - X-Small (6 1/2) - \$6.15
- VG-3004 - Medium (7 1/2) - \$6.15
- VG-3005 - Large (8 1/2) - \$6.15
- VG-3006 - X Large (9) - \$6.15
- 100 Per Box



### Powder-Free Nitrile Textured Exam Gloves

Nitrile Protection Meets Superior Comfort  
Our Powder Free Nitrile Textured Exam Gloves are lightweight, yet highly puncture resistant, have great tactile strength and dexterity and are a great alternative for individuals sensitive to latex materials.

**ORDER CODE**

- NG-2003 - Small (6 1/2) - \$8.50
- NG-2004 - Medium (7 1/2) - \$8.50
- NG-2005 - Large (8 1/2) - \$8.50
- NG-2006 - X Large (9) - \$8.50
- 100 Per Box



### Powder-Free Micro Textured Latex Exam Gloves

Our Powder-Free Latex Exam Gloves are thin for extra tactile sensitivity and have a micro-textured finish. The gloves are white and Polymer coated for easy on/off comfort.

**ORDER CODE**

- LG-2102 - X-Small (6) - \$9.30
- LG-2103 - Small (6 1/2) - \$9.30
- LG-2104 - Medium (7) - \$9.30
- LG-2105 - Large (8 1/2) - \$9.30
- LG-2106 - X-Large (9) - \$9.30
- 100 Per Box



### Powdered Professional Latex Exam Gloves

Our lightly powdered professional latex exam gloves are soft, strong, low modulus gloves which are low protein and contain minimal powder. The boxes are color coded according to size.

**ORDER CODE**

- LG-1102 - X-Small (6) - \$8.05
- LG-1103 - Small (6 1/2) - \$8.05
- LG-1104 - Medium (7 1/2) - \$8.05
- LG-1105 - Large (8 1/2) - \$8.05
- 100 Per Box



### Powder-Free Textured Nitrile Exam Gloves (LITE)

Our Powder-Free Textured Nitrile Exam Gloves (LITE), purple in color, are packed 200 per box to provide maximum value. These latex free gloves are highly puncture resistant and provide great tactile strength and dexterity.

**ORDER CODE**

- NG-3002 - X-Small (6) - \$6.15
- NG-3003 - Small (6 1/2) - \$6.15
- NG-3004 - Medium (7 1/2) - \$6.15
- NG-3005 - Large (8 1/2) - \$6.15
- NG-3006 - X-Large (9) - \$6.15
- 100 Per Box



# Apparel & Eyewear



## DEFEND+PLUS Latex-Free Disposable Jackets

DEFEND+PLUS latex-free disposable jackets provide comfort and style! Made of a soft, breathable fabric, these fluid-resistant jackets meet OSHA requirements for blood-borne pathogens. The autoclavable, washable DEFEND+PLUS jackets are available in four vibrant colors. Rated 4 1/2+ out of 5 from the Dental Advisor.

These jackets also provide the latest technology in protection, with 28% better filtration (Nelson Laboratories 1/2011), higher quality stitching, cuffs, and collar compared to the leading competition.

- Static Free - No Clinging
- Soft 3-Layer SMS Fabrics
- Fluid Resistant, Breathable
- Snap Front, Knit Cuffs & Collar
- 1 Breast Pocket & 2 Hip Pockets
- Autoclavable and Washable
- Latex Free

### ORDER CODE

|         | Ceil Blue   | Sky Blue    | Light Pink  | Purple      |
|---------|-------------|-------------|-------------|-------------|
| Small   | SJ-3630CBS  | SJ-3630SBS  | SJ-3630LPS  | SJ-3630PPS  |
| Medium  | SJ-3630CBM  | SJ-3630SBM  | SJ-3630LPM  | SJ-3630PPM  |
| Large   | SJ-3630CBL  | SJ-3630SBL  | SJ-3630LPL  | SJ-3630PPL  |
| X-Large | SJ-3630CBXL | SJ-3630SBXL | SJ-3630LPXL | SJ-3630PPXL |

10 Per Bag - \$19.95



## U.V. Protective Eyewear

These economical curing/bonding protective eyewear, filter out over 99% of U.V. and dangerous blue light (up to 540nm). One piece lens construction with protective side shields are durable and able to be worn over most eyeglasses.

ORDER CODE GL-2020

Each - \$20.55



## Clear Protective Eyewear

These economical eyewear offer excellent protection in the operatory or in the lab. The one piece lens is made of a high impact, scratch resistant clear polycarbonate material. They provide non-glare and optically clear protection.

ORDER CODE GL-2019

Each - \$5.10



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

## Apparel



### Disposable Lab Coats

DEFEND has improved the quality of its lab coats to make them more comfortable for dental professionals. The new coats are now more breathable, while retaining the same fluid-resistance and the knit collar and cuffs they have always had. The nonflammable, disposable/reusable coats feature a front pocket and easy-to-use snap closures. The blue coats come designated for both men and women, sized small, medium, large and extra-large.

#### ORDER CODE

|                            |                              |
|----------------------------|------------------------------|
| SG-9003 - Small - \$20.40  | SG-9006 - XLarge - \$20.40   |
| SG-9004 - Medium - \$20.40 | SG-9007 - XXLarge - \$21.15  |
| SG-9005 - Large - \$20.40  | SG-9008 - XXXLarge - \$21.15 |

10 Per Bag



### Disposable Tie-Back Protective Gown

These full length gowns are made from a fluid resistant, cool and comfortable material to provide both protection and comfort. They are an inexpensive alternative to lab coats while having key features such as waist ties and elastic cuffs. The blue gowns are designed for both men and women and come in 2 sizes.

#### ORDER CODE

SG-1000 - Medium (40" length) - \$8.45  
 SG-1001 - Large (48" length) - \$8.45  
 10 Per Bag



### Non-Sterile Bouffant Caps

Non-Woven Polypropylene Bouffant caps are ideal for clean room environments where a light weight head covering is necessary. These 100% polypropylene caps are highly breathable and low linting.

#### ORDER CODE

BC-3501 - 21" (53cm) - \$15.95  
 BC-3503 - 24" (61cm) - \$15.95  
 100/bag



### Non-Sterile Shoe Covers

Elastic pulls shoe cover in around top of shoe for better fit. Made from polypropylene. These are very generously sized shoe covers and will fit most sized shoes and smaller work boots. Opening stretches to approximately 16" heel to toe. Non-Skid Thread, Non-Conductive, One Size Fits All.

#### ORDER CODE

SC-5000 Sc-5000 Shoe Covers - \$15.95  
 100/box



Call Toll Free at 1-800-888-4435



# Disposables



## Cotton Tipped Applicators

Non-sterile, made from all wood dowels with very tight 100% cotton bud to minimize linting and prevent shredding. Comes sealed in autoclavable bag with indicator and are available in either 3" or 6" sizes.

**ORDER CODE**

CS-3000 (3" dowel) - \$5.60

CS-6000 (6" dowel) - \$8.20

1000 Per Box



## Adult Tongue Depressors

These tongue depressors are made of 100% wood and are smooth and splinter-free. They measure 6" in length and are non-sterile.

**ORDER CODE TD-9000 - \$7.20**

500 Box



## Green Seal® Certified Paper Products

Green Seal Certified 100% recycled paper products are the perfect addition to any eco-friendly practice! Certified, recycled white C-Fold Towels meet Green Seal's environmental standard for bleaching, de-inking and packaging. These soft, ultra absorbent towels contain 100% recovered paper fiber and are chlorine-free, with no added pigments, inks, dyes or fragrances. The 100% recycled 2-ply white facial tissue is soft yet strong, while the 2-ply, textured toilet tissue offers softness and durability, and is suitable for plumbing and septic systems. Green Seal Certified paper products are available per case only.

**ORDER CODE**

C-Fold Towel CF-1000 (10 1/8" x 13", 200 sheets/pkg) 12 per case - \$29.95

Facial Tissue FT-1000 (8" x 8", 100 count) 30 per case - \$40.95

Toilet Tissue BT-1000 (500 sheets per roll) 48 per case - \$41.95



## Patient Towels

The unique technique used to bond the layers of tissue and poly securely together eliminates separation of the layers and the reinforced edges provide added strength and durability. 2-ply tissue and 1-ply poly for absorption and patient protection.

**ORDER CODE**

PB-8000 White PB-8003 Green

PB-8001 Blue PB-8004 Lavendar

PB-8002 Mauve (Pink)  
13" x 18" - 500 Per Case - \$23.85



## 2" x 2" Cotton Filled Gauze Sponges

DEFEND has improved the performance of its cotton filled gauze sponges by adding 20% more cotton. Now more absorbent than before, these pure white, 8 ply sponges are shred-resistant and made of 100% cotton filled USP gauze. The 2" x 2" pads are available either sterile or non-sterile.

**ORDER CODE**

CS-0100 Non-Sterile - \$41.50

CS-0101 Sterile - \$57.45

5000 Per Case



## Non-Woven Sponges

A great alternative to cotton filled sponges, these are made from a 4-ply, 32 gram rayon/polyester blend which provides for maximum absorbency and little adhesion to wounds. These are soft, smooth, strong, non-linting and non-sterile.

**ORDER CODE**

NW-0200 (2" x 2" - 5000/cs) - \$29.75

NW-0300 (3" x 3" - 4000/cs) - \$42.35

NW-0400 (4" x 4" - 2000/cs) - \$38.90



## Cotton Rolls

These quality cotton rolls are made from 100% cotton, are non-linting, non-sterile, soft texture, pliable and 20% more absorbent than the leading brand. The virgin cotton holds together better, maintains its shape when wet and will not adhere to mucous membranes.

**ORDER CODE**

CS-0200 (#2, Medium 1.5" x .375") (3/8") - \$15.20  
2000 Box



## Headrest Covers

Made with soft and comfortable textured tissue outer covering and a poly backing to help prevent moisture soak through. The headrest covers are fluid-resistant and disposable.

**ORDER CODE**

HC-2001 (10" x 10") - \$31.90

HC-1001 (10" x 13") - \$35.95

White - 500 Per Case



## Disposable Cups

DEFEND disposable drinking cups have a special embossed grip that provides greater strength and makes dispensing easier. These recyclable cups have a lustrous, high-gloss finish and are available in five designer colors.

**ORDER CODE**

DC-7000 White

DC-7001 Blue

DC-7002 Mauve

5oz/148 ml - 1000 Per Case - \$26.80



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Procedure Accessories



### Pre-Bent Dispensing Tips

These dispensing tips are pre-bent to the proper angle for placing etchants, sealants, cements and other flowable materials. German made & final step washed and blown dry to avoid sediment within the metal tubes.

**ORDER CODE**

- PN-1227 Pink, Liners & Bases (18 gauge) - \$11.30
- PN-1228 Yellow, Sealant & Composites (20 gauge) - \$11.30
- PN-1229 Black, Flow Composites & Cements (19 gauge) - \$11.30
- PN-1231 Blue, Etchants (25 gauge) - \$11.30
- 100 Per Bag




### DEFEND Micro Applicators

DEFEND Micro Applicators are ideal for the application of bonding agents, etchants, hemostatic solutions, cavity liners and disclosing solutions. Its bendable neck "stays put" to facilitate the placement of materials into difficult-to-reach areas.

Non-absorbent, non-linting fibers suspend solutions as small as 1/8 of a drop, without dripping or spilling. DEFEND Micro Applicators are available in multiple colors to easily identify multi-step procedures.

**ORDER CODE**

- |                         |                         |
|-------------------------|-------------------------|
| <b>Superfine:</b>       | <b>Regular:</b>         |
| MB-1000 White - \$6.15  | MB-1100Blue - \$6.15    |
| <b>Fine:</b>            | MB-1101 Green - \$6.15  |
| MB-1001 Mauve - \$6.15  | MB-1102 Orange - \$6.15 |
| MB-1002 Yellow - \$6.15 | MB-1103 Purple - \$6.15 |
| 100 Per Tube            |                         |



### Surgical Aspirator Tips

Made with smooth edges for maximum patient comfort. The surgical aspirator is molded at a 30 degree angle and fits 11mm suction holders. This item is also autoclavable up to 135° C or may be disposed after single use. Available in 1/8" or 1/4" diameters.

**ORDER CODE**

- ST-1020 - Blue 1/16" Dia. - \$6.15
- ST-1021 - White 1/8" Dia. - \$6.15
- ST-1023 - Green 1/4" Dia. - \$6.15
- 25 Per Bag



### DEFEND Air/Water 3-Way Syringe Tips

DEFEND Air/Water 3-Way Syringe Tips are an excellent alternative to typical metal autoclavable tips. The disposable tips have an inner core water chamber that seats into the syringe, eliminating water contamination in the air chamber. The unique quad core design ensures instant dry air, prevents cross contamination, and delivers superior air/water spray.

**ORDER CODE**

- AW-1000 Clear (250/bag) - \$47.90
- AW-2000 White (250/bag) - \$47.90
- AW-3000 Blue (250/bag) - \$47.90
- AW-3500 Blue (1500/bag) - \$249.95



### Standard Saliva Ejectors

Fit all standard saliva ejector hose ends. Pliable, holds their shape and the tip minimizes tissue impingement.

**ORDER CODE**

- SE-7000 - White - \$4.10
- SE-7001 - Clear - \$4.10
- SE-7002 - Assorted Colors - \$4.10
- 100 Per Bag



### Vented Oral Evacuation Tips

These evacuators fit all standard evacuator hose ends and contain smooth edges for maximum patient comfort. The straight handle allows for increased control and are convenient single-use disposable tips. One end is vented and the other is non-vented (dual use). Non sterile, 1/8" diameter. Available in standard and short (3") lengths.

**ORDER CODE**

- VE-1121 - 6" Length - \$5.70
- VE-1123 - 3" Length - \$5.70
- 100 Per Bag

- DC-7003 Green
- DC-7004 Lavendar



## Scaler Inserts



### IF-50 Universal Series Insert

The IF-50 series is our universal sized tip with a plastic grip. This tip is best for general scaling procedures.

**TM0004-292 25K ART IF-50 #10 streamline, internal flow, plastic handle, ultrasonic insert (Blue grip) - \$69.00**

**TM0004-362 30K ART IF-50 #10 streamline, internal flow, plastic handle, ultrasonic insert (White grip) - \$69.00**



### IF-100 Slim Series Insert

The IF-100 series insert is a slim, internal flow, plastic grip insert. The slim tip is good for perio scaling of calculus on sensitive areas of the mouth.

**TM0004-452 ART 25K slim line series, internal flow, plastic handle, ultrasonic insert (Blue grip) - \$79.00**

**TM0004-462 ART 30K slim line series, internal flow, plastic handle, ultrasonic insert (White grip) - \$79.00**



### IF-5 BeaverTail Series

The IF-5 series insert is a beavertail, spatula-shaped insert. This tip is best for removing hidden calculus.

**TM0004-872 ART 25K IF-5 beavertail series 3, internal flow, plastic handle, ultrasonic insert (Blue grip) - \$69.00**

**TM0004-422 ART 30K IF-5 beavertail series 3, internal flow, plastic handle, ultrasonic insert (White grip) - \$69.00**



### P-10 Series Insert

The P-10 series tip is a universal sized tip with a metal grip. Water flows from an external water tube above the tip. This tip is best for general scaling procedures.

**TM0006-162 ART 25K P10 cavitron style, universal straight angled, external flow, metal handle - \$69.00**

**TM0007-062 ART 30K P10 cavitron style, universal straight angled, external flow, metal handle - \$69.00**



### P-100 Series Insert

The P-100 series is a slim, external flow insert with a metal grip. The slim tip is best for perio scaling of calculus on sensitive areas of the mouth. Water flows from an external tube above the tip. Straight (P-100), left (P-100L) and right-angled (P-100R) models available.

**TM0103-032 ART 25K P100 slim line furcation series, external flow, universal straight angled, metal handle - \$79.00**

**TM003-032 ART 25K P100L slim line furcation series, external flow, universal left angled, metal handle - \$79.00**

**TM003-012 ART 25K P100R slim line furcation series, external flow, universal right angled, metal handle - \$79.00**

**TM0103-052 ART 30K P100 slim line furcation series, external flow, universal straight angled, metal handle - \$79.00**

**TM0003-052 ART 30K P100 slim line furcation series, external flow, universal left angled, metal handle - \$79.00**

**TM0003-042 ART 30K P100 slim line furcation series, external flow, universal right angled, metal handle - \$79.00**

## Sabra Prophy Powder

### Sabra Prophy Powder

Sabra's Prophy polishing powder is totally compatible with your Prophy Jet as well as any other air-slurry polisher on the Market.

- No waste with simple screw top pour spout
- Great spearmint flavor
- No Clogging
- Save tons of money!

**ORDER CODE PP-03 - 3 Lb. Sabra Prophy Powder - \$47.30**

**ORDER CODE PP-06 - 6 Lb. Sabra Prophy Powder - \$83.60**

**ORDER CODE PP-12 - 12 Lb. Sabra Prophy Powder - \$156.00**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.



# Impression Products



### HP Mixing Tips

High Performance Mixing Tips for all HP style 50ml impression material cartridges. Tips are color coded according to size.

**ORDER CODE**

- VP-8105 Yellow (4.2mm)
- VP-8106 Pink (5.4mm)
- VP-8104 Teal (6.5mm)
- VP-8107 Purple (7.5mm)
- 48 Per Bag - \$35.50



### Clear Intra-Oral Tips

Tips provide precise placement interproximally. Fits with universal mixing tips, or with HP Mixing Tip #VP-8104.

**ORDER CODE**

- VP-8002 Large
- VP-8101SM Small
- 100 Per Bag - \$26.20



### Yellow Intra-Oral Tips

Tips provide precise placement interproximally. Fits 4.2mm mixing tips and to be used with HP Mixing Tip.

**ORDER CODE**

- VP-115 50 Per Bag - \$12.99




### Universal Mixing Tips

For Regular, Monophase and Heavy Body. Use with standard style cartridge systems, clear color.

**ORDER CODE**

- VP-8100 48 Per Bag - \$28.90
- ORDER CODE
- VP-8100 100 Per Bag - \$52.45



### Universal Monophase Mixing Tips

Designed for Vinyl Polysiloxane Impression Materials. Use with standard style cartridge systems, clear with purple mixer.

**ORDER CODE**

- VP-8120 48 Per Bag - \$28.90



### Ribbon Tips

For dispensing a ribbon-like layer of bite registration material on the occlusal surface. For use with #VP-8104. Easy snap on to Teal (green) HP Mixing Tip and to the Universal Mixing Tip.

**ORDER CODE**

- VP-8102 24 Per Bag - \$14.45



### Dynamic Mixing Tips

Ultra high quality mixing tips for use with automatic impression material mixing systems.

**ORDER CODE**

- VP-134 50 Per Bag - \$39.95



### Crown & Bridge Mixing Tips

Tips are high performance and for temporary crowns, bridges and acrylics. Available in 10:1, 4:1 and 1:1 ratios

**ORDER CODE**

- VP-8108 (10:1 & 4:1 ratio)
- VP-8111 (1:1 ratio)
- 25 Per Bag - \$19.90



### Temporary Cement Mixing Tips

Tips are high performance for Temporary cement. Available in 1:1 ratio

**ORDER CODE**

- VP-124 50 Per Bag - \$29.95

### Core Mixing Tip with Intra Nail Tip

**ORDER CODE**

- VP-125 50 Per Bag - \$32.95



### Dynamic Yellow Mixing Tips

Ultra high quality mixing tips for use with DEFEND Auto Mix System #VP-6001 & VP-6002.

**ORDER CODE**

- VP-8112 50 Tips Per Bag - \$84.95



### Bayonet Ring

Replacement pack for DEFEND Auto Mix System #VP-6001 & VP-6002.

**ORDER CODE**

- VP-8113 2 Rings Per Bag - \$13.95



### Dispensing Gun

The dispensing gun is used to extrude impression materials and bite registration materials from cartridges. It consists of a release lever and a plunger, which are coordinated simultaneously to load and remove the cartridge. The gun also consists of a trigger to dispense the material. Works with 1 to 1 and 2 to 1 cartridges. Available in two styles for HP cartridges and old style cartridges.

**ORDER CODE**


- VP-8202 For HP Cartridge (50ml) - \$66.65
- VP-8222 For Old Style Universal Cartridge (48ml) - \$47.20



Call Toll Free at 1-800-888-4435

# Impression Products

Empty or fully-loaded.



The BFC Syringe

The BFC Complete

**The choice is yours.**

## BFC Syringe

**Introducing the single use, disposable**

**Better Faster Cheaper Syringe**

The material wasted in a standard mixing tip (1.6cc) is enough to do 3 crowns with the BFC syringe. The small yet sensibly angled tip of the BFC gives easy access, and allows the material to be applied right at the sulcus minimizing bubbles and retakes. Save up to \$2.50 per impression using a tool that is easier and better than the existing technology. Simple.

- Pre-load and have them ready when you need them
- Better angle – narrow tip
- Lower cost
- No messy clean ups
- Minimizes bubbles
- Lower waste

**ORDER CODE BFC - \$72.99**

100 pcs

## BFC's Loaded Syringe

Pre-loaded with Vaccu-Sil, Light Body Impression material.

Compatible with all fast set heavy body material

**ORDER CODE BFC-PL - \$99.99**

100 per box

## Lock N'Reload™

**Reload your impression cartridges quickly and easily for double-barrel benefits.**

Save valuable time and up to 30% on the cost of VPS impression material by refilling partially used cartridges with the innovative Lock N'Reload™ coupler. Use one coupler and an extruder gun to refill dozens of cartridges per session. It's fast, easy to use and the revolutionary coupling design makes it impossible to cross mix the impression material.

- Fast & Easy. Anyone in the office can do it.
- Save over 30% by eliminating leftover waste.
- Eliminate wasted time spent switching cartridges on the fly.
- Use less mixing tips.

**ORDER CODE 72L - \$24.95**

20 Pk

## Occlusal Sweep

Disposable impression gun tip. Allows for smooth even disbursement material over the tooth, forcing material down into the gingival fold, eliminating bubbles.

**ORDER CODE 910CCL - \$18.95**

20 pcs

## Dexiter Tray Adhesive

Comparable to Bosworth. A paint on tray adhesive for use with alginates.

It helps to keep alginate firm in the tray during removal from the mouth.

Comes with convenient brush inside the cap.

**ORDER CODE DTA-01 - \$24.95**

1 oz.



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Impression Products



### The 5th Hand

The 5th Hand Lip and cheek retractor. Perfect tool for taking perfect impressions, as it keeps lips, cheeks, mustaches, etc completely out of the way. Retractor provides eschismic compression on salivary ducts, to cut down saliva production during impression taking. Also great retractor for photography, black surface does not reflect light and provides great contrast to the images.

#### ORDER CODE

- CR01B** - Full Arch lip and cheek retractor, 25 per bag, Small - **\$25.99**
- CR02B** - Full Arch lip and cheek retractor, 25 per bag, Medium - **\$25.99**
- CR03B** - Full Arch lip and cheek retractor, 25 per bag, Large - **\$25.99**
- CRSB** - Full Arch lip and cheek retractor, 25 per bag, Sampler Pack, 5 sm., 15 med., 5 lg. - **\$25.99**



### Cheek Retractor/Bite Blocks

Now available in bags of ten and molded in resilient plastic, our disposable lip/cheek retractor provides superior, patient-assisted retraction. Ideal access back to the second molar. The cheek retractor also doubles as a disposable bite block.

#### ORDER CODE

- 51SC** - Disposable Cheek Retractor/Bite Block, Small 10 per bag - **\$8.15**
- 51MC** - Disposable Cheek Retractor/Bite Block, Medium 10 per bag - **\$8.15**
- 51LC** - Disposable Cheek Retractor/Bite Block, Large 10 per bag - **\$8.15**




### Alginate Tray Cleaner

Remove stubborn alginate residues. Safe powder formula stores easily. Indispensable, efficient cleaner for plastic and metal trays. Easy to use. Just add water and soak.

#### ORDER CODE

- T734-1** - Alginate tray cleaner, 2 lbs powder - **\$28.39**
- T734-7** - Alginate tray cleaner, 5 lbs powder - **\$52.23**



### Posterior Dam & Periphery Extension Membrane -“Damit”

The Damit Post Dam Membrane is one of Sabra's most popular items. As the name suggests it closes off the back of the dental impression tray preventing the impression material from running down the back of the throat. It also forces the material up around the teeth and into the peripheral fold. Damit is a micropore membrane which becomes part of the dental impression. Placing dental wax across the back of an impression tray is a speed bump compared to damit which creates a fence. Use Damit around your 3-1, triple tray impressions to create support for the impression material. Damit is compatible with alginate, VPS and Ether impression materials. For pennies per application you will save on impression material, take more accurate impressions and be more comfortable for the patient.

- Captures second and third molars.
- Prevent gagging! Never have material run down the patients' throat again.
- Allow tray side removal with flexible extension to accommodate dental & lingual tori.
- Ideal for all viscosities of all impression materials.
- Saves impression material
- Increases accuracy of triple tray impressions
- Saves Money

#### ORDER CODE

- DamitDF**- 1 Roll with dispenser - **\$14.95**



# Impression Products & Organizers

## Sabra Bite Trays

Plastic disposable Bite trays with

Available in 5 sizes: varying quantities - please see below - **\$14.99 per bag**

### ORDER CODE


**SIT-EP01** - Posterior Bite Tray, 50 Trays/50 ID stickers

**SIT-EP02** - Sideless Posterior Bite Tray, 50 Trays/50 ID stickers

**SIT-EP03** - Full Arch Bite Tray, 25 Trays/25 ID stickers

**SIT-EP04** - Anterior Bite Tray, 30 Trays/30 ID stickers

**SIT-EP05** - Quadrant Bite Tray, 30 Trays/30 ID stickers



## Crystal Tray

Rigid Crystal disposable trays that are easily molded with any heat source and trimmed for a custom fit.

Available in packs of 12 perforated or solid - **\$7.15 per bag**

### ORDER CODE

| Solid       |                        | Perforated  |
|-------------|------------------------|-------------|
| <b>01S9</b> | Lower Right Quad       | <b>01P9</b> |
| <b>02S9</b> | Lower Left Quadrant    | <b>02P9</b> |
| <b>03S9</b> | Universal Anterior     | <b>03P9</b> |
| <b>04S9</b> | Small Upper Full Arch  | <b>04P9</b> |
| <b>05S9</b> | Medium Upper Full Arch | <b>05P9</b> |
| <b>06S9</b> | Large Upper Full Arch  | <b>06P9</b> |
| <b>07S9</b> | Small Lower Full Arch  | <b>07P9</b> |
| <b>08S9</b> | Medium Lower Full Arch | <b>08P9</b> |
| <b>09S9</b> | Large Lower Full Arch  | <b>09P9</b> |



## Crystal Bite Trays

3-1 bite registration trays are 100% injection molded, for superior asepsis and stability.

### ORDER CODE

**101B9** - Universal Anterior bite tray

**103B9** - Posterior Quadrant bite Tray

**105B9** - Full Arch Bite Tray

Packs of 12 - **\$7.15**



## Impression Organizer

Impression Tray Organizer ( Fits all Crystal Trays ) -

### ORDER CODE

**13GE** - Empty Grey Impression Organizer- **\$23.99**



## Composite Caddy

Join Organized dentistry today with the clear composite organizer. Each caddy holds 28 syringes and has 12 compartments to store all your loose items.

### ORDER CODE

**# 80PC** - Affordable Dental Composite Caddy- **\$29.95**



## Sabra Carousel

Holds up to 34 syringes, with multiple compartments, all on a " lazy Susan" type turn table - \$69.95 or FREE with the purchase of over \$150 worth of Sabra Composite Materials- Call 1-800-888-4435 to find out details.

### ORDER CODE

**SCC-01** - Sabra Composite/materials Carousel - **\$69.95**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
©2013 Sabra Dental Products, All Rights Reserved.



# Endodontics



For Best access to Apical Surgery - The Sabra OMS 45 II handpiece - See page 3. Also - See page 79 for Sabra Amalgam. Great 20:1 reduction handpiece for endo procedures - The SRC-Endo Handpiece - see page 5. For best perpendicular access to posteriorly treated canals - The Sabra 105 handpiece see page 2. See pages - 19 for the sabra 958SU endo bur.



### Root Canal Filling Material *Comparable to AH26/Dentsply*

Root canal filling and sealing material is used for permanent obturation of root canals of teeth of the secondary dentition with or without the aid of obturation points. Root canal powder and resin are mixed to produce a silver-free material with excellent sealing properties. As a result of its close application to the root canal walls and its very small contraction during setting, it allows for tightly sealed root canal obturation with excellent tissue compatibility.


**ORDER CODE DEX-022-010** - Root Canal Filling Material (Includes: 1 bottle 8g root canal filling powder, 1 tube 10g root canal filling resin, scoop, and pad) - **\$59.99**



### Hema Desensitizer *Comparable to Gluma Desensitizer/ Kulzar*

Contains hema and benzalkonium chloride. Fast, effective way of sealing off dentinal tubules by combining with proteins and forming a physiologic barrier of protection. Benzalkonium chloride adds an antimicrobial action to eliminate bacteria. Used as a normal part of resin bonding. Also effective for cervical sensitivity, root planning, and scaling.


**ORDER CODE HEMA** Hema Desensitizer 7ml - **\$39.99**



### Gutta Percha *Made from the finest Gutta Percha*

Uniformly hand rolled, Stiff but flexible to provide ideal working qualities. Uniform and non-distorting crimp, Excellent radiopacity, long shelf life and freshness assured. Suitable for vertical and lateral condensation, Solvent procedures, Variety in the selection and sizes of the packaging. Available in ISO standardized & non standardized. X-fine, fine-fine, medium-fine, fine, fine-medium, medium, medium-large, large, x-large.

**ORDER CODE 8-80** Gutta Percha - .04 Taper, .06 Taper - **\$7.99** Reg 120pts - **\$4.99**



### Absorbent Paper Points

Rigid enough to insert in canals without curling or bending, Rolled firmly to hold shape when withdrawn from canal, Highly absorbent, Sterilized individual packages. Variety in the selection of sizes and packaging available.

**ORDER CODE 8-80** Absorbent paper points - .04 Taper, .06 Taper - **\$7.99** Reg 200pts - **\$3.99**



### Formo Creosol

Used to devitalize and sterilize infected pulp chambers - 1oz.

**ORDER CODE DEX12-660** Formo Creosol - **\$24.99**



### Eugenol USP Grade

Has analgesic and antiseptic qualities

**ORDER CODE DEX12-651** Eugenol USP 1oz - **\$7.99**



### Endo/Implant Ruler

This limited-use 40mm clear Endodontic Implant ruler is available in bags of 10. Each bag contains ten 40mm rulers and finger clips. Ideal for measuring intraoral lesions, spaces and images on X-Rays. Because it's clear, you can lay it over your object or image for greater accuracy.

**ORDER CODE 40mm Clear Disposable Endo/Implant Rulers 10per** - **\$9.65**

*Please see [www.sabradent.com](http://www.sabradent.com) for the entire line of EDS products at great prices.*



Call Toll Free at 1-800-888-4435

# Amalgams



## Nogama Dispersed Phase Alloy 69.5% AG (Similar to Dispersalloy)

Nogama is an admix of spherical, copper-silver eutectic particles and lathe cut conventional alloy

### ORDER CODE

- 110471 **Silver Shining surface**, Nogama Box of 50, (1 Spill) - **\$56.00**
- 110472 **Supreme marginal adaptability**, Nogama Box of 50, (2 Spill) - **\$76.00**
- 110473 **Excellent durability**, Nogama Box of 50, (3 Spill) - **\$99.00**
- 110491 **Firm pack & smooth packing**, Nogama Box of 500, (1 Spill) - **\$550.00**
- 110492 **Nogama Box of 500**, (2 Spill) - **\$750.00**
- 110493 **Nogama Box of 500**, (3 Spill) - **\$999.00**

### Physical Properties

Composition . . . . . Ag 69.5%, Sn 19.5%, Cu 10.5%, ZN 0.5%  
 Particle shape . . . . . 80% lathe cut 20% spherical  
 Creep (7 days) . . . . . 0.35%  
 Dimensional Change . . . . . +0.80%  
 Compressive Strength 1 hr. . . . . > 164 MPa  
 Compressive Strength 24 hr. . . . . > 470 MPa



## Spherodon-M Admix Alloy 45.5% AG (Similar to Contour)

Spherodon M provides excellent stability, enhanced corrosion resistance and reduced ditching

- 120471 **Outstanding dimensional**, Spherodon-M Box of 50, (1 Spill) - **\$42.00**
- 120472 **Stability**, Spherodon-M Box of 50, (2 Spill) - **\$58.00**
- 120473 **Improved marginal integrity**, Spherodon-M Box of 50, (3 Spill) - **\$73.00**
- 120491 **More body and easy**, Spherodon-M Box of 500, (1 Spill) - **\$399.00**
- 120492 **Condensation versatility**, Spherodon-M Box of 500, (2 Spill) - **\$575.00**
- 120493 **Spherodon-M Box of 500**, (3 Spill) - **\$725.00**

### Physical Properties

Composition . . . . . Ag 45.5%, Sn 31.5%, Cu 23%  
 Creep (7 days) . . . . . 0.35%  
 Dimensional Change . . . . . +0.08%  
 Compressive Strength 1 hr. . . . . > 180 MPa  
 Compressive Strength 24 hr. . . . . > 460 MPa



## Spherodon+ 100% spherical Alloy 60% AG (Similar to Tytin)

Spherodon+ is a single composition ultra fine spherical ternary alloy

- 130471 **Superb compressive**, Spherodon+ Box of 50, (1 Spill) - **\$49.00**
- 130472 **Tensile strength**, Spherodon+ Box of 50, (2 Spill) - **\$72.00**
- 130473 **Decreased micro leakage**, Spherodon+ Box of 50, (3 Spill) - **\$94.00**
- 130491 **Extra low mercury alloy ratio**, Spherodon+ Box of 500, (1 Spill) - **\$499.00**
- 130492 **Spherodon+ Box of 500**, (2 Spill) - **\$699.00**
- 130493 **Spherodon+ Box of 500**, (3 Spill) - **\$939.00**

### Physical Properties

Composition . . . . . Ag 60.5%, Sn 28%, Cu 11.5%  
 Creep (7 days) . . . . . 0.62%  
 Dimensional Change . . . . . +0.01%  
 Compressive Strength 1 hr. . . . . > 171 MPa  
 Compressive Strength 24 hr. . . . . > 485 MPa



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Cosmetic Dentistry



### Light cure core build up material

Comparable to Clearfill photocure/Kuraray

A highly filled composite specifically formulated for core build ups. Available in blue or natural A2 shade. The distinctive blue color aids in distinguishing the core material from tooth structure. A 6-7mm depth of cure can be obtained. The blue core is best suited under metal crowns, the natural shade for all other applications. Radiopaque for easy identification in radiographs and has high viscosity for easy build up with no slumping.

VLC core material 5 gram syringe natural shade or contrast blue shade.

**ORDER CODE**

DEX-003-022 (Natural) - \$32.99

DEX-003-021 (Blue) - \$32.99



### Chemical cure core build up material

Comparable to compcore/premier or encore/centrix

A highly radiopaque, high compressive strength material. The non slumping consistency allows the material to be syringed directly into the preparation. The viscosity is ideal for pre fabricated post cementation and crown and bridge relining. Easily distinguished from tooth structure by its distinctive white, blue or natural color.

14 gram catalyst and base core build up material.

**ORDER CODE**

DEX-003-018 - \$59.99



### Anesthetic Gel

ORDER CODE - \$5.99 per

STA-B Topical Anesthetic, 1oz, Bubble Gum

STA-C Topical Anesthetic, 1oz, Cherry

STA-M Topical Anesthetic, 1oz, Mint

STA-P Topical Anesthetic, 1oz, Pina Colada

STA-S Topical Anesthetic, 1oz, Strawberry



### Profil Temp

ORDER CODE

PFJar-01 PROFIL TEMP, Temporary Cavity Filling Material, 28g jar - \$14.99

PFTube-01 PROFIL TEMP, Temporary Cavity Filling Material, 4 x 7g tubes - \$19.99



### ProShield

Scientifically designed to block the harmful effects from ultra-violet during composite restoration procedures. The ProShield is the only protective gear in the market that really protect your eyes. It also blocks the LED head light from prematurely curing the composite material. The ProShield Handheld is large enough to block off all harmful light for both the dentist, the assistant and the patient. The specially designed angled ProShield adapts to the operating field. The only shield in the market that really protects your eyes.

ORDER CODE:

#SPS-01 - \$34.95



# Handpiece Maintenance

## Sabra-Care Plus Lubricant Kit



12oz. with quick connect nozzle for Sabra handpieces  
 Performance rating 5-10 times better than standard lubricants  
 (red straw included for 4 hole handpieces)

Forms a protective film around metal parts, Lowers operating temperature, reducing bearing wear, Cleans as it lubricates,  
 Preventing contaminant build up, Fewer repairs,  
 Save Time - Save Money

**ORDER CODE SLK01 - \$38.00**

**ORDER CODE SLK02 Sabra-Care Plus lubricant refill 12fl oz. ( no nozzle included ) - \$32.00**



## Sabra Oil Kit


1x 4 oz bottle, 2 10cc bottles  
**ORDER CODE SLL-01 - \$34.00**



Lube nozzle for Sabra quick connect handpieces  
**ORDER CODE SLK-03 - \$9.00**



Sabra Oil, 5cc  
**ORDER CODE SLO01 - \$6.00**



## E-Type Spray Nozzle

**ORDER CODE SLK-04 Sabra E- Type Lubricant Nozzle - \$9.00**

## Handpiece Maintenance Purge Station



- ORDER CODE SC-1875 Sabra handpiece purge station with Hepa filters - \$229.00**
- ORDER CODE SC1883 Sabra adaptor for Sabra handpiece purge system (Also fits NSK, Kinetics, Viper, and Healthco handpieces ) - \$49.95**
- ORDER CODE SC-1880 Midwest 4 hole adaptor for Sabra handpiece purge system - \$39.95**
- ORDER CODE SC-1881 Star adaptor for Sabra handpiece purge system - \$39.95**
- ORDER CODE SC-1882 Kavo adaptor for Sabra handpiece purge system - \$54.95**
- ORDER CODE SC-1885 2 hole borden adaptor for Sabra handpiece purge system - \$39.95**
- ORDER CODE SC-1887 Midwest XGT adaptor for Sabra handpiece purge system - \$54.95**
- ORDER CODE SC-1888 Lares adaptor for Sabra handpiece purge system - \$42.99**
- ORDER CODE SC-1889 W&H adaptor for Sabra handpiece purge system - \$42.95**
- ORDER CODE SC-1K Replacement bottle-less cap - \$10.95**
- ORDER CODE SC-2690 Replacement Hepa filter for Sabra purge system - \$10.95**



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)  
 ©2013 Sabra Dental Products, All Rights Reserved.



## Turbines and Replacement Parts



### Replacement Turbines

#### ORDER CODE

- ST-II R turbine for Sabra OMS 45-II, Sabra OMS 45-FO, Sabra 105, Sabra 105 FO - **\$209.00**
- ST-12 Replacement Turbine for Sabra 45-LS - **\$209.00**
- ST-ETPB Sabra Econotorque Push Button Replacement Turbine- **\$109.00**
- ST-ET Sabra Econotorque Replacement Turbine, Jacobs- **\$89.00**
- STOMSD Sabra OMS-45 Replacement Turbine (For original OMS-45, handpiece discontinued) - **\$209.00**
- STBD Sabra Black Diamond/Topair Replacement Turbine - **\$209.00**

### Replacement Back Caps

#### ORDER CODE

- SP-106 Head Cap for Sabra OMS-45 II/FO/105 - **\$42.99**
- SP-107 Head Cap, Sabra Econotorque - **\$14.99**
- SP-108 Replacement Head Cap, Sabra Econotorque Push Button - **\$24.99**



### Replacement O-Rings

#### ORDER CODE

- SP-105 O-Ring Kit - Sabra OMS-45 II/FO/105 (Set of 5) - **\$9.00**
- SP-QDJ Set "O" Rings for QD-J M4/ QD-J B2 Coupling- **\$9.00**



### Misc Parts

#### ORDER CODE

- SP-108 Head Cap Wrench, Sabra Handpieces ( Series OMS, 105 ) - **\$12.00**
- SP-109 Bur Wrench, Sabra Econotorque - **\$15.95**
- SP-099 Standard 2 Hole Gasket - **\$4.00**
- SP-102 Standard 4 Hole Gasket - **\$4.00**
- SP-103 Standard 5 Hole Gasket - **\$4.00**
- SP-104 Standard 6 Hole Gasket - **\$4.00**
- SP-107 Anti-Retraction Valve for Sabra C-4, C-5, C-6 Coupler - **\$42.99**



# Sabra Rapid Repair Services



## Sabra Rapid Repair Service

Sabra Dental is a full service repair facility for your handpiece repair needs. We offer quick turn-around times and competitive pricing on high-speeds, low speeds, electric units, lab motors, surgical drills etc.

Sabra Dental is the Authorized repair facility for the New York County Dental Society.

Contact us at **1-800-888-4435** for a Free postage paid handpiece mailing box or simply go to [www.sabradent.com](http://www.sabradent.com) and print a mailing label. Place it on your own mailer and send it directly to us --- and we'll pick up the shipping costs.

**Sabra repairs all makes and models including:** Sabra, NSK, Kavo, Midwest, Star, Kinetics, W&H, Sirona, Schein and other dealer house brands, MKDent, Micro-Motors, Lares, Palisades, Dabi Alante, Lynx, MTI, Athena/Champion, Medidenta, Aseptico/Anthogyr, Sybron, Hall, Stryker, most small equipment (scalars, and prophy units), and many brands of implant/surgical units.

**We Offer!** 24 hour/same day repairs on most handpieces.

**4x Sabra now offers:**  
4x month warranty on high speed rebuilds  
6 months warranty on low speed overhauls

### Sabra Small Equipment Repair Take a closer look at what we fix!

#### Curing Lights

Demetron, Patterson, Healthco,  
Degussa, Dentsply/Caulk, 3M

#### Plasma Arc Curing Lights

DMD, ADT, Wave Light, Arc Light, Rolence

#### Ultrasonic Cleaners

L&R, Healthco, Patterson, Health Sonic

#### Electrosurge Units

Coles, Parkel, Healthco, Macan, Strobex, Hampton Research, Elman

#### Ultrasonic Scalars

Dentsply, Sonatron, Parkell, Sonus, Vetrason, Microson

#### Amalgamators

Caulk, Crescent, Degussa, SDI, Adec, Kerr, Espe

### CALL NOW!

for a postage paid handpiece repair mailer

# 1-800-888-4435

or VISIT OUR WEBSITE

to print a FREE return repair label

## [www.sabradent.com](http://www.sabradent.com)

289 Suburban Ave. Suite E  
Deer Park, N.Y. 11729



Toll Free: 800.888.4435  
[www.sabradent.com](http://www.sabradent.com)



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# “Top 10” OSHA Compliance Products



### OSHA Compliance System

The comprehensive HPTC OSHA Compliance Program covers all aspects of the OSHA Act.

#### Complete Program Includes:

- Written training manual •Implementation DVD or Video •Hazardous Materials Wall Chart
  - Implementation instructions •Compliance label packet •6 CEU's for three individuals following exam submission
  - Exposure control plan •Toll free Technical Hot Line •TB & Ergonomics section
  - Upgrade service (1 year) • Chemical MSDS
- Upgrade Available. Call for more information.

**ORDER CODE DCPE - \$499.00**



### Master Spill Kit

This all purpose kit combines absorbents for all of the most common spills that occur in a healthcare facility. **Contents include:**

- Biological spill powder, 1 pair nitrile gloves
- MercoSorb, 1 Mercury magnet, Chemabsorb, 1 Mercury scrubber, 1 Scooper & pan, Biohazard disposal bags, 1 Pair safety goggles, Biohazard labels

**ORDER CODE MSK - \$136.95**



### Compliance Label Kit

Various sized chemical, biohazard and radiation caution labels are combined with a hazardous materials wall chart and labeling instructions to assist your office in OSHA requirements for labeling. Toll free technical hot line service included.

**ORDER CODE CLK - \$129.95**



### In-Office OSHA Assessment Video Program IOA/VP w/ DVD

An annual “Certified Written Hazard Assessment” is required by OSHA regulations for every workplace.

Our complete program contains everything you will need to achieve compliance, including:

- OSHA regulations
- compliance checklist
- instructional video
- TB risk assessment
- Certificate of Assessment
- HPTC Toll-free Technical Hotline Service

Completion of training exams yields 2 CEU's for 3 individuals.

**ORDER CODE IOA/VP - \$86.50**



### Bloodborne Compliance Wallchart

Designed to educate and encourage employee compliance with OSHA'S Bloodborne Disease Standard, this chart employs pictorials and written text in a clear and concise format. This laminated poster should be displayed in all workplaces to detail measures taken to protect your employees.

17" x 31" poster

**ORDER CODE BCWC - \$39.95**



### Labor Law Poster

Federal law mandates that the Labor Law Poster must be displayed in the workplace in view of all employees. Employee information notices required are:

- Fair Labor Standard Act
- Employee Polygraph Protection
- Family Medical Leave Act
- Age Discrimination in Employment Act
- Title VII of the Civil Rights Act
- OSHA Form 3165
- American Disability Act

This comprehensive all-in-one laminated poster includes all required information.

**ORDER CODE LLP - \$44.95**

See [www.sabradent.com](http://www.sabradent.com) for the complete line of Sabra OSHA products



Call Toll Free at 1-800-888-4435

# OSHA Compliance



## Emergency Eyewash Station

This unit fits any faucet equipped with a removable aerator, and features an automatic shut-off valve.

With Eliminator Valve

ORDER CODE EESE - \$116.95



## Scalpel Blade Remover

This "Blade Mate" remover is designed to eliminate scalpel blade cuts. The blades drop safely into the sealed reservoir (stores 15-20 blades). The reservoir acts as an autoclavable sharps container.

6+ Units ea.

ORDER CODE SBR6 - \$32.95



## X-Ray Monitoring Badge

Experts agree that there is no safe level of radiation. Even the smallest dose may cause cancer and genetic damage. Any healthcare facility may be exposing staff to unsafe levels of radiation. That's why accurately measuring and recording radiation exposure over the long term is important for your people and your practice. The Landauer Luxel Plus with state-of-the-art Optically Stimulated Luminescence (OSL) technology, has the highest sensitivity available today -10 times the sensitivity of film badges!

### Features include:

- Fastest reporting in the industry- five day average turnaround
- A lifetime exposure record is archived for safekeeping
- Personalized badges
- Web-based reporting

Quarterly Badge Service  
ORDER CODE XMBQ - \$99.99



## Emergency Signs-Glow in The Dark-Emergency Exit Only

OSHA requires that all "emergency exits," "exits," and "non-exits" be marked with clear, visible, illuminated signs. These light weight signs glow in the dark for up to 8 hours and are 10"x7". The easy mount fire extinguisher sign measures 4"x18". All signs meet OSHA & NFPA requirements.

'Emergency Exit Only'

ORDER CODE EOS - \$53.95



## Emergency Signs-Glow in The Dark-Not and Exit

OSHA requires that all "emergency exits," "exits," and "non exits" be marked with clear, visible, illuminated signs. These light weight signs glow in the dark for up to 8 hours and are 10"x7". The easy mount fire extinguisher sign measures 4"x18". All signs meet OSHA & NFPA requirements.

'Not An Exit'

ORDER CODE NXS - \$53.95



## Emergency Signs-Glow in The Dark-Exit

OSHA requires that all "emergency exits," "exits," and "non exits" be marked with clear, visible, illuminated signs. These light weight signs glow in the dark for up to 8 hours and are 10"x7". The easy mount fire extinguisher sign measures 4"x18". All signs meet OSHA & NFPA requirements.

'Exit'

ORDER CODE EXS - \$53.95



## Emergency Signs-Glow in The Dark-Fire Extinguisher

OSHA requires that all "emergency exits," "exits," and "non exits" be marked with clear, visible, illuminated signs. These light weight signs glow in the dark for up to 8 hours and are 10"x7". The easy mount fire extinguisher sign measures 4"x18". All signs meet OSHA & NFPA requirements.

'Fire Extinguisher'

ORDER CODE FES - \$59.95

See [www.sabradent.com](http://www.sabradent.com) for the complete line of Sabra OSHA products



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.



# Notes



# Notes



Visit Our Web Site at [www.sabradent.com](http://www.sabradent.com)

©2013 Sabra Dental Products, All Rights Reserved.

# Notes





[WWW.SABRADENT.COM](http://WWW.SABRADENT.COM)  
1-800-888-4435